

## Product datasheet for RC232610

### KCNJ15 (NM\_001276437) Human Tagged ORF Clone

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | Expression Plasmids   |
| Product Name:             | KCNJ15 (NM_001276437) Human Tagged ORF Clone                      |
| Tag:                      | Myc-DDK   |
| Symbol:                   | KCNJ15  |
| Synonyms:                 | IRKK; KIR1.3; KIR4.2  |
| Mammalian Cell Selection: | Neomycin  |
| Vector:                   | pCMV6-Entry (PS100001)  |
| E. coli Selection:        | Kanamycin (25 ug/mL)  |
| ORF Nucleotide Sequence:  | >RC232610 ORF sequence<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGATGCCATTCACATCGGCATGTCCAGCACCCCTGGTGAAGCACACTGCTGGGGCTGGGCTCAAGG  
CCAACAGACCCCGCTCATGTCCAAGAGTGGGCACAGCAACGTGAGAATTGACAAAGTGGATGGCATATA  
CCTACTCTACCTGCAAGACCTGTGGACCACAGTTATCGACATGAAGTGGAGATACAACTCACCTGTTC  
GCTGCCACTTTTGTGATGACCTGGTTCCTTTTGGAGTCATCTACTATGCCATCGCGTTTATTCATGGGG  
ACTTAGAACCCGGTGAGCCATTTCAAATCATACCCCTGCATCATGAAAGTGGACTCTCTCACTGGGGC  
GTTTCTCTTTCCCTGGAATCCCAGACAACATTGGCTATGGAGTCCGTTCCATCACAGAGGAATGTCT  
CATGCCATCTTCTGTTGGTTGCTCAGTTGGTCATCACGACCTTGATTGAGATCTTCATCACCGGAACCT  
TCTGGCCAAAATCGCCAGACCCAAAAGCGGGCTGAGACCATCAAGTTCAGCCACTGTGCAGTCATCAC  
CAAGCAGAATGGGAAGCTGTGCTTGGTATTAGGTAGCCAATATGAGGAAGAGCCTCTTGATTCAAGTGC  
CAGCTCTCTGGCAAGCTCCTGCAGACCCACGTACCAAGGAGGGGGAGCGGATTCTCTCAACCAAGCCA  
CTGTCAAATCCACGTGGACTCCTCCTGAGAGCCCCCTCTCATTCTGCCATGACATTCTACCATGT  
GCTGGATGAGACGAGCCCCCTGAGAGACCTCACACCCAAAACCTAAAGGAGAAGGAGTTTGAGCTTGTG  
GTCCTCCTCAATGCCACTGTGGAATCCACCAGCGCTGTCTGCCAGAGCCGAACATCTTATATCCAGAG  
AAATCTACTGGGTTTTGAGTTTGTGCCTGTGGTATCTCTCTCCAAAAATGGAAAATATGTGGCTGATT  
CAGTCAGTTTGAACAGATTTCGAAAAGCCAGATTGCACATTTTACTGTGCAGATTCTGAGAAAACAGCAA  
CTCGAGGAGAAGTACAGGCAGGAGGATCAGAGGAAAAGAGAAGTGAAGACTTTTATTACAACAGAGCA  
ATGTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC232610 protein sequence  
Red=Cloning site Green=Tags(s)

MDAIIHIGMSSTPLVKHTAGAGLKANRPRVMSKSGHSNVRIDKVDGIYLLYLQDLWTTVIDMKWRYKLTFL  
 AATFVMTWFLFGVIYYAIAFIHGDLPEGEPISNHTPCIMKVDSL TGAFLFSLESQTTIGYGVRSITEECP  
 HAIFLLVAQLVITLIEIFITGTFLAKIARPKKRAETIKFSHCAVITKQNGKLC LVIQVANMRKSLLIQC  
 QLSGKLLQTHVTKEGERILLNQATVKFHVDSSESPFLILPMTFYHVLDETSPLRDLTPQNLKEKEFELV  
 VLLNATVESTSAVCQSRYSYIPEEIIYWGFEFVPVVSLSKNGKYVADF SQFEQIRKSPDCTFYCADSEKQQ  
 LEEKYRQEDQRERELRLLLLQQSNV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk6575\\_c12.zip](https://cdn.origene.com/chromatograms/mk6575_c12.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_001276437

**ORF Size:** 1125 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001276437.2](#)

**RefSeq Size:** 2903 bp

**RefSeq ORF:** 1128 bp

**Locus ID:** 3772

**UniProt ID:** [Q99712](#)


**Cytogenetics:** 21q22.13-q22.2

**Protein Families:** Druggable Genome, Ion Channels: Potassium, Transmembrane

**MW:** 42.6 kDa

**Gene Summary:** Potassium channels are present in most mammalian cells, where they participate in a wide range of physiologic responses. The protein encoded by this gene is an integral membrane protein and inward-rectifier type potassium channel. The encoded protein has a greater tendency to allow potassium to flow into a cell rather than out of a cell. Eight transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Feb 2013]

### Product images:



Circular map for RC232610