

Product datasheet for **RC232638**

Transcription termination factor 1 (TTF1) (NM_001205296) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Transcription termination factor 1 (TTF1) (NM_001205296) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TTF1
Synonyms:	TTF-1; TTF-I
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC232638 representing NM_001205296 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTACCGGGACGACTTGAACGGTTTAAAGGAATTTAAAGCACAAAGGTGTCGCTATTAATTTGGCAAGT
TTTCTGTAAGGAAAATAAGCAGTTAGAGAAAAATGTGGAAGACTTTCTAGCCCTGACAGGCATTGAGAG
TGCAGACAAGCTGCTGTACACGGACAGATATCCTGAGGAAAAATCTGTGATCACCAACTAAAAAGGAGA
TACTCGTTTAGATTACACATTGGTAGGAACATTGCCCGGCCCTGGAACTTATATACTATCGAGCAAAGA
AGATGTTTCGATGTCAACAATTACAAAGGCAGGTATAGCGAAGGAGATACTGAGAAGTTAAAGATGTACCA
TTCTCTCCTTGGGAATGACTGGAAGACGATTGGTGAGATGGTGGCCCGAAGTAGCCTCTCCGTGGCCCTC
AAGTTCACAGATCAGCAGTCAAAGAAATCGTGGTGCTTGGAGTAAGTCTGAAACCCGGAACTAATCA
AGGCTGTCGAAGAAGTGATTCTGAAGAAGATGTCTCCCCAGGAGTTAAAAGAGGTGGATTCCAACTCCA
AGAAAACTCTGAAAGTTGCCTATCAATTGTTCCGGAAAACTCTACAAGGCATATCTTGGGTAGAAGTA
GAAGCTAAAGTGCAAACCAGAAATTGGATGCAGTGTAAAAGTAAGTGGACAGAAATCTAACCAAGAGGA
TGACTAATGGTCGGCGTATCTACTATGGCATGAATGCCCTGCGGGCCAAGGTCAGCCTTATTGAAAGGTT
GTATGAAATAAATGTGGAAGATACTAATGAAATAGACTGGGAAGATCTTGCTAGTGCCATAGGTGATGTT
CCTCCATCTTACGTTCAAATAAATTTCTAGGCTGAAAGCTGTCTATGTTCCATTTTGGCAGAAAAAGA
CTTTTCCAGAGATCATCGACTACCTTTATGAGACGACTCTACCTTTGCTGAAGGAAAAGTTAGAAAAAT
GATGGAGAAAAAGGCACTAAAAATCCAGACTCCTGCAGCACCCAAGCAAGTTTTCCCATTTCCGAGACATC
TTTTATTATGAAGACGATAGTGAAGGAGAGGACATAGAAAAAGAAAGCGAAGGCCAGGCCATGCATGG
CTCACGCTGTAATTCAGTACTTTGGGAGCCAAGGCCGGTGGATCATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC232638 representing NM_001205296
Red=Cloning site Green=Tags(s)

MYRDDLERFKEFKAQGVAIKFGKFSVKENKQLEKNVEDFLALTGIESADKLLYTDRYPEEKSVITNLKRR
 YSFRLHIGRNIARPWKLIYYRAKKMFDVNNYKGRYSEGDETEKLMYHSLLGNDWKTIGEMVARSSLSVAL
 KFSQISSQRNRGAWSKSETRKLKAVEEVILKKMSPQELKEVDSKLQENPESCLSIIVREKLYKGISWVEV
 EAKVQTRNWMQCKSKWTEILTKRMTNGRRIYYGMNALRAKVSLIERLYEINVEDTNEIDWEDLASAIGDV
 PPSYVQTKFSRLKAVYVPFWQKKTFPEIIDYLYETTLLPLLEKLEKMEKKGTKIQTPAAPKQVFPFRDI
 FYYEDDSEGEDIEKESEGQAPCAHACNSSTLGGQGRWII

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001205296

ORF Size:	1170 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001205296.1 , NP_001192225.1
RefSeq Size:	1780 bp
RefSeq ORF:	1173 bp
Locus ID:	7270
UniProt ID:	Q15361
Cytogenetics:	9q34.13
MW:	46.1 kDa
Gene Summary:	This gene encodes a transcription termination factor that is localized to the nucleolus and plays a critical role in ribosomal gene transcription. The encoded protein mediates the termination of RNA polymerase I transcription by binding to Sal box terminator elements downstream of pre-rRNA coding regions. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. This gene shares the symbol/alias 'TFF1' with another gene, NK2 homeobox 1, also known as thyroid transcription factor 1, which plays a role in the regulation of thyroid-specific gene expression. [provided by RefSeq, Apr 2011]