

Product datasheet for **RC232644**

ANKS1B (NM_001204067) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKS1B (NM_001204067) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ANKS1B
Synonyms:	AIDA; AIDA-1; ANKS2; cajalin-2; EB-1; EB1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC232644 representing NM_001204067 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGACCCAGATGCCCTGTCCAACAGTGGGACAATGGTTGGAAAGCATTGGGCTACCTCAGTACGAGA
ACCACCTGATGGCTAATGGATTTGACAATGTGCAGTTTATGGGAAGCAATGTTATGGAAGATCAGGATTT
GTTGGAAATTGGAATCCTTAATTCTGGGCACAGACAAAGAATTCTACAGGCAATCCAGCTCCTTCCAAG
ATGAGACCCATTGGGCATGATGGCTACCATCCACCTCTGTAGCTGAGTGGCTGGATTCCATTGAACTGG
GCGACTACACCAAGCCTTTCTAATTAATGGCTACACTTCGATGGACCTGTTGAAAAAATCTGGGAGGT
TGAACCTATTAATGTTTTAAAAATCAATTTGATTGGCCACAGGAAACGTATTTGGCATCTCTGGGAGAC
AGGCTGCACGACGATCCCCACAGAAGCCCCCTCGGTCCATCACCTCAGGACAGGAGACTGGGGAGAAC
CTTCCATTACCTTGCACCTCCGAATGAAGCCACAGCCTCTACCCCGGTACAGTACTGGCAGCATCACCC
AGAAAAGCTTATCTTCCAGTCGTGTGATTACAAAGCTTTTTATTTAGGTTCTATGCTGATAAAAGAGCTT
AGGGGGACAGAATCAACCAAGATGCTTGTGCAAAAATGCGGGCTAACTGTCAGAAGTCTACAGAGCAA
TGAAGAAGTCCCTACTATTATTTCTTCTCATATAAAGGAGTCAAATTTATTGATCAACAATAA
GAACATAATTGCTGAGCATGAAATTCGTAATATCTCCTGTGCTGCCAGGACCCAGAAGACCTCTCAACA
TTTGCCTATATCACAAAAGATTTGAAGTCTAATCACCCTACTGTGATGTTTACTGCCTTTGATGTGA
ATTTAGCCTATGAAATCATCCTAACCTGGGACAGGCATTGGAAGTCGCTTACCAGCTAGCACTACAAGC
AAGAAAAGGGGACACTCCTCCACACTTCCAGAAAGCTTTGAAAACAAACCCTCCAAACCATCCCCAAG
CCCCGCTTAGCATTGCAAGTCCGTGATCGACCATCTGAGCAAAGACTCTGGCCAATCTACCGTGA
TTGTGGAGCCGGCCAAGAAGCCAAGAGGGCATTAAATACCAAGTATGAAACCAGATTTTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC232644 representing NM_001204067
 Red=Cloning site Green=Tags(s)

MGPRCPVQTVGQWLESIGLPQYENHLMANGFDNVQFMGSNVMEDQDLLEIGILNSGHRQRILQAIQLLPK
 MRPIGHGTYHPTSAEWLDSIELGDYTKAFLINGYTSMDLLKKIWEVELINVLKINLIGHRKRILASLGD
 RLHDDPPQKPPRSITLRTGDWGEPSITLRPPNEATASPVOYQWHHPEKLIHQSCDYKAFYLGSMLIKEL
 RGTESTQDACAKMRANCQKSTEQMKKVPITIIILSVSYKGVKFFIDATNKNIIEAHEIRNISCAAQDPEDLST
 FAYITKDLKSNHHYCHVFTAFDNLAYEIIILTLGQAFEVAYQLALQARKGGHSSSTLPESFENKPSKPIPK
 PRVSIRKSVIDPSEQKTLANLPWIVEPGQEAARGINTKYETTIF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001204067

ORF Size:	1182 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001204067.1 , NP_001190996.1
RefSeq Size:	3157 bp
RefSeq ORF:	1185 bp
Locus ID:	56899
UniProt ID:	Q7Z6G8
Cytogenetics:	12q23.1
MW:	45 kDa
Gene Summary:	This gene encodes a multi-domain protein that is predominantly expressed in brain and testis. This protein interacts with amyloid beta protein precursor (AbetaPP) and may have a role in normal brain development, and in the pathogenesis of Alzheimer's disease. Expression of this gene has been shown to be elevated in patients with pre-B cell acute lymphocytic leukemia associated with t(1;19) translocation. Alternatively spliced transcript variants encoding different isoforms (some with different subcellular localization, PMID:15004329) have been described for this gene. [provided by RefSeq, Aug 2011]