

Product datasheet for **RC232709**

ZNF410 (NM_001242926) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF410 (NM_001242926) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF410
Synonyms:	APA-1; APA1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC232709 representing NM_001242926
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGTTATCAGATGAGTTAGAATCCAAACCAGAGCTCCTGGTACAGTTTGTTCAGAATACGTCCATCCCAT
 TGGGACAGGGGCTTGTAGAATCAGAAGCTAAAGATATTACTTGCTTGCCCTCCTTCCCGTGACTGAAGC
 CTCAGAATGCAGTCGGCTAATGTTACCAGATGATACTACAAATCATTCTAACTCCTCCAAGGAGGTCCTT
 TCCTCAGCTGTTTTGAGAAGCCTTCGGGTGAATGTGGTCCAGACGGAGAGGAGACGAGAGCTCAGACTG
 TACAGAAATCCCGGAGTTTTTGTCCACTTCAGAGTCTTCTAGCTTGTGCAAGATCTACAGCCAAGTGA
 TAGCACTTCTTTTATTCTTCTTAACCTAACAAGAGCAGGTCTGGGCTCTTCAGCTGAGCACTTAGTGTTT
 GTACAGGATGAGGCAGAAGATTCAGGGAATGATTTCTCTCCAGTGAGAGCACAGACAGTAGCATTCCAT
 GGTTCTCCGGGTTCCAGAGTTGGCCCATGACAGTTTGATTGCTGCTACTCGTGCAACATGGCAAAGAA
 TGCAAAAACCAGCAGCAATGGAGAAAATGTCCACCTTGGTTCTGGTGATGGGCAGTCAAAAGATTCTGGG
 CCCCTTCTCAAGTGGAAGAAAGCTCAAGTGTACAGTTGAAGGTTGTGACCGGACATTTGTATGGCCAG
 CTCACCTTAAATACCACCTCAAGACTCATCGAAATGACCGCTCCTTCATCTGTCTGCAGAAGGTTGTGG
 GAAAAGCTTCTATGTGCTGCAGAGGCTGAAGGTGCACATGAGGACCCACAATGGAGAGAAGCCCTTTATG
 TGCCATGAGTCTGGCTGTGGTAAGCAGTTTACTACAGCTGGAAACCTGAAGAACCACCGGCGCATCCACA
 CAGGAGAGAAACCTTTCCCTTTGTGAAGCCCAAGGATGTGGCCGTTCCCTTTGCTGAGTATTCTAGCCTCCG
 AAAACATCTGGTGGTTCCTCAGGAGAGAAGCCTCATCAGTGCCAAGTCTGTGGGAAGACCTTCTCTCAG
 AGTGGAAGCAGGAATGTGCATATGAGAAAGCATCACCTGCAGCTGGGAGCAGCTGGGAGTCAAGAGCAGG
 AGCAAACTGAGGTGCTTGCTGAAGGATCCCCACGTTCCCTGTCTTCAGTGCCTGATGTGACACATCACCT
 GGTGACCATGCAGTCAGGAGGCAATCATATGAAGTTTCTGTCTTAAGTGCAGTAAATCCACAAGAGTTA
 CTAACCAAGGAGATTTAACTGAAAGACGGACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC232709 representing NM_001242926
 Red=Cloning site Green=Tags(s)

MLSDELESKPELLVQFVQNTSIPLGQGLVESEAKDITCLSLLPVTEASECSRLMLPDDTTNHSNSSKEVP
 SSAVLRSLRVNVPDGEETRAQTVQKSPEFLSTSESSLLQDLQPSDSTSFILLNLTRAGLGSSAEHLVF
 VQDEAEDSGNDFLSSESTDSSIPWFLRVQELAHDSLIAATRAQLAKNAKTSNNGENVHLGSGDGQSKDSG
 PLPQVEKKLKCTVEGCDRTFVWPAHFYHLKTHRNDRSFCPAEGCGKSFYVLQRLKVHMRTHNGEKPFM
 CHESGCGKQFTTAGNLKNHRIHTGEKPFLEAQGCGRSFAEYSSLRKHLVVHSGEKP HQCQVCGKTF SQ
 SGRNVHMRKHHLQLGAAGSQEQTEVLAEGSPRSLSSVPDVTHHLVTMQSGRQSYEVSVLTAVNPQEL
 LNQGDLTERRT

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001242926

ORF Size: 1293 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001242926.2](#)


RefSeq Size: 2504 bp

RefSeq ORF: 1296 bp

Locus ID: 57862

UniProt ID:	<u>Q86VK4</u>
Cytogenetics:	14q24.3
Protein Families:	Transcription Factors
MW:	47.8 kDa
Gene Summary:	Transcription factor that activates transcription of matrix-remodeling genes such as MMP1 during fibroblast senescence.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC232709