

## Product datasheet for **RC232746**

### **MFSD11 (NM\_001242534) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	MFSD11 (NM_001242534) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MFSD11
Synonyms:	ET
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RC232746 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTC**CCCCGAATCTAAAAAGCTTTTCAACATCATTATTTTAGGAGTTGCCTTTATGTTTATGTTCACTG**  
**CCTTTCAAAC**TTGTGGAAATGTGGCGCAA**ACTGTCATCAGGAGCTTAAATAGGACAGATTTTCACGGCAG**  
**TGGATATACCAGCATGGCTATTATCTATGGAGTGTCTCTGCTTCAAATTTGATTACACCGTCAGTGGTT**  
**GCCATTGTAGGACCTCAACTCTCTATGTTTGCCAGTGGTTATTTTACAGCATGTACATTGCCGTTTTCA**  
**TCCAGCCTTTCCCGTGGTCTTCTACACAGCCTCTGTTTTATTGGAATTGCTGCTGTGCTTTGGAC**  
**AGCACAAGGAACTGCCTGACAATCAATTCGGATGAGCACAGCATTGGGAGAAACAGTGGGATTTCTGG**  
**GCATTCTGCAGTCTAGCTTGTCTTTGGAAATCTCTACATATATTTGCCTGGCAAGGGAAAACTCAGA**  
**TATCAGAGAGTGACCGAAGAACAGTGTATTGCCCTAACGGTGATTAGCCTTGTGGGACAGTTCTATT**  
**CTTTCTCATTGCGAAACAGATTCTGAAAATGTCTAGGAGAAGATGAGTCTTCTGATGACCAGGACATG**  
**GAAGTCAACGAGTCTGCCAGAACAATCTGACAAAGGCAGTAGATGCTTTTAAAAAGTCTTTAAGTTAT**  
**GTGTCACCAAGGAGATGCTCCTTCTTAGTATTACAACCTGTTATACAGGTCTGGAATTAACCTTCTCTC**  
**TGGTGTATATGGAACCTGTATTGGTGTACAAATAAATTTGGAGCAGAAGAGAAAAGCCTTATTGGACTT**  
**TCTGGCATTTCATCGGCATTGGAGAAATTTAGGTGGAAGCCTCTTCGGCCTGCTGAGCAAGAACAATC**  
**GTTTTGGTAGAAATCCAGTTGTGCTGTTGGGCATCCTGGTGCACCTCATAGCTTTTTATCTAATATTTCT**  
**CAACATGCCTGGAGATGCCCGATTGCTCCTGTTAAAGGAACTGACAGCAGTGTACATCAAATCCAGC**  
**AAAGAAGTTGCCATTCTCTGCAGTTTTCTGTTGGCCCTGGAGACAGTCTTTAATACCCAGCTGCTTA**  
**GTATCTTGGGCTTCTGTATTCTGAAGACAGCGCCCAAGCATTGGCCATTTCAAGTTTGTTCAGTCTAT**  
**TTGCGCAGCCGTGGCATTTTTCTACAGCAACTACCTTCTCCTTCACTGGCAACTCCTGGTCATGGTGATA**  
**TTTGGTTTTTTGGAACAATTTCTTTCTCACTGTGGAATGGGAAGCTGCCGCTTTGTAGCCCGCGCT**  
**CTGACTACCGAAGTATC**

**ACGCGT**ACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RC232746 protein sequence  
 Red=Cloning site Green=Tags(s)

MSPE**SKKLFNII**ILGVAFMFMTAFQTCG**NVAQTVIRSLNRTDFHGSGYTSMAIIYGVFSASNLITPSV**  
**AIVGPQLSMFASGLFYSMYIAVFIQPPW**SFYTASVFIGIAAAVLWTAQGNCLTINSDEHSIGRNSGIFW  
 ALLQSS**LLFFGNLYIYFAWQK**TQISESDRRTVFI**ALTVISLVGT**VLFFLIRKPDSE**NVLGEDESSDDQDM**  
 EVNES**AQNNLTKAVDAFKKSF**KLCVTKEM**LLSITTAYTGLEL**TFPSGVYGT**CIGATNKFGAEEKSLIGL**  
 SGIFIGI**GEILGGSLFLLSKNNR**FGRNPV**LLGILVHFIAFYLI**FLNMPGD**APIAPVKGT**DSSAYIKSS  
 KEVA**ILCSFLLGLGDS**FN**TQLLSILGFLY**SEDSAPAF**AIFKFVQ**SICAAVAFFY**SNYLLLHWQLL**VMVI  
 FGFFGTIS**FFTVEW**EAAAFVARGSDYRSI

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Chromatograms:**

[https://cdn.origene.com/chromatograms/mk8112\\_b10.zip](https://cdn.origene.com/chromatograms/mk8112_b10.zip)

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_001242534

**ORF Size:** 1347 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001242534.2](#)

**RefSeq Size:** 2850 bp

**RefSeq ORF:** 1350 bp

**Locus ID:** 79157


**UniProt ID:** [O43934](#)

**Cytogenetics:** 17q25.1-q25.2

Protein Families: Transmembrane

MW: 49.2 kDa

Product images:



Circular map for RC232746