

Product datasheet for **RC232795**

AKT3 (NM_001206729) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AKT3 (NM_001206729) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	AKT3
Synonyms:	MPPH; MPPH2; PKB-GAMMA; PKBG; PRKBG; RAC-gamma; RAC-PK-gamma; STK-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC232795 representing NM_001206729
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAGCGATGTTACCATTGTGAAAGAAGGTTGGGTTTCAGAAGAGGGGAGAATATATAAAAACTGGAGGC
CAAGATACTTCTTTTGAAGACAGATGGCTCATTCATAGGATATAAAGAGAAACCTCAAGATGTGGATTT
ACCTTATCCCCTCAACAACCTTTTCAGTGGCAAAATGCCAGTAAATGAAAAAGAACGACCAAGCCAAAC
ACATTTATAATCAGATGTCTCCAGTGGACTACTGTTATAGAGAGAACATTTTCATGTAGATACTCCAGAGG
AAAGGGAAGAATGGACAGAAGCTATCCAGGCTGTAGCAGACAGACTGCAGAGGCAAGAAGAGGAGAGAAT
GAATTGTAGTCCAACCTCACAAATTGATAATATAGGAGAGGAAGAGATGGATGCCTCTACAACCCATCAT
AAAAGAAAGACAATGAATGATTTTGACTATTTGAAACTACTAGGTAAGGCATTTTGGGAAAGTTATTT
TGGTTCGAGAGAAGGCAAGTGGAAAATACTATGCTATGAAGATTCTGAAGAAAGAAGTCATTATTGCAA
GGATGAAGTGGCACACACTCTAACTGAAAGCAGAGTATTAAGAACAAGTACATCCCTTTTAAACATCC
TTGAAATATTCTTCCAGACAAAAGACCGTTTGTGTTTTGTGATGGAATATGTTAATGGGGCGAGCTGT
TTTTCCATTTGTGCGAGAGAGCGGGTCTCTGAGGACCGCACACGTTTCTATGGTGCAGAAATTGTCTC
TGCTTGGACTATCTACATTCGGAAAGATTGTGTACCGTGTACTCAAGTTGGAGAATCTAATGCGGGAC
AAAGATGGCCACATAAAAATTACAGATTTTGGACTTTGCAAGAAGGGATCACAGATGCAGCCACCATGA
AGACATTCTGTGGCACTCCAGAATATCTGGCACCAGAGGTGTTAGAAGATAATGACTATGGCCGAGCAGT
AGACTGGTGGGCGCTAGGGTGTGTCATGTATGAAATGATGTGTGGGAGGTTACCTTTCTACAACCAGGAC
CATGAGAACTTTTGAATTAATATTAATGGAAGACATTAATTTCTCGAACACTCTCTCAGATGCAA
AATCATTGCTTTTCAGGGCTCTTGATAAAGGATCCAAATAAACGCCTTGGTGGAGGACAGATGATGCAAA
AGAAATTATGAGACACAGTTTCTCTCTGAGTAAACTGGCAAGATGTATATGATAAAAAGCTTGTACCT
CCTTTTAAACCTCAAGTAAACATCTGAGACAGATACTAGATATTTTGTGAAAGATTTACAGCTCAGACTA
TTACAATAACACCACCTGAAAAATGTCAGCAATCAGATTGTGGCATGCTGGGTAAGTGGAAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC232795 representing NM_001206729
Red=Cloning site Green=Tags(s)

MSDVTIVKEGWVQKRGEYIKNWRPRYFLLKTDGSFIGYKEKPQDVL PYPLNNSVAKCQLMKTERPKPN
TFIIRCLQWTTVIERTFHVDTPEREWEATEIQAVADRLQRQEERMNCSPTSQIDNIGEEEMDASTTHH
KRKTMNDFDYLKLLGKGTFGKVLVREKASGKYAMKILKKEVIAKDEVAHTL TESRVLKNTRHPFLTS
LKYSFQTKDRLCFVMEYVNGGELFFHL SRERVF SEDRTRFYGAEIVSALDYLHSGKIVYRDLKLENLMRD
KDGHIKITDFGLCKEGITDAATMKTF CGTPEYLAPEVLEDNDYGRAVDWWGLGVVMEYEMMGRLPFYNQD
HEKLFELILMEDIKFPRTLSSDAKSLLSGLLIKDPNKRLGGPDDAKEIMRHSFFSGVNWQDVYDKKLV
PFKPQVSETDTRYFDEEFTAQTITITPPEKQQSDCGMLGNWKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6563_b02.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_001206729

ORF Size: 1395 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001206729.1](#), [NP_001193658.1](#)

RefSeq Size: 1584 bp

RefSeq ORF: 1398 bp

Locus ID: 10000

UniProt ID: [Q9Y243](#)

Cytogenetics: 1q43-q44


Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase

Protein Pathways: Acute myeloid leukemia, Adipocytokine signaling pathway, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Glioma, Insulin signaling pathway, Jak-STAT signaling pathway, MAPK signaling pathway, Melanoma, mTOR signaling pathway, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer, Renal cell carcinoma, Small cell lung cancer, T cell receptor signaling pathway, Tight junction, Toll-like receptor signaling pathway, VEGF signaling pathway

MW: 54.1 kDa

Gene Summary: The protein encoded by this gene is a member of the AKT, also called PKB, serine/threonine protein kinase family. AKT kinases are known to be regulators of cell signaling in response to insulin and growth factors. They are involved in a wide variety of biological processes including cell proliferation, differentiation, apoptosis, tumorigenesis, as well as glycogen synthesis and glucose uptake. This kinase has been shown to be stimulated by platelet-derived growth factor (PDGF), insulin, and insulin-like growth factor 1 (IGF1). Alternatively splice transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC232795