

Product datasheet for **RC232992**

CEP76 (NM_001271989) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CEP76 (NM_001271989) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CEP76
Synonyms:	C18orf9; HsT1705
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC232992 representing NM_001271989
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCGCTGCCTCCGGAGAAAGCCTCCGAGCTGAAGCAGCTCATCCACCAGCAGCTGAGCAAGATGGATG
 TCCATGGTAGAATAAGAGAAATCCTTGCTGAGACTATACGGGAAGAAATGGCACCTGATCAACAGCATTT
 ATCAACAGAAAGATTTGATCAAAGCCCTTAGACGTCGAGGAATCATTGACGATGTGATGAAAGAATTAAT
 TTTGTTACTGACAGTGTGAGCAAGAACTCCCTTCTCTCCAAAAACAACCTATTTGTTTTGATAGACAAT
 CGACATTAATAAAAAAGTATGGAAGTGAATGGCTGATTAACAACAATGTTATCAATAAGTATCAAT
 TCATATGGTGCTAATCAAACAGACATATTTGGTGAGACGACTTAGTAGCATCATATTTCTGGAATGG
 CGATCGGTTTTGGGCTCAGAAAATGGAGTGACCAGTCTGACTGTGGAATTATGGGTGTAGGCACAGAAT
 CAAAAGTTTCTGTGGGAATTTAAATATAAACTTGAATGTATCCACCCTCAATCAAACGTTATCTCA
 AGAAGTAGTGAACACACAGCTTGCTTTGGAACGTCAGAAAAGAGCGGATTATTTCTTGTA
 TATGCTAAGCAGTGGTGGAGAGAATATTTGCAAATTCGACCCTCACACAACCTCACGACTGGTTAAGATT
 TTGCACAGGATGAAAATGGGATAAATAGACCAGTCTGTTCCATGTTAAACCCTTCGAGCTGGACGGCT
 TCTTGATACTCCAAGGCAAGCAGCAAGATTTGTTAATGTCTTGGTTATGAACGAGCCCTGTTATTGGA
 GGAGGAGGTAACAGGAGCAGTGGTGCCTCTGCTGGCCTTTCTCTGTAGAAAACAAGGGTACTGTGAAG
 ATCACGCTAACCTTCTGTGCAGCCTTCTTCTGGATATGGATTAGAAGCCTTTGTTTGTGTTGGGACCA
 GGCAAAAGGAGTACCTCATGCATGGGTTATGACTGTGGAAGTATGAGGAGGATCACTTTTTGGGAGAGT
 TTAACAGGACACAGGTACATCCATAAACCTACCAATCCTGATGAACCTCCAGTTGCTGAACAGCCAAAC
 CACTGTACCCATATCGAACAATGGTTGTTTTCAACCATCAGATGTTCTGGGAAATGTCAACCTC
 TGATGCAGTAGAAACCTGTGATTTGATTTGAACGATGAATCCAAATGGAAACCATGAGTGAGGAAGCA
 ATTAATCTGTGTGTGCTCCTGGAGCTACAACATCCCTTCTCCCTTTCCACCTCTGTGTGCATCCACAA
 TTGACGCGTCAGTAACAAGTAATGAAATGAAATGCAGCTGAGGCTCCTGGTGTGAGAACACAGGAAGGA
 TCTTGGCCTCACTACTGTTTGGGAAGACCAGCTCTCTACCTTTTATCACAGCTTTGGCTTCTTATGAA
 TTTGAGCGTACAACAAGTATATCAGCAGGCAATGAAGAATTTCAAGATGCCATAAGAAGGGCTGTACCTG
 ATGGTCACACATTTAAGGGTCCCAATACATTTTGTGTATAGAAATGCAAGACGTGCATTTGCCACATG
 TCTTCGATCTCCTTTCTGTGAAGAAATAATCTGTTGCCGTGGAGACCAAGTGCAGCTGGCAGTTCGTGTC
 CGAGTATTACTTACCCTGAATCTGCATGTGCTGTTGGATCATGTTGCTTGTAAATATCGCTCGGTAT
 TA

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC232992 representing NM_001271989
 Red=Cloning site Green=Tags(s)


MSLPPEKASELKQLIHQQLSKMDVHGRIREILAETIREELAPDQQHLSTEDLIKALRRRGIIDVVKELN
 FVTDSVEQELPSSPKQPICFDRQSTLKKSDGTRMADSTTMLSIDPIHMLIKTDIFGETTLVASVFLEW
 RSVLGSSENGVTSLTVELMGVGTESKVSIGILNIKLEMYPLNQTLSQEVVNTQLALERQKTAEKERLFLV
 YAKQWREYLQIRPSHNSRLVKIFAQDENGINRPVCSYVKPLRAGRLDTPRQAARFVNVLYGERAPVIG
 GGGKQEQWCTLLAFLCRNKGDCEHDHANLLCSLLLGYGLEAFVVCVGTAKAGVPHAWVMTCGTDGAIWFES
 LTGHRYIHKPTNPDEPPVAEQPKPLYPYRTIGCVFNHQMFLGNCQPSDAVETCVFDLNDESKWKPMSEEA
 IKSVCAPGATTSPPFPPLCASTIDASVTSNEIEMQLRLLVSEHRKDLGLTTVWEDQLSYLLSPALASYE
 FERTTISIAGNEEFQDAIRRAVPDGHFTKGFPIHFVYRNARRAFATCLRSPFCEEIICCRGDQVRLAVRV
 RVFTYPESACAVWIMFACKYRSVL

SGPTRTRPLEQKLI**SEEDLAANDILDYKDDDDKV**


Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001271989

ORF Size: 1752 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001271989.1 , NP_001258918.1
RefSeq Size:	2722 bp
RefSeq ORF:	1755 bp
Locus ID:	79959
UniProt ID:	Q8TAP6
Cytogenetics:	18p11.21
MW:	66.2 kDa
Gene Summary:	This gene encodes a centrosomal protein which regulates centriole amplification by limiting centriole duplication to once per cell cycle. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012]