

Product datasheet for **RC233019**

MELK (NM_001256689) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MELK (NM_001256689) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MELK
Synonyms:	HPK38
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC233019 representing NM_001256689
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGATGAACTTCTCAAATATTATGAATTACATGAAACTATTGGGACAGAGTGATTGGCCCGGATCAAAA
 CGGAGATTGAGGCCTTGAAGAACCTGAGACATCAGCATATATGTCAACTCTACCATGTGCTAGAGACAGC
 CAACAAAATATTCATGGTTCTTGAGTACTGCCCTGGAGGAGAGCTGTTTGACTATATAATTTCCAGGAT
 CGCCTGTCAGAAGAGGAGACCCGGTGTCTCCGTCAGATAGTATCTGCTGTTGCTTATGTGCACAGCC
 AGGGCTATGCTCACAGGGACCTCAAGCCAGAAAATTTGCTGTTTGTGAATATCATAAATTAAGCTGAT
 TGACTTTGGTCTCTGTGCAAAACCAAGGTAACAAGGATTACCATCTACAGACATGCTGTGGGAGTCTG
 GCTTATGCAGCACCTGAGTTAATACAAGGCAAATCATATCTTGGATCAGAGGCAGATGTTGGAGCATGG
 GCATACTGTTATATGTTCTTATGTGGATTTCTACCATTGATGATGATAATGTAATGGCTTTACAA
 GAAGATTATGAGAGGAAAATATGATGTTCCAAGTGGCTCTCCAGTAGCATTCTGCTCTTCAACAA
 ATGCTGCAGGTGGACCCAAAGAAACGGATTTCTATGAAAAATCTATTGAACCATCCCTGGATCATGCAAG
 ATTACAATATCCTGTTGAGTGGCAAAGCAAGAATCCTTTTATTACCTCGATGATGATTGCGTAACAGA
 ACTTTCTGTACATCACAGAAACAACAGGCAACAATGGAGGATTAATTTCACTGTGGCAGTATGATCAC
 CTCACGGCTACCTATCTTCTGCTTCTAGCCAAGAAGGCTCGGGGAAAACAGTTCGTTTAAAGCTTTCTT
 CTTTCTCTGTGGACAAGCCAGTGTACCCATTCACAGACATCAAGTCAAATAATTGGAGTCTGGAAGA
 TGTGACCGCAAGTGATAAAAAATTATGTGGCGGGATTAATAGACTATGATTGGTGTGAAGATGATTTATCA
 ACAGGTGCTGCTACTCCCGAACATCACAGTTTACCAAGTACTGGACAGAATCAAATGGGTGGAATCTA
 AATCATTAACTCCAGCCTTATGCAGAACACCTGCAAATAAATTAAGAACAAGAAAATGTATATACTCC
 TAAGTCTGCTGTAAGAATGAAGAGTACTTTATGTTTCCCTGAGCCAAAGACTCCAGTTAATAAGAACCAG
 CATAAGAGAGAAAATACTCACTACGCCAAATCGTTACACTACCCCTCAAAGCTAGAAAACAGTGCCTGA
 AAGAACTCCAATTAATAACAGTAAATCAACAGGAACAGACAAGTTAATGACAGGTGTCATTAGCCC
 TGAGAGGCGGTGCCGCTCAGTGGAAATGGATCTCAACCAAGCACATATGGAGGAGACTCCAAAAAGAAAG
 GGAGCCAAAGTGTGGGAGCCTTGAAGGGGGTTGGATAAGGTTACTCTGTGCTACCAGGAGCAAAA
 GGAAGGGTTCTGCCAGAGACGGGCCAGAAGACTAAAGCTTCACTATAACGTGACTACAAGTATAGT
 GAATCCAGATCACTGTTGAATGAAATAATGTCTATTCTTCAAAGAAGCATGTTGACTTTGTACAAAAG
 GGTATACACTGAAGTGTCAAACACAGTCAGATTTGGGAAAGTGACAATGCAATTTGAATTAGAAGTGT
 GCCAGCTTCAAAAACCCGATGTGGTGGTATCAGGAGGCAGCGGTTAAGGGCGATGCCTGGGTTTACAA
 AAGATTAGTGAAGACATCCTATCTAGCTGCAAGGTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTAA

Protein Sequence:


>RC233019 representing NM_001256689
 Red=Cloning site Green=Tags(s)

MMNFSNIMNYMKLLGQSDLPRIKTEIEALKNLRHQHICQLYHVLETANKIFMVLEYCPGGELFDYIISQD
 RLSEETRVRVFRQIVSAVAYVHSQGYAHRDLKPENLLFDEYHKLKIDFGLCAKPKGNKDYHLQTCGSL
 AYAAPELIQGKSYLGSEADVWSMGILLVLMCGFLPFDDDNVMALYKIMRGKYDVPKWLSPSSILLQQ
 MLQVDPKKRISMKNLLNHPWIMQDYNYPVEWQSKNPF IHLDDCVTEL SVHHRNNRQTMEDLISLWQYDH
 LTATYLLLLAKKARGKPVRLRLSSFSCGQASATPFDIKSNWSLEDVTSKDKNYVAGLIDYDWCEDDLS
 TGAATPRTSQFTKYWTESNGVESKSLTPALCRTPANKLKNKENVYTPKSAVKNEEYFMFPEPKTPVNKNQ
 HKREILTPNRYTTPSKARNQCLKETPIKIPVNSTGTDKLMTGVI SPERRCRSVELDLNQAHMEETPKRK
 GAKVFGSLERGLDKVITVLRTRSKRKGSAARDGPRRLKLHYNVTTTRLVNPQQLNEIMSILPKKHVDFVQK
 GYTLKQCQTQSDFGKVTMQFELEVQQLQKPDVVGIRRLKGDWVYKRLVEDILSSCKV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001256689

ORF Size: 1857 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256689.1 , NP_001243618.1
RefSeq Size:	2400 bp
RefSeq ORF:	1860 bp
Locus ID:	9833
UniProt ID:	Q14680
Cytogenetics:	9p13.2
Protein Families:	Druggable Genome, Protein Kinase
MW:	71.6 kDa
Gene Summary:	Serine/threonine-protein kinase involved in various processes such as cell cycle regulation, self-renewal of stem cells, apoptosis and splicing regulation. Has a broad substrate specificity; phosphorylates BCL2L14, CDC25B, MAP3K5/ASK1 and ZNF622. Acts as an activator of apoptosis by phosphorylating and activating MAP3K5/ASK1. Acts as a regulator of cell cycle, notably by mediating phosphorylation of CDC25B, promoting localization of CDC25B to the centrosome and the spindle poles during mitosis. Plays a key role in cell proliferation and carcinogenesis. Required for proliferation of embryonic and postnatal multipotent neural progenitors. Phosphorylates and inhibits BCL2L14, possibly leading to affect mammary carcinogenesis by mediating inhibition of the pro-apoptotic function of BCL2L14. Also involved in the inhibition of spliceosome assembly during mitosis by phosphorylating ZNF622, thereby contributing to its redirection to the nucleus. May also play a role in primitive hematopoiesis.[UniProtKB/Swiss-Prot Function]