

Product datasheet for **RC233042**

NOL10 (NM_001261394) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NOL10 (NM_001261394) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NOL10
Synonyms:	PQBP5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC233042 representing NM_001261394
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCAGGTCTCCAGCCTCAATGAGGTGAAGATTTACAGCCTCAGCTGCGGCAAGTCCCTTCTGAGTGGC
 TTTCTGATAGGAAGAAGAGAGCGCTACAGAAGAAAGATGTAGATGTCCGTAGGAGAATTGAACTTATTCA
 GGACTTTGAAATGCCTACTGTGTGTACCACTATTAAGGTGTCAAAAGATGGACAGTACATTTTAGCAACT
 GGAACATATAAACCTCGGGTTCGATGTTATGACACCTATCAATTATCCTTGAAGTTTGAAAGGTGTTTAG
 ATTCAGAAGTTGTCACCTTTGAAATTTGTCTGATGACTACTCAAAGATTGTCTTCTACATAATGATAG
 ATACATTGAATTTCAATCGCAATCAGGTTTTTACTACAAAACCAGAATACCAAAGTTTGGGAGAGATTTT
 TCTTACCACTATCCATCCTGTGACTTGTACTTTGTTGGTCAAGTTCTGAAGTTTATAGGTTAACTTAG
 AACAGGACGATACCTGAATCCTCTACAACTGATGCTGCGGAGAATAATGTTTGTGACATAAATTCAGT
 GCATGGCTTGTGGCACAGGAACCATAGAGGGTAGAGTGGAGTCTGGGACCCAAGAACTCGAAACAGA
 GTTGGCCTGTAGACTGCGCCTTAAACAGTGTACAGCAGATTACAGATAAACAGTTTACCAACAACTCT
 CTGCTTTGAAATTTAATGGTGCCTTGACCATGGCAGTTGGAACAACACAGGGCAGGGAAAAATATTAC
 TTCCTTGAGCCAGAGCATGACCTTAATGATGTTTGTCTCTACCCAACTCAGGCATGCTTCTGACGGCC
 AATGAAACCCCAAGATGGGCATCTATTACATTCCAGTTTTGGGTCTGCTCCTCGGTGGTGTCTCTTCT
 TAGACAACTTGACCGAAGAATTAGAAGAGAATCCAGAAAGCACAGTCTATGATGATTATAAATTTGTCAC
 CAAGAAAGACCTTGAAAAATTTAGGGCTCACCCACCTCATTGGATCTCCTTTCTCCGGGCATATATGCAT
 GGGTTTTTCATGGATATAAGACTCTATCACAAGGTGAAACTGATGGTAAATCCATTTGCTTATGAAGAAT
 ATAGGAAAGATAAAATACGACAGAAAAATAGAAGAAACACGTGCACAGAGAGTCCAGTTAAAGAAATTGCC
 AAAAGTTAACAAAGAGCTGGCACTTAAATTAATTGAGGAAGAAGAGGAGAAGCAGAAATCTACATGGAAA
 AAGAAAGTTAAGAGTCTTCTAATATTCTCACCGATGATCGATTTAAAGTTATGTTTGAGAACCTGACT
 TCCAAGTAGATGAAGAGAGTGAAGAATTTAGGCTTCTGAATCCACTTGTTTCAAAAATTAGTAAAAAAG
 GAAGAAGAACTAAGACTCTTAGAGCAACAAGAACTTCGTGAAAAAGAAGAGGAGGAAGAGCCGGAAGGA
 AAACCAAGTGATGCAGAAAGTTCGGAGAGTTCAGATGATGAAAAAGCCTGGGTGAAGAGGTGAGGAAGC
 AACGCAGACTCCTCCAGCAGGAGGAAAAAGTGAAGCGGCAGGAACGACTCAAGGAGGACCAGCAGACAGT
 CCTAAAGCCCCAGTTTTATGAGATCAAAGCAGGAGAAGAATTTAGAAGCTTCAAAGATTCTGCCACAAAG
 CAGAACTGATGAACAAAACCTTGAAGATCGTTTGAAAATTGAAGCAAAAAATGGGACATTGAGTGTAT
 CCGACACCACCGTTGGCAGCAACAATTGACATTACGTTAAAGAGGTCTGAACAGCAGAAGAAGCAACA
 GGAGGCTGAGAACTGCATCGACAAGAAAGGAAAGACTCCGTCGTTTCGGCCGGACACCTGAAGTCAAGA
 CACAAAAGAGGACGGTCGTTTCAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC233042 representing NM_001261394
 Red=Cloning site Green=Tags(s)

MQVSSLNEVKIYSLSCGKSLPEWLSDRKKRALQKKDQVRRRIELIQDFEMPTVCTTIKVSQKDGQYILAT
 GTYKPRVRCYDITYQLSLKFERCLDSEVVTFEILSDDYSKIVFLHNDRIEFHSQSGFYKTRIPKFGDRF
 SYHYPSCDLYFVGASSEVYRLNLEQGRYLNPLQTDAAENNVCDINSVHGLFATGTIEGRVECWDPRTRNR
 VGLLDCALNSVTADSEINSLPTISALKFNGALTMVGTGQKIFTSLEPEHDLNDVCLYPNSGMLLTA
 NETPKMGIYYIPVLGPAPRWCSFLDNLTEELEENPESTVYDDYKFVTKKDLENLGLTHLIGSPFLRAYMH
 GFFMDIRLYHKVKLMVNPFAEYEEYRKDKIRQKIEETRAQVRVQLKKLPKVNKELALKLIEEEEEKQKSTWK
 KKVKSLPNILTDDRKFVMFENPDFQVDEESEEFRLNPLVSKISEKRKKKLRLLEQQELREKEEEEEPEP
 KPSDAESSESDDEKAWVEVRQRLLQEEKVKRQERLKEDQQTVLKPQFYEIKAGEEFRSFKDSATK
 QKLMNKTLEDRLIEAKNGTLSVSDTTVGSKQLTFTLKRSEQKKQAEKLRQERKLRRSAGHLKSR
 HKRGRSFH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Cloning Scheme:



ACCN: NM_001261394

ORF Size: 1914 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001261394.1, NP_001248323.1

RefSeq Size: 3365 bp

RefSeq ORF: 1917 bp


Locus ID: 79954

UniProt ID: Q9BSC4

Cytogenetics: 2p25.1

MW: 75 kDa

Product images:



Circular map for RC233042