

## Product datasheet for **RC233065**

### ATRIP (NM\_001271022) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ATRIP (NM_001271022) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ATRIP
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RC233065 representing NM\_001271022  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGAAAGTAAATGGAAGAAGAAGTTCTCATTAAAGATGGAGAAATTTAAATTTTTCGAGACTCACTACATC  
AGACGGAATCCGTTCTAGAGGAACAGAGAAGATCACATTTTCTTCTTGAGCAAGAGAAAACCAAGCACT  
CAGTGACAAGGAAAAGGAATTCTCCAAAAGCTCCAATCATTGCACTGAACTCCAGTTTAAAGATGCA  
GAGATGAATGAATTAAGGACAAAGCTCCAGACCAGTGAACGAGCAATAAACTGGCTGCTCCCTCTGTTT  
CCCATGTCAGTCTAGGAAAAACCTTCTGTGGTTATAAAGCCAGAAGCATGTTCTCCACAATTTGGAAA  
AACATCTTTTCTACAAAGGAGTCTTTTAGTGCTAACATGTCCCTTCCCCACCCCTGCCAGACGGAGTCA  
GGATAAAGCCTCTGGTGGGCAGAGAGGATAGTAAGCCCCACAGTCTGAGAGGTGACTCCATAAAACAAG  
AAGAGGCCAGAAAAGCTTTGTTGACAGCTGGAGACAGAGATCAAACACTCAAGGTTCCATTTTGATAAA  
CCTGCTCTGAAGCAGCCTTTGATCCCAGGGTCATCCCTAAGCCTTGGCCACCTCCTGAGTAGTGTCT  
GAGTCTCCTGCTGGCACCCTTGCAGCCACCAGGGTTTGGCAGTACCTTGGCTGGAATGTCAGGCCTCA  
GGACCACAGGTTCTTATGATGGGTCATTTCCCTCTCAGCCCTGAGAGAAGCACAGAACCTGGCATTAC  
TGGACTGAATCTGGTTGCCCGGAATGAGTGCTCACGTGATGGAGACCCAGCAGAGGGAGGCAGAAGGGCC  
TTCCCACTCTGCCAGTCTCTGGAGCCGTGATTTCTCCCTTGTACAGTTCTTCATCGGCTTACACT  
GCCAGGCCTGCAGGACTTGGCAGCTGCTAAGAGAAGCGGAGCACCTGGGACTCACCGACATTCCTC  
CTGCGTGAGCTCTGGGTAGAGACCAACCTGAGGACTCAGTGTGCATCCTGGAAGGCTTCTCTGTGACT  
GCATTAGCATTCTCAGCACCTGGTGTCCACAGCGGAGCAGTCGTCTCCCTATTACTGTCAGGAGTGG  
GGCAGATTTCTGCTGCTGGGAAGGAAACAGGAGCCTGGTTCACAGGCTTAGTGATGGAGATATGACCTC  
AGCCCTAAGGGGGTTGCTGATGACCAAGGACAGCACCCACTGTTGAAGATGCTTCTCACCTGTTGGCT  
TTCTTCTGTCAGCAACAGGTCACCTTCAAGCCAGTGTCTGACCCAGTGCCTAAGGTTTTGGTAAAT  
TAGCCGAAAACACTTCTGTGATTTCTTGCCAGGTTCCAGTGTGTGTTCCAAGTCTGCCAAAGTGCCT  
CAGCCCAGAGACACCCCTGCCTAGCGTGTCTGGCTGTTGAGCTCCTCTCCCTGCTGGCGGACCACGAC  
CAGCTGGCACCTCAGCTCTGTTCCCACTCAGAAGGCTGCCTCCTGCTGCTGTACATGTACATCACAT  
CACGGCCTGACAGAGTGGCCTTGGAGACACAATGGCTCCAGCTGGAACAAGAGGTGGTGTGGCTCCTGGC  
TAAGCTTGGTGTGCAGAGCCCTTGCCCCAGTCACTGGCTCCAAGTCCAGTGTAAATGTGGAGGTGGTC  
AGAGCGCTCACGGTGTGTTGCACAGACAGTGGCTGACAGTGGGAGGGCAGGGGACCCCAAGGACCG  
ACCAGCAGAGGGCGACAGTGCCTGTCTGCGGGACAGGTGCTGCTGCTGCACGGCTATCGCAGAAGGA  
CAAGCTTTCATGATGCACTGCGTGGAGGTCCTGCATCAGTTTGACCAGGTGATGCCGGGGTGCAGCATG  
CTCATCCGAGGGCTTCTGATGTGACGGACTGTGAAGAGGCAGCCCTGGATGACCTCTGTGCCGCGGAAA  
CCGATGTGAAGACCCCGAGGTGGAGTGTGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC233065 representing NM\_001271022  
 Red=Cloning site Green=Tags(s)


MKVMEEEVLIKNGEIKILRDSLHQTESVLEEQRSSHFLLEQEKQALSDKEKEFSKKLQSLQSELQFKDA  
 EMNELRTLKQTSERANKLAAPSVSHVSPRKNPSVVIKPEACSPQFGKTSFPTKESFSANMSLPHPCQTES  
 GYKPLVGREDSKPHSLRGDSIKQEAAQKSFVDSWRQRSNTQGSILINLLKQPLIPGSSLSLCHLLSSSS  
 ESPAGTPLQPPGFGSTLAGMSGLRTTGSYDGSFSLSALREAQNLAFTGLNLVARNECSRDPGPAEGGRRRA  
 FPLCQLPGAVHFLPLVQFFIGLHCQALQDLAAAKRSGAPGDSPTHSSCVSSGVETNPEDSVCILEGFSVT  
 ALSILQHLVCHSGAVVSLLLSGVGADSAAGEGNRSLVHRLSDGDMTSALRGVADDQGHPLLKMLLHLLA  
 FSSAATGHLQASVLTQCLKVLVLAENTSCDFLPRFQCVFQVLPKCLSPETPLPSVLLAVELLSSLADHD  
 QLAPQLCSHSEGCLLLLLYMYITSRPDRVALETQWLQLEQEVVWLLAKLGVQSPLPPVTGSNCQCNEVV  
 RALTVMLHRQWLTVRRAGGPPRTDQQRRTVRCLRDTVLLLHGLSQKDKLFMMHCVEVLHQFDQVMPGVSM  
 LIRGLPDVTDCEEALDDLCAEATDVEDPEVECG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



## Plasmid Map:



ACCN: NM\_001271022

ORF Size: 1992 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001271022.1](#), [NP\\_001257951.1](#)

RefSeq Size: 2419 bp

RefSeq ORF: 1995 bp

Locus ID: 84126

UniProt ID: [Q8WXE1](#)

Cytogenetics: 3p21.31

MW: 72.8 kDa

**Gene Summary:** This gene encodes an essential component of the DNA damage checkpoint. The encoded protein binds to single-stranded DNA coated with replication protein A. The protein also interacts with the ataxia telangiectasia and Rad3 related protein kinase, resulting in its accumulation at intranuclear foci induced by DNA damage. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2012]