

Product datasheet for **RC233073**

IL31RA (NM_001242637) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IL31RA (NM_001242637) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	IL31RA
Synonyms:	CRL; CRL3; GLM-R; GLMR; GPL; hGLM-R; IL-31RA; PLCA2; PRO21384; zcytoR17
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC233073 representing NM_001242637
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGCATCAGGCAACTCAAGTTTTCCACACGGCATGTGTCTGTGAATGTCGCAAACATTCTCTCTC
CCCAGCCTTCATGTGTTAACCTGGGGATGATGTGGACCTGGGCACTGTGGATGCTCCCTCACTCTGCAA
ATTCAGCCTGGCAGCTCTGCCAGCTAAGCCTGAGAACATTTCTGTGTCTACTACTATAGAAAAATTTA
ACCTGCACTTGGAGTCCAGGAAAGGAAACCAGTTATACCCAGTACACAGTTAAGAGAACTTACGCTTTTG
GAGAAAAACATGATAATTGTACAACCAATGTTCTACAAGTAAAATCGTGCTTCGTGCTCTTTTTCT
TCCAAGAATAACGATCCCAGATAATTATACCATTGAGGTGGAAGCTGAAAATGGAGATGGTGAATTA
TCTCATATGACATACTGGAGATTAGAGAACATAGCGAAAACCTGAACACCTAAGATTTCCGTGTGAAAC
CAGTTTTGGGCATCAAACGAATGATTCAAATTGAATGGATAAAGCCTGAGTTGGCGCCTGTTTCATCTGA
TTTAAAAATACACACTTCGATTCAGGACAGTCAACAGTACCAGCTGGATGGAAGTCAACTTCGCTAAGAAC
CGTAAGGATAAAAACCAACGTACAACCTCACGGGGCTGCAGCCTTTTACAGAATATGTCATAGCTCTGC
GATGTGCGGTCAAGGAGTCAAAGTCTGGAGTGACTGGAGCCAAGAAAAAATGGGAATGACTGAGGAAGA
AGCTCCATGTGGCCTGGAAGTGTGGAGAGTCTGAAACCAAGCTGAGGCGGATGGAAGAAGGCCAGTGCGG
TTGTTATGGAAGAAGGCAAGAGGAGCCCAAGTCTAGAGAAAACACTTGGCTACAACATATGGTACTATC
CAGAAAGCAACTAACCTCACAGAAAACATGAACACTACTAACCAGCAGCTTGAAGTGCATCTGGGAGG
CGAGAGCTTTGGGTGTCTATGATTTCTTATAATTCTCTGGGAAGTCTCCAGTGGCCACCCTGAGGATT
CCAGCTATTCAAGAAAAATCATTTCAGTGCATTGAGGTGATGCAGGCCTGCGTTGCTGAGGACCAGCTAG
TGGTGAAGTGGCAAAGCTCTGCTTAGACGTGAACACTTGGATGATTGAATGGTTCCGGATGTGGACTC
AGAGCCCACCACCTTTCTGGGAATCTGTGTCTCAGGCCACGAACCTGGACGATCCAGCAAGATAAATTA
AAACCTTTCTGGTGTATAACATCTCTGTGTATCCAATGTTGCATGACAAAAGTTGGCGAGCCATATTCCA
TCCAGGCTTATGCCAAAAGAAGGCGTTCCATCAGAAGTCTGAGACCAAGGTGGAGAACATTGGCGTGAA
GACGGTACGATCACATGGAAAGAGATTCCAAGAGTGAAGAAAAGGATATCATCTGCAACTACACCATC
TTTTACCAAGCTGAAGGTGGAAAAGGATTCTCCAAGACAGTCAATTCCAGCATCTTGCACTACGGCCTGG
AGTCCCTGAAACGAAAGACCTCTTACATTGTTCCAGTCTGAGGATATCCTCATAACTTCTCTGATTGGT
GGAGGCCCTTCTTATTCTCATTATCCTGACAGTGGCATATGGTCTCAAAAAACCAACAAATTGACTCATC
TGTGTTGGCCACCGTTCCCAACCCTGCTGAAAGTAGTATAGCCACATGGCATGGAGATGATTTCAAGGA
TAAGCTAAACCTGAAGGAGTCTGATGACTCTGTGAACACAGAAGACAGGATCTTAAACCATGTTCCACC
CCCAGTGACAAGTTGGTATTGACAAGTTGGTGGTGAACCTTTGGGAATGTTCTGCAAGAAATTTTACAG
ATGAAGCCAGAACGGTTCAGGAAAACAATTTAGGAGGGGAAAAGAATGGGACTAGAATTCTGTCTTCTG
CCCAACTCAATA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC233073 representing NM_001242637
 Red=Cloning site Green=Tags(s)


MCIRQLKFFTTACVCECPQNILSPQPSCVNLGMMWTWALWMLPSLCKFSLAALPAKPENISCVYYYRKNL
 TCTWSPGKETSQYQYTVKRTYAFGEKHDNCTTNSSTSENRASCSFFLPRITIPDNYTIEVEAENGDGVK
 SHMTYWRLENIAKTEPPKIFRVKPVLGIKRMIQIEWIKPELAPVSSDLKYTLRFRTVNSTSWMEVNF
 AKNRKDKNQTYNLTGLQPFTEYVIALRCAVKESKFWSDWSQEKMGMTTEEEAPCGLELWRVLPKPAEADGRRPVR
 LLWKKARGAPVLEKTLGYNIWYYPESNTNLTETMNTNQQLELHLGGESFWVSMISYNSLGKSPVATLRI
 PAIQEKSFQCIQEVMAQACVAEDQLVVKWQSSALDVNTWMIWFPDVPDSEPTLSWESVSQATNWTIQDKL
 KPFWCYNI SVY PMLHDKVGEPIYIQAAYAKEGVPSEGPETKVENIGVKTVTITWKEIPKSERKGIICNYTI
 FYQAEQGGKFSKTVNSSILQYGLESLKRKTSYIVQVMASTAGGTNGTSINFKTLFSVFEEIILITSLIG
 GGLLILIIILTVAYGLKPKNKLTHLCWPTVPNPAESSIATWHGDDFKDKLNLKESDDSVNTEDRILKPCST
 PSDKLVIDKLVVNFVNLQEIFTDEARTGQENNLGGEKNGTRILSSCPTS I


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001242637

ORF Size: 2043 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242637.2](#)

RefSeq Size: 2672 bp

RefSeq ORF: 2046 bp

Locus ID: 133396

UniProt ID: [Q8NI17](#)

Cytogenetics: 5q11.2

Protein Families: Druggable Genome, Transmembrane

MW: 77.5 kDa

Gene Summary: The protein encoded by this gene belongs to the type I cytokine receptor family. This receptor, with homology to gp130, is expressed on monocytes, and is involved in IL-31 signaling via activation of STAT-3 and STAT-5. It functions either as a monomer, or as part of a receptor complex with oncostatin M receptor (OSMR). Several alternatively spliced transcript variants encoding different isoforms have been noted for this gene.[provided by RefSeq, Jun 2011]