

Product datasheet for **RC233089**

NUP93 (NM_001242795) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NUP93 (NM_001242795) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NUP93
Synonyms:	NIC96
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC233089 representing NM_001242795
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTGAGGAGTACCATCGGGAGTCAATGTTGGTTGAGTGGGAGCAAGTGAAACAGCGAATTCGCACA
 CACTGCTGGCATCAGGAGAAGACGCCCTTGACTTTACTCAAGAAAGCGAGCCAAGCTACATCAGTGATGT
 GGGACCCCTGGTCGAAGCTCTCTGGATAACATCGAGATGGCCTATGCGCGGCAATTTATATCTATAAT
 GAGAAAATTGTAATGGACACCTGCAGCCTAACCTGGTGGACCTTTGTGCTTCCGTCGCAGAGCTCGATG
 ATAAGAGCATTTCCGACATGTGGACCATGGTAAACAAATGACAGACGTGTTGTTGACACCGGCAACGGA
 TGCCCTGAAGAACCAGCAGCAGCTGGAAGTGGCATGGAGTTTGTGAGGAGGCTTGGCGTACCTTGAG
 CAGAGTTATAAGAATTACACCTTGTGACTGTCTTTGAAATTTGCATCAGGCCAGCTGGGCGGGGTGC
 CTGGGACTTACCAATTGGTTGGAAGTTTCTGAACATTAAGTCCAGCTCCCTTGCCTGGACTACAGGA
 TGGAGAGGTGGAAGGCCATCCTGTGTGGCGCTAATTTACTACTGCATGCGCTGTGGAGACCTGCTTGCC
 GCTTCACAGGTAGTTAATCGAGGCCAGCACAGCTGGGAGAGTTTAAACCTGGTTCCAGGAGTACATGA
 ACAGCAAGGACAGAAGATTGTCCCCAGCTACGGAACAAAGCTCCGGCTGCATTACCGTAGGGCCCTCAG
 GAACAATACAGATCCCTACAAGCGGGCCGTGACTGTATCATTGGCAGATGTGACGTACCGACAACAG
 AGTGAAGTGGCGGACAAACTGAGGATTACCTGTGGTGAAAGTTGAACCAAGTGTGTTTTGACGACGATG
 GCACAGCTCCCAACAAGACAGGCTCACTCTCTCAGATTCAGAAAGCAGTTGTTGGAAGACTATGGCGA
 GTCCCACTTTACGGTGAACAGCAACCTTCTCTACTTCCAAGTCTGTTCTGACAGCGCAGTTTGAA
 GCAGCAGTTGCCCTTTCTTTCCGATGGAGCGGCTGCGCTGCCATGCTGTCCATGTAGCACTGGTGCTGT
 TTGAGCTGAAGCTGCTTTTAAAGTCTCTGGACAGAGTGCTCAGCTCCTCAGCCACGAGCCTGGTGACCC
 TCCTTGCTTGGCGCGGCTGAACCTTCGTGCGGCTCCTCATGCTGTACACCCGGAAGTTTGAGTCCACGGAC
 CCAAGGGAGGCCCTCCAGTACTTCTATTTCTCAGGGATGAGAAAGATAGTCAAGGAGAAAAATGTTTC
 TGCGCTGTGTGAGTGAGCTTGTGATTGAAAGCCGAGAGTTCGATATGATTCTTGGGAACTAGAGAATGA
 CGGAAGTAGAAAGCCTGGAGTCATAGATAAGTTTACTAGTGACACAAAGCCTATTATCAACAAAGTTGCT
 TCTGTGGCAGAAAATAAAGGACTGTTTGAAGAGGCAGCAAGCTGTATGACCTTGCCAAGAATGCTGACA
 AGGTACTGGAGCTGATGAACAACTGCTGAGCCCTGTCGTCCTCCAGATCAGTGCCCGCAATCCAACAA
 GGAGAGGCTGAAGAACATGGCACTCTCCATTGCCGAACGGTATAGGGCTCAAGGAATAAGCGCAATAAA
 TTTGTGGACTCCACGTTCTATCTCTTTTGGACTTGATCACCTTTTTTGACGAGTATCATAGTGGTCATA
 TTGATAGAGCTTTTGATATCATTGAGCGCTTGAAGCTGGTGCCCTGAATCAGGAAAGTGTGGAAGAGAG
 AGTGGCTGCCTTCAGAAATTTAGTGATGAAATCAGGCACAACCTCTCAGAAAGTGCTTCTTGCCACCATG
 AACATCTTGTTACACAGTTTAAAGAGGCTCAAGGGGACAAGTCCATCCTCGTCATCCAGGCCCGCAGGAG
 TCATCGAGGACCGCACTCTCAACTCCGAAGTCAAGCCCGCACTCTGATTACCTTTGCTGGAATGATACC
 ATACCGAACGTCTGGGGACCAATGCGAGGCTGGTGCAGATGGAGGTCTCATGAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC233089 representing NM_001242795
 Red=Cloning site Green=Tags(s)


MAEEYHRESMLVEWEQVKQRILHTLLASGEDALDFTQESEPSYISDVGPPGRSSLDNIEMAYARQIYIYN
 EKIVNGHLQPNLVDLCASVAELDDKSISDMWTMVKQMTDVLTPATDALKNRSSVEVRMEFVRQALAYLE
 QSYKNYTLVTVFGNLHQAQLGGVPGTYQLVRSFLNIKLPAPLPGLQDGEVEGHPVWALIYYCMRCGDLLA
 ASQVVNRAQHQLGEFKTWFEYMNKDRRLSPATENKLRLHYRRALRNNTDPYKRAVYCIIGRCDVTDNQ
 SEVADKTEDYLVKLNQVCFDDGTSSPDRLTLSTQFQKQLLEDYGESHTVNNQPPFLYFQVFLTAQFE
 AAVAFLLFRMERLRCHAVHVALVLFELKLLKSSGQSAQLLSHEPGDPPCLRRLLNFVRLMLYTRKFESTD
 PREALQYFYFLRDEKDSQGENMFLRCVSELVIESREFDMILGKLENDGSRKPGVIDKFTSDTKPIINKVA
 SVAENKGLFEEAAKLYDLAKNADKYLEMKNLLSPVVPQISAPQSNKERLKNMALIAERYAQQISANK
 FVDSTFYLLDLITTFDEYHSGHIDRAFDIERLKLVLNQESVEERVAEFRNFSDIEIRHNLSEVLLATM
 NILFTQFKRLKGTSPSSSRPQRVIEDRDSQLRSQARTLITFAGMIPYRTSGDTNARLVQMEVLMN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001242795

ORF Size: 2088 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001242795.2](#)

RefSeq Size: 2927 bp

RefSeq ORF: 2091 bp

Locus ID: 9688


UniProt ID: [Q8N1F7](#)

Cytogenetics: 16q13

MW: 80.3 kDa

Gene Summary: The nuclear pore complex is a massive structure that extends across the nuclear envelope, forming a gateway that regulates the flow of macromolecules between the nucleus and the cytoplasm. Nucleoporins are the main components of the nuclear pore complex in eukaryotic cells. This gene encodes a nucleoporin protein that localizes both to the basket of the pore and to the nuclear entry of the central gated channel of the pore. The encoded protein is a target of caspase cysteine proteases that play a central role in programmed cell death by apoptosis. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2016]

Product images:



Circular map for RC233089