

Product datasheet for **RC233266A1V**

Human Insulin (INS) (NM_001185098) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Insulin (INS) (NM_001185098) AAV Particle
Tag: Myc-DDK
Symbol: Insulin
Synonyms: IDDM; IDDM1; IDDM2; ILPR; IRDN; MODY10; PNDM4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC233266 representing NM_001185098
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCCTGTGGATGCGCCTCCTGCCCTGCTGGCGCTGCTGGCCCTCTGGGGACCTGACCCAGCCGCAG
CCTTTGTGAACCAACACCTGTGCGGCTCACACCTGGTGAAGCTCTACCTAGTGTGCGGGGAACGAGG
CTTCTTCTACACACCAAGACCCGCGGGAGGCAGAGGACCTGCAGGTGGGCAGGTGGAGCTGGCGGG
GGCCCTGGTGCAGGCAGCCTGCAGCCCTTGGCCCTGGAGGGTCCCTGCAGAAGCGTGGCATTGTGGAAC
AATGCTGTACCAGCATCTGCTCCCTCTACCAGCTGGAGAACTACTGCAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC233266 representing NM_001185098
Red=Cloning site Green=Tags(s)

MALWMRLPLLLALLALWGPDPAAAFVNQHLGSHLVEALYLVCGERGFFYTPKTRREAEDLQVQVELGG
GPGAGSLQPLALEGSLQKRGIVEQCCTSICSLYQLENYCN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001185098
ORF Size: 330 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001185098.1</u>
RefSeq Size:	648 bp
RefSeq ORF:	333 bp
Locus ID:	3630
UniProt ID:	<u>P01308</u>
Cytogenetics:	11p15.5
MW:	12 kDa