

## Product datasheet for **RC233275**

### **TMEM218 (NM\_001258247) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TMEM218 (NM\_001258247) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** TMEM218  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC233275 representing NM\_001258247  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTGGCACTGTGCTCGGAGTCGGTGCGGGCGTGTTTCATCTTAGCCCTGCTCTGGGTGGCAGTGCTGC  
TGCTGTGTGTGCTGCTGTCCAGAGCCTCCGGGGCGGCGAGGTTCTCTGTCATTTTTTTATTCTTCGGTGC  
TGTGATCATCACATCAGTTCTGTTGCTTTCCCGCAGCTGGTGAATCCAGCCCCAGAAGTGGAAAGTT  
AAGATTGGATGACTTTTTTCATTGGCCGCTATGCTCTGCTGGCTTTCCTTAGTGCCATCTTCCTTGGAG  
GCCTCTTCTGGTTTTAATCCATTATGTTCTGGAGCCGATCTATGCCAAACCACTGCACCTCTAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC233275 representing NM\_001258247  
**Red=Cloning site Green=Tags(s)**

MAGTVLGVGAGVFILALLWVAVLLLCVLLSRASGAARFSVIFLFFGAVIITSVLLLFPRAGEFPAPEVEV  
KIVDDFFIGRYVLLAFLSAIFLGLFLVLIHYVLEPIYAKPLHSY

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001258247

**ORF Size:** 345 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001258247.2](#)

**RefSeq Size:** 3891 bp

RefSeq ORF: 348 bp  
 Locus ID: 219854  
 UniProt ID: [A2RU14](#)  
 Cytogenetics: 11q24.2  
 Protein Families: Transmembrane  
 MW: 12.9 kDa  
 Gene Summary: May be involved in ciliary biogenesis or function.[UniProtKB/Swiss-Prot Function]

**Product images:**


Circular map for RC233275