

Product datasheet for RC233315

JPH3 (NM_001271605) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: JPH3 (NM_001271605) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: JPH3
Synonyms: CAGL237; HDL2; JP-3; JP3; TNRC22
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC233315 representing NM_001271605
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCCAGTGGGGCAGGTTTAATTTTGACGACGGAGGGTCTACTGTGGAGGCTGGGAGGACGGCAAGG
 CGCACGGCCATGGCGTCTGCACCGGCCCAAGGGCCAAGGCGAATACACCGGCTCGTGGAGCCACGGCTT
 CGAGGTGCTGGGCGTCTACACCTGGCCAGCGGCAACACGTACCAGGGCACCTGGCGCAGGGCAAGCGC
 CACGGCATCGGCCTGGAGAGCAAGGGGAAGTGGGTGTACAAGGGCGAGTGGACGCACGGATTCAAGGGG
 GCTACGGGTGCGGGAGTGCGCGGGCAACGGGGCAAATACGAAGGGACCTGGAGCAACGGGCTGCAGGA
 CGGCTACGGGACCGAGACCTACTCGGACGGAGGTGGTTGGACGGGACTAGCGCTGTCAAGCAGGCTGC
 TGGAGGGAGGAGGGAGGCCCATCTGCTCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC233315 representing NM_001271605
 Red=Cloning site Green=Tags(s)

MSSGGRFNFDDGGSYCGGWEDGKAHGHVCTGPKGQGEYTGWSHGFVGLGVYTWPSGNTYQGTWAQGKR
 HGIGLESKWKWYKGEWTHGFKGRYGVRECAGNGAKYEGTWSNGLQDGYGTETYSDDGGWRTSAVEAGC
 WREEGGPSAQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001271605

ORF Size: 450 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271605.2](#)

RefSeq Size: 1238 bp

RefSeq ORF: 453 bp

Locus ID: 57338


Cytogenetics: 16q24.2

Protein Families: Druggable Genome, Transmembrane

MW: 16.5 kDa

Gene Summary: Junctional complexes between the plasma membrane and endoplasmic/sarcoplasmic reticulum are a common feature of all excitable cell types and mediate cross talk between cell surface and intracellular ion channels. The protein encoded by this gene is a component of junctional complexes and is composed of a C-terminal hydrophobic segment spanning the endoplasmic/sarcoplasmic reticulum membrane and a remaining cytoplasmic domain that shows specific affinity for the plasma membrane. CAG/CTG repeat expansion from normally 6-28 repeats to 40-59 repeats in the 3' UTR of this gene have been associated with Huntington disease-like 2 (HDL2). This gene is a member of the junctophilin gene family. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Jul 2016]

Product images:



Circular map for RC233315