

Product datasheet for RC233317

C1orf144 (SZRD1) (NM_001271869) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C1orf144 (SZRD1) (NM_001271869) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C1orf144
Synonyms: C1orf144
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC233317 representing NM_001271869
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAAGATGAGGAGGTCGCTGAGAGCTGGGAAGAGGCGGCAGACAGCGGGAAATAGACAGACGGTTGG
 AAAAAAAAACTGAAGATCACACAAAAGAGAGGAAATCAAATCTCTCCCAAAGTCCCCATTGTGATTCA
 GGACGATAGCCTTCCCGCGGGCCCCCTCCACAGATCCGCATCTCAAGAGGCCACCAGCAACGGTGTG
 GTCAGCAGCCCCAACTCCACCAGCAGGCCACCCTTCCAGTCAAGTCCCTAGCACAGCGAGAGGCCGAGT
 ACGCCGAGGCCCGAAGCGGATCTGGGCAGCGCCAGCCCGAGGAGGAGCAGGAGAAACCCATCTCGA
 CAGGCCAACAGGATCTCCAACCCGAAGACAGCAGGCAGCCCAATAATGTGATCAGACAGCCTTTGGGT
 CCTGATGGGTCTCAAGGCTTCAAACAGCGCAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC233317 representing NM_001271869
 Red=Cloning site Green=Tags(s)
 MEDEEVAESWEEAADSGEIDRRLEKCLKITQKERKSKSPPKVPIVIQDDSLPAGPPPQIRILKRPTSNV
 VSSPNSTRPTLPVKSLAQREAEYAEARKRILGSASPEEEQEKPLDRPTRISQPEDSRQPNNVIRQLPG
 PDGSQGFKQRR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:

ACCN: NM_001271869

ORF Size: 453 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)
OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271869.2](#)
RefSeq Size: 3692 bp

RefSeq ORF: 456 bp


Locus ID: 26099

UniProt ID: [Q7Z422](#)

Cytogenetics: 1p36.13

MW: 17.4 kDa

Product images:



Circular map for RC233317