

Product datasheet for **RC233332**

DBNDD2 (NM_001197140) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: DBNDD2 (NM_001197140) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: DBNDD2
Synonyms: C20orf35; CK1BP; HSMNP1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC233332 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

RCATGGACCCAAATCCTCGGGCCGCCCTGGAGCGCCAGCAGCTCCGCCTTCGGGAGCGGCAAAAATTCTT
CGAGGACATTTTACAGCCAGAGACAGAGTTTGTCTTCTCTGTCCCATCTGCATCTCGAGTCGCAGAGA
CCCCCATAGGTAGTATCTCATCCATGGAAGTGAATGTGGACACACTGGAGCAAGTAGAACTTATTGACC
TTGGGGACCCGGATGCAGCAGATGTGTTCTTGCCTTGCGAAGATCTCCACCAACCCCCAGTCGCTGG
GGTGGACAACCATTTGGAGGAGCTGAGCCTGCCGGTGCCTACATCAGACAGGACCACATCTAGGACCTCC
TCCTCCTCCTCCGACTCCTCCACCAACCTGCATAGCCCAAATCCAAGTGATGATGGAGCAGATACGC
CCTTGGCACAGTCGGATGAAGAGGAGGAAAGGGTGATGGAGGGGCAGAGCCTGGAGCCTGCAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC233332 protein sequence
Red=Cloning site Green=Tags(s)

XWTQILGPPWSASSAFGSGKNSRTFYRQSLSFCLPICISSRRDPP*VVSHPWK*MWTHWSK*NLTT
LGTMRMQMCSCLAKILHQPPSRLGWTTIWR*ACRCLHQTGPLHPPPPPPPTPPPTCIAIQVMMEQIR
PWHSRMKRRKGVMEQSLEPA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6384_e10.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001197140

ORF Size: 483 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001197140.1](#), [NP_001184069.1](#)


RefSeq Size: 1044 bp

RefSeq ORF: 486 bp

Locus ID: 55861

UniProt ID:	<u>Q9BQY9</u>
Cytogenetics:	20q13.12
MW:	17.5 kDa
Gene Summary:	May modulate the activity of casein kinase-1. Inhibits CSNK1D autophosphorylation (in vitro). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC233332