

Product datasheet for RC233337

EIF1AD (NM_001242482) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: EIF1AD (NM_001242482) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: EIF1AD
Synonyms: haponin; OBELIX
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC233337 representing NM_001242482.
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTCAGGCCACCAAGAGGAAGCATGTGGTGAAGGAGGTGCTAGGGGAGCACATAGTGCCTCCAAC
CAGCAGCAGATTGTCAGGGTACTCAGGACCCAGGGAACAATCTGCATGAGGTGGAGACAGCCCAAGGG
CAGCGCTTCTGGTGAGCATGCCCTCCAATACCGCAAGAACATCTGGATCAAGAGAGGGGACTTTCTC
ATTGTTGACCCATTGAAGAGGGAGAAAAGGTGAAGGCTGAAATCTCGTTTGTGCTCTGCAAGGACCAC
GTGCGCTCTCTGCAGAAGGAGGGTTTTGGCTGAGGCCTTCTCTGAAGTGGCTGAGAAACACAACAAC
AGGAACAGACAACTCAACCAGAACTCCAGCTGAGCCACAGTTATCAGGAGAGGAGTCCAGCTCAGAA
GATGATTCGACCTGTTTGTAAACAAAACCGCAGACAGTATCATGAGAGTGAGGAGGAGAGTGAAGAG
GAGGAGGCAGCC
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Protein Sequence: >Peptide sequence encoded by RC233337
Blue=ORF Red=Cloning site Green=Tag(s)

```
MSQATKRKHVVKEVLGEHIVPSNQQQIVRVLRTPGNNLHEVETAQGQRFLVSMPSKYRKNIWIKRGDFL
IVDPIEEGKVKAEISFVLCKDHVRSQKEGFWPEAFSEVAEKHNNRNRQTQPELPAEPQLSGEESSE
DDSDLFVNTNRRQYHESEEESEEEEA
TRTRPLEQKLISEEDLAANDILDYKDDDDKV
```

Chromatograms: https://cdn.origene.com/chromatograms/mk6405_e07.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001242482

ORF Size: 495 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq Size: 2901 bp

RefSeq ORF: 498 bp

Locus ID: 84285


UniProt ID: [Q8N9N8](#)

Cytogenetics: 11q13.1

MW: 19.1 kDa

Gene Summary: Plays a role into cellular response to oxidative stress. Decreases cell proliferation.
[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC233337