

Product datasheet for **RC233361**

Ankyrin repeat domain containing protein 65 (ANKRD65) (NM_001243535) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ankyrin repeat domain containing protein 65 (ANKRD65) (NM_001243535) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ankyrin repeat domain containing protein 65
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233361 representing NM_001243535 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACTCCCAGAGGCCTGAGCCCAGAGAGGAGGAGGAGGAACAGGAAGTCCGGTGGATGGAGCTGG
ACTCCGAAGAGGCCCTGGGAACAGGACAGAGGGGCTAGTGTTCAGGGCTGGGGGCACCTGCTCCA
GGCCGTGTGGAGGGCCCTGCAGGCCTGGTGACGACGTGCTGCGCAAGGTGCCAGCGTGGAGGAGAGG
ACATTGAGGTGCTGCTGGGCCACGGGCAGACCCAGGCATCAGGGACAGGCATGGCCGCTCTGCGCTGCA
CAGGGCTGCCGCCGAGGACACCTGCTTGCCGTCCAGTTGCTGGTCACCCAGGGGGCCGAGGTGGATGCG
CGGGACACCCTGGGCCTCACACCCTGCATCACGCCTCTCGGGAAGGCCACGTGGAGGTTGCCGGCTGCC
TGCTGGACAGGGGTGCCAGGTGGATGTACCGCTGGCTCCGAAAGACCCCTACACCTGGCTGCAGA
GCGAGGGCATGGGCCTACCGTGGGGCTTCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC233361 representing NM_001243535 Red=Cloning site Green=Tags(s)
-------------------	---

MDSQRPEPREEEEEQELRWELDSEELGTRTEGPSVVQGWGHLLQAVWRGPAGLVTLQLLRQASVEER
TLRCCWATGQTQASGTGMAALRCTGLPPEDTCLPSSCWSRGPWRMRGTPWASHPCITPLGKATWRLPAA
CWTGVPRWMLPAGSERPPYTWLQSEGMGLPWGFC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM 001243535

ORF Size: 522 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001243535.2](#)

RefSeq Size: 1664 bp

RefSeq ORF: 525 bp

Locus ID: 441869
UniProt ID: [E5RIM6](#)
Cytogenetics: 1p36.33
MW: 19.8 kDa

Product images:



Circular map for RC233361