

Product datasheet for RC233399

C1orf149 (MEAF6) (NM_001270875) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C1orf149 (MEAF6) (NM_001270875) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C1orf149
Synonyms:	C1orf149; CENP-28; EAF6; NY-SAR-91
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233399 representing NM_001270875 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGCGATGCACAACAAGGCGGCGCCGCCGAGATCCCGACACCCGGCGGGAGCTGGCGGAGCTCGTGA
 AGCGGAAGCAGGAGCTGGCGGAAACATTGGCAAATTTGGAGCGACAGATCTATGCTTTTGAGGGAAGCTA
 CCTGGAAGACACTCAGATGTATGGCAATATTATTCGTGGCTGGGATCGGTATCTGACCAACCAAAAAAC
 TCCAATAGCAAAAATGATCGAAGGAACCGAAGTTTAAGGAAGCTGAGCGGCTCTTCAGTAAATCCTCGG
 TTACCTCAGCAGCTGCAGTAAGTGCATTGGCAGGAGTTCAGGACCAGCTCATTGAAAAGAGGGAGCCAGG
 AAGTGGGACGGAAGTGACACTTCTCCAGACTCCACAATCAGGAAAATGAGCCAGCCAGGAGGACCCT
 GAGGATCTGGATGGATCTGTGCAGGAGTGAAACCTCAGAAGGCTGCTTCTTCTACTTCTCAGGGAGTC
 ACCACAGCAGCCATAAAAAGCGAAAGAATAAAACCGGCACAGGATTGATCTGAAGTTAAACAAAAAAC
 ACGAGCTGACTAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence:	>RC233399 representing NM_001270875 Red=Cloning site Green=Tags(s)
	MAMHNKAAPPQIPDTRRELAELVKRKQELAETLANLERQIYAFEGSYLEDQMYGNIIRGWDRYL TNQKN SNSKNDRNRKFKEAERLFSKSSVTSAAAVSALAGVQDQLIEKREPGSGTESDTSPDFHNQENEPSQEDP EDLDGSVQGVKPKAASSTSSGSHSHSHKKRKNKNRHRIDLKLNKKPRADY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_001270875

ORF Size: 573 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001270875.3](#)

RefSeq Size: 4788 bp

RefSeq ORF: 576 bp

Locus ID: 64769


UniProt ID: [Q9HAF1](#)

Cytogenetics: 1p34.3

MW: 22.1 kDa

Gene Summary: This gene encodes a nuclear protein involved in transcriptional activation. The encoded protein may form a component of several different histone acetyltransferase complexes. There is a pseudogene for this gene on chromosome 2. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012]

Product images:



Circular map for RC233399