

Product datasheet for **RC233421**

ZFP41 (NM_001271156) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZFP41 (NM_001271156) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ZFP41
Synonyms: zfp-41; ZNF753
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC233421 representing NM_001271156
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGAAGCCTGCAGGCAGAAAAAGAAGACGCCGACCCCAAGGGAGGAGGCAGACGTGCAGAAGAGTG
CGCTCAGAGAGGAGAAGGTGTCCGGGGACAGAAAGCCACCTGAGAGGCCCACTGTGCCAGGAAGCCCCG
CACAGAGCCCTGCCTGAGTCCTGAAGACGAAGAGCACGTCTTTGATGCCTTCGACGCTTCATTTAAAGAT
GACTTTGAGGGGTTCCCGTGTTCACTCTTTTCAGAGGAAGAAACCCTATGAGTGCAGTGAGTGTGGGC
GGATCTTTAAGCACAAGACAGACCACATTCGCCATCAGAGGGTCCACACTGGAGAGAAGCCCTTCAAGTG
TGCAGTGCAGGGAAGGCCTTCGGGCACAGCTCTGACGTACCAAACACCAGAGGACTCACACGGGAGAG
AAGCCCTTCAAATGCGGGGAGTGCGGGAAGCCTTTAACTGCGGCTCAAATCTCCTGAAACATCAGAAGA
CGCACACCGGGGAGAAGCCCTACGAATGCACGCACACTGTGGGAAAGCCTTTGCCTACAGCTCCTGTCTCAT
CCGCCATCAGAAACGCCACCCTCGGAAGAAGCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC233421 representing NM_001271156
Red=Cloning site Green=Tags(s)

MEKPAGRKKKTPTPREADVQKSALREEKVSQDRKPPRPTVPRKPRTEPCLSPEDDEHVFDADFASFKD
DFEGVPVFIPIFQRKKPYECSEGRIFKHKTDIRHQRVHTGEKPFKCAQCGKAFRHSDDVTKHQRTHTGE
KPFKCGECGKAFNCGSNLLKHQKTHTEKPYECTHCGKAFAYSSCLIRHQKRHPKKP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001271156

ORF Size: 594 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001271156.3](#)

RefSeq Size: 3298 bp

RefSeq ORF: 597 bp

Locus ID:	286128
UniProt ID:	Q8N8Y5
Cytogenetics:	8q24.3
MW:	23.3 kDa
Gene Summary:	A putative DNA-binding regulatory protein associated with meiosis in spermatogenesis. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC233421