

Product datasheet for RC233424

STX10 (NM_001271611) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	STX10 (NM_001271611) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	STX10
Synonyms:	hsyn10; SYN10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233424 representing NM_001271611 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTCTCGAAGACCCCTTTTTGTAGTCCGAGGCGAGGTGCAGAAGGCGGTGAACACGGCCCGCGGGC
TGTACCAGCGCTGGTGCAGCTCCTGCAGGAAAGCGCGGCGTCCGACGCGAGGAGCTGGACTGGACGAC
CAATGAGCTGCGGAATGGCTGCGCAGCATCGAGTGGACCTCGAGGACCTGGAAGAGACCATCGGTATA
GTGGAAGCCAACCCAGGCAAGCCAGCTGCCAGAAGTCAACAGCGACCTGCTGGATGCCAGCGCAGTCT
CGGCCACATCTCGCTACATCGAGGAGCAGCAGGCCACACAGCAGCTGATCATGGATGAACAGGATCAACA
GCTGGAGATGGTGTCTGGGAGCATCCAGTTCTGAAGCACATGTCCGGCCGCGTTGGAGAAGAGCTGGAC
GAGCAGGGCATCATGCTGGATGCCTTCGCCAAGAGATGGACCACACCCAGTCCCGCATGGACGGGTCC
TCAGGAAGTTGGCAAAGTATCCACATGACGAGTGACCGCCGACAGTGGTGTGCCATCGCCGTGCTAGT
GGGGGTGCTTCTCCTCGTTCTCATCTTACTATTCTCTCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC233424 representing NM_001271611 Red=Cloning site Green=Tags(s)
-------------------	---


MSLEDPFVVRGEVQKAVNTARGLYQRWCELLQESAAGVREELDWTNENLRNGLRSIEWDLEDLEETIGI
VEANPGKPAAQKSPSLLDASAVSATSRYIEEQATQQLIMDEQDQQLMVSQSIQVLKHMSSGRVGEELD
EQGIMLDAFAQEMDHTQSRMDGVLRLAKVSHMTSDRRQWCAIAVLVGVLLLVLILLFSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_001271611

ORF Size: 600 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001271611.2](#)

RefSeq Size: 1173 bp

RefSeq ORF: 603 bp

Locus ID:	8677
UniProt ID:	O60499
Cytogenetics:	19p13.13
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	SNARE interactions in vesicular transport
MW:	22.9 kDa
Gene Summary:	This gene belongs to the syntaxin family and encodes a soluble N-ethylmaleimide sensitive factor attachment protein receptor (SNARE). The encoded protein is involved in docking and fusion events at the Golgi apparatus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2012]

Product images:



Circular map for RC233424