

Product datasheet for RC233427

RhoGDI (ARHGDI) (NM_001185077) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RhoGDI (ARHGDI) (NM_001185077) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: RhoGDI
Synonyms: GDIA1; HEL-S-47e; NPHS8; RHOGDI; RHOGDI-1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC233427 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGCTGAGCAGGAGCCACAGCCGAGCAGCTGGCCAGATTGCAGCGGAGAACGAGGAGGATGAGCACT
 CGGTCAACTACAAGCCCCGGCCAGAAGAGCATCCAGGAGATCCAGGAGCTGGACAAGGACGACGAGAG
 CCTGCGAAAGTACAAGGAGGCCCTGCTGGGCCGCGTGGCCGTTTCCGAGACCCCAACGTCCCAACGTC
 GTGGTACTGGCCTGACCCTGGTGTGCAGCTCGGCCCGGGCCCCCTGGAGCTGGACCTGACGGGCGACC
 TGGAGAGCTTCAAGAAGCAGTCGTTTGTGCTGAAGGAGGGTGTGGAGTACCGGATAAAAAATCTCTTCCG
 GGTTAACCAGAGATAGTGTCCGCATGAAGTACATCCAGCATACGTACAGGAAAGGCGTCAAGATTGAC
 AAGACTGACTACATGGTAGGCAGCTATGGGCCCGGGCCGAGGAGTACGAGTTCCTGACCCCGTGGAGG
 AGGCACCAAGGGTATGCTGGCCCGGGGACGCTACAGCATCAAGTCCCGCTTACAGACGACGACAAGAC
 CGACCACCTGTCTGGGAGTGGAAATCTACCATCAAGAAGGACTGGAAGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC233427 protein sequence
 Red=Cloning site Green=Tags(s)

MAEQEPTAEQLAQIAAENEDEHSVNYKPPAQKSIQEIQLDKDDESLRKYKEALLGRVAVSADPNV
 VVTGLTLVCSSAPGPLELDLTGDLESFVKQSFVLEKEGVEYRIKISFRVNREIVSGMKYIQHTYRKGVKID
 KTDYMGVSYGPRAEYEFLLTPVEEAPKMLARGSYSIKSRFTDDDKTDHLSWEWNLTIKDKWDK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Chromatograms: https://cdn.origene.com/chromatograms/mk6085_d09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001185077

ORF Size: 612 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001185077.3](#)

RefSeq Size: 2150 bp


RefSeq ORF: 615 bp

Locus ID: 396

UniProt ID: [P52565](#)
Cytogenetics: 17q25.3
Protein Families: Druggable Genome
Protein Pathways: Neurotrophin signaling pathway
MW: 23.2 kDa

Gene Summary: This gene encodes a protein that plays a key role in the regulation of signaling through Rho GTPases. The encoded protein inhibits the disassociation of Rho family members from GDP (guanine diphosphate), thereby maintaining these factors in an inactive state. Activity of this protein is important in a variety of cellular processes, and expression of this gene may be altered in tumors. Mutations in this gene have been found in individuals with nephrotic syndrome, type 8. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]

Product images:



Circular map for RC233427