

Product datasheet for **RC233452**

PAQR6 (NM_001272111) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PAQR6 (NM_001272111) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PAQR6
Synonyms:	PRdelta
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233452 representing NM_001272111 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTACCATCTCTTCTGCGCGCTGCTCACTGGCTTCTCTTCGCCTCCACCTGCCTGAAAGGCTGGC
ACCAGGACGCTTTGATTACATCGGCCACAGCCACCAAGTTATTCCACATCTGTGCAGTGCTGGGCACCCAC
TTCCAGCTGGAGGCAGTGCTGGCTGATATGGGATCACGCAGAGCCTGGCTGGCCACACAGGAACCTGCC
TGGGCTGGCAGGCACAGTGCCCACTGGTCTTGGCTGCAGCTGGGAACCTACTCATTATTGCTGCTTT
CACAGCCACCCTGCTTCGGGCCCCAGTACATGCCCTCTGCTGCAGGGTGGCCCACTGGAGGGGGTACC
CAGGCCAAACAACAGTGAGGCCCATCCCTGACCCTGTCCTGGAGGGGGCAGAGGCCAGGCCCACTGCT
GACGAGGAGCCCAGATTTGGGCTAATCAGGTGGGACGCATCTCAGCCTGGAACCAACAGGGGCTGAGG
AGAGAGGGCAGAGGAGAGGGCAGAGAAGAGGAGGGGTGTCTAGGGGGACTGGCAGAGTGTGAGAGGGA
CCGTGAGGGGGCTCTTGATGGGAGTGGAAGAAGTGTGAGGGTCTGAGAGGGGAGATGCATGCGTGTCCA
GGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC233452 representing NM_001272111

Red=Cloning site Green=Tags(s)

MATISSARCSLASSPPTCLKGWHQDALITSATATSYSTSVQCWAPTSSWRQCWLIDWHAEPGWPHRNL
 WAWQAQWPHWSWLQLGTYSLLLSQPPCFGPPVHALCCRVAHWRGVPRPNNSEAPSLTLWRGQRPGPSA
 DEEPRFGPNQVGTHLSLEPTGAEEGRHRREGREEEGCLGGLAECERDREGALDGSGRSAEGLRGEMHACP
 G



TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM 001272111

ORF Size: 633 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001272111.2</u>
RefSeq Size:	1902 bp
RefSeq ORF:	636 bp
Locus ID:	79957
UniProt ID:	<u>Q6TCH4</u>
Cytogenetics:	1q22
Protein Families:	Druggable Genome, Transmembrane
MW:	23.7 kDa
Gene Summary:	Plasma membrane progesterone (P4) receptor coupled to G proteins (PubMed:23763432, PubMed:23161870). Seems to act through a G(s) mediated pathway (PubMed:23161870). Involved in neurosteroid inhibition of apoptosis (PubMed:23161870). May be involved in regulating rapid P4 signaling in the nervous system (PubMed:23763432). Also binds dehydroepiandrosterone (DHEA), pregnanolone, pregnenolone and allopregnanolone (PubMed:23763432, PubMed:23161870).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC233452