

Product datasheet for RC233455

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

EMA (MUC1) (NM_001204296) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: EMA (MUC1) (NM 001204296) Human Tagged ORF Clone

Tag: Myc-DDK

Symbol: EMA

Synonyms: ADMCKD; ADMCKD1; ADTKD2; CA 15-3; CD227; EMA; H23AG; KL-6; MAM6; MCD; MCKD;

MCKD1; MUC-1; MUC-1/SEC; MUC-1/X; MUC1/ZD; PEM; PEMT; PUM

Mammalian Cell

Selection:

Neomycin

Vector:pCMV6-Entry (PS100001)E. coli Selection:Kanamycin (25 ug/mL)

ORF Nucleotide >RC233455 representing NM_001204296

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA





Protein Sequence: >RC233455 representing NM_001204296

Red=Cloning site Green=Tags(s)

MTPGTQSPFFLLLLLTVLTATTAPKPATVVTGSGHASSTPGGEKETSATQRSSVPSSTEKNAIPAPTTTK SCRETFLKCFCRFINKGVFWASPILSSVSDVPFPFSAQSGAGVPGWGIALLVLVCVLVALAIVYLIALAV CQCRRKNYGQLDIFPARDTYHPMSEYPTYHTHGRYVPPSSTDRSPYEKVSAGNGGSSLSYTNPAVAATSA NL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:



CTATAGGGCGGCCG	EcoR I GGAATTCGTC	BamHI GACTGGATOO	<u> </u>	RBS_ AGGAGAT		Kozac onsensu Sg	f/	C ATG		
ORF	NNN		Not a	CCG CTC	GAG CA	G AAA K	Myc.1	Tag ATC TC I S	A GAA E	GAG E
GAT CTG GCA GC		RV ATC CTG G I L <u>I</u>	AT TAC A	ag.Tag AG GAT K D	GAC GAC	GAT D		Pme I	ACGG	se I cagaac

^{*} The last codon before the Stop codon of the ORF

ACCN: NM_001204296

ORF Size: 636 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM 001204296.2

 RefSeq Size:
 1037 bp

 RefSeq ORF:
 639 bp

 Locus ID:
 4582

 UniProt ID:
 P15941

 Cytogenetics:
 1q22

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

MW: 22.8 kDa


Gene Summary: This gene encodes a membrane-bound protein that is a member of the mucin family. Mucins

are O-glycosylated proteins that play an essential role in forming protective mucous barriers on epithelial surfaces. These proteins also play a role in intracellular signaling. This protein is expressed on the apical surface of epithelial cells that line the mucosal surfaces of many different tissues including lung, breast stomach and pancreas. This protein is proteolytically cleaved into alpha and beta subunits that form a heterodimeric complex. The N-terminal alpha subunit functions in cell-adhesion and the C-terminal beta subunit is involved in cell signaling. Overexpression, aberrant intracellular localization, and changes in glycosylation of this protein have been associated with carcinomas. This gene is known to contain a highly polymorphic variable number tandem repeats (VNTR) domain. Alternate splicing results in

multiple transcript variants.[provided by RefSeq, Feb 2011]



Product images:



Circular map for RC233455