

## Product datasheet for **RC233463**

### Shugoshin (SGO1) (NM\_001199257) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Shugoshin (SGO1) (NM_001199257) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Shugoshin
Synonyms:	CAID; NY-BR-85; SGO; SGOL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233463 representing NM_001199257 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCAAGGAAAGATGCCTGAAAAAGTCCTTTCAAGATAGTCTTGAAGACATAAAGAAGCGAATGAAAG  
AGAAAAGGAATAAAAACCTTGGCAGAGATTGGCAAACGCAGGTCTTTTATAGCTGCACCATGCCAAATAAT  
CACCAACACTTCTACACTGCTGAAAAATTACCAAGACAACAACAAAATGTTAGTTTTAGCTTTGGAAAAT  
GAAAAATCCAAAGTGAAAGAAGCCCAAGATATCATCCTACAGCTGAGAAAAGAATGTTACTATCTCACAT  
GTCAGCTATATGCATTGAAAGGAAAACCTACATCACAACAACAGTAGAACCTGCTCAGAACCAGGAAAT  
ATGTTCTCTGGAATGGACCCCAATAGTGATGACAGCTCCAGAAATTTATTTGTGAAGGATTTACCGCAA  
ATTCTCTTGAAGAAAACCTGAACTTCCAGGACAAGGAGAATCATTTCAAATAGAAGATCAGATACCTACTA  
TTCCTCAAGACACACTGGGAGTTGATTTTGATTGAGGAGCCCTGGAGGTATCACCTGCCAAAGAAGCAAT  
TTTTATTTTATATTATGTTTCGAGAATTTGTTTCGAGATTTCCAGACTGTAGGAAATGTAACCTGAAACC  
CACATCTGCTTGAGG

**ACGCGT**ACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001199257.3](#)

**RefSeq Size:** 1454 bp

**RefSeq ORF:** 648 bp

**Locus ID:** 151648

**UniProt ID:** [Q5FBB7](#)


**Cytogenetics:** 3p24.3

**Protein Pathways:** Oocyte meiosis

**MW:** 25.1 kDa

**Gene Summary:** The protein encoded by this gene is a member of the shugoshin family of proteins. This protein is thought to protect centromeric cohesin from cleavage during mitotic prophase by preventing phosphorylation of a cohesin subunit. Reduced expression of this gene leads to the premature loss of centromeric cohesion, mis-segregation of sister chromatids, and mitotic arrest. Evidence suggests that this protein also protects a small subset of cohesin found along the length of the chromosome arms during mitotic prophase. An isoform lacking exon 6 has been shown to play a role in the cohesion of centrioles (PMID: 16582621 and PMID:18331714). Mutations in this gene have been associated with Chronic Atrial and Intestinal Dysrhythmia (CAID) syndrome, characterized by the co-occurrence of Sick Sinus Syndrome (SSS) and Chronic Intestinal Pseudo-obstruction (CIPO) within the first four decades of life (PMID:25282101). Fibroblast cells from CAID patients exhibited both increased cell proliferation and higher rates of senescence. Pseudogenes of this gene have been found on chromosomes 1 and 7. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2015]

Product images:



Circular map for RC233463