

Product datasheet for **RC233572**

TMEM99 (NM_001195387) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	TMEM99
Synonyms:	MGC21518
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC233572 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGTTGGCATCCTCCCACTCTGTTGCTCCGGCTGTGTCCCCTCGCTCTGTTGTTCCAGCTATGTGCCCT
CTGTTGCTCCAAGTGCAGCTCATTCTGTTAGAGTTCCTCATTAGCTGGTCACTGTGGCCAGAGGGTGT
GGCCTGCTCCCTTCCTCAAGTATTCTTAAAGCCATGGATTTTGTGGAGCATTTTCTTCTGGCTCTCC
CTTGAGTTATTTCTTTCTTCGCTATCTTGGGACTCTTCTTGTGCTTGGGTCATCGGTTGAGAGAAG
GACTTCTTCTTCTTGTCTCCTTGGTGTGGCTCGTGGTTGCTCTTCAACAAGTGGAGTGGAGGCTCTTG
GTTTTCTTCTATCTTCAACAAGTCACTCTCTCAAGGTCTCACGTTGCAGCATTCTTACCAGAGGCC
ATTGGGCTGGAGTTCAGTTCAGTGTCTGGAGAGTCCACCTCAGCTCAGCAATCTCATGCCGTTGGC
AATTGTCAGCAGAAGCCGATGCCTGCCATCAGTCTTTACTCTGAGGTGTAGAGTGAATAAAATAT
AAATACTTATACTAGTTTTCATGACTTCTGCTTAATATTGGGTATTTTTGTTTTGTTTTGTTGGCG
GTGATAGGCTTACCTTACATTAACAGGCCTTAGCCTTTCTGTGGCTTTGTTATGGCAAAGCCTCATAT
TACTCTCTAGTCTGGTTCAGCAGGACAGTCAGGTCCACCTGGGGCTGTTGTTTTCTACGTTTACCTC
AACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>RC233572 protein sequence
 Red=Cloning site Green=Tags(s)

MVGILPLCCSGCVPSLCCSSYVPSVAPTAHHSVRVPHSAGHCGQRLACSLPQVFLKPWIFVEHFSSWLS
 LELFSFLRYLGTLLCACGHRLEGLLLPCLLGVGSWLLFNNWTGGSWFLHLQQVSLSQGSHVAAFLPEA
 IGPGVPVPVSGESTSAQQSHAGWQLSAEADACPSVLYSEVLEWKNKINTYTSFHDCLILGIFLCFVLA
 VIGLPYIKPGLSLSVALLWQSLILLSSLVQQDSQVHTWGCLFSTFTST

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk8012_d06.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001195387

ORF Size:

774 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001195387.1](#), [NP_001182316.1](#)

RefSeq Size: 2067 bp

RefSeq ORF: 777 bp


Locus ID: 147184

Cytogenetics: 17q21.2

Protein Families: Secreted Protein, Transmembrane

MW: 28 kDa

Product images:



Circular map for RC233572