

Product datasheet for **RC233587**

FRMD6 (NM_001267047) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FRMD6 (NM_001267047) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	FRMD6
Synonyms:	C14orf31; c14_5320; EX1; Willin
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233587 representing NM_001267047 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGGACCAGCTGGAAAAACGGTCGCGGGCCAGCGGGAGCAGTGCGGGCAGCATGAAACACAAGCGCCTGT
 CCCGTATTCCACCGCCAGCCACAGCAGTTCACACCTCGGGCATTGAGGCAGACCAAGCCCCGGGA
 CACGGGGCCAGAAGACAGCTACTCCAGCAGTCCATCCACGCAAGCTGAAAACCTGCAGCTCAATGACC
 AGTCATGGCAGCTCCACACCTCAGGGGTGGAGAGTGGCGGCAAGACCGGCTGGAAGAGGACTTACAGG
 ACGATGAAATAGAGATGTTGGTTGATGACCCCGGGATCTGGAGCAGATGAATGAAGAGTCTCTGGAAGT
 CAGCCCAGACATGTGCATCTACATCACAGAGGACATGCTCATGTGCGGGAAGCTGAATGGACACTCTGGG
 TTGATTGTGAAAGAAATTGGGTCTTCCACCTCGAGCTCTTCAGAAACAGTTGTTAAGCTTCGTGGCCAGA
 GTACTGATTCTTCCACAGACTATATGTCGAAACCAAGACCTCCACTGATCGACACAGCTTGAGCCT
 CGATGACATCAGACTTTACCAGAAAGACTTCTGCGCATTGCAGGTCTGTGTCAGGACACTGCTCAGAGT
 TACACCTTTGGATGTGGCCATGAACTGGATGAGGAAGGCCTCTATTGCAACAGTTGCTTGGCCCAGCAGT
 GCATCAACATCCAAGATGCTTTCCAGTCAAAAGAACCAGCAAATACTTTCTCTGGATCTCACTCATGA
 TGAAGTCCAGAGTTTGTGTG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)


Protein Sequence: >RC233587 representing NM_001267047
 Red=Cloning site Green=Tags(s)

MDQLEKRSRAGSSAGSMKHKRLSRHSTASHSSSHTSGIEADTKPRDTGPEDSYSSSAIHRKLKTCSSMT
 SHGSSHTSGVESGGKDRLEEDLQDDEIEMLVDDPRDLEQMNEESLEVSPDMCIYITEDMLMSRKLNGHSG
 LIVKEIGSSTSSSSETVVKLRGQSTDLPQTICKPKTSTDHSLSLDDIRLYQKDFLRIAGLCQDTAQS
 YTFGCGHELDEEGLYCNSCLAQQCINIQDAFPVKRTSKYFSLDLTHDEVPEFVV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001267047

ORF Size: 792 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001267047.1, NP_001253976.1

RefSeq Size: 3623 bp

RefSeq ORF: 795 bp


Locus ID: 122786

UniProt ID: Q96NE9

Cytogenetics: 14q22.1

MW: 29.7 kDa

Product images:



Circular map for RC233587