

Product datasheet for **RC233592**

CCER2 (NM_001243212) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCER2 (NM_001243212) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CCER2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233592 representing NM_001243212 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCGCCCCGAGGGCCAGCCTCTGAGCTGCTGCTGCTGCGGCTGCTCCTGCTGGGGCGGCCACCGCTG
CTCCCTTGGCACCAGACCCCTCCAAGGAGGAGCTGACCCGCTGTCTGGCAGAGGTGGTCACAGAGGTGCT
GACCGTGGGCCAGGTCCAGAGAGGACCCCTGCACTGCTCTTCTCCACAAGGAGTTGTGCGGGACAGAGCCC
CACGGCTGTGCGTCCACCGAGGAGAAAGCCTGCTGCTTGGGGATTTCAGAAGCAGGAGGCTGGGAAGA
TGAGGTCCAGCCAGGAGGTGAGGGATGAGGAAGAGGAGGAGGTAGCAGAGAGGACCCACAAGTCTGAGGT
CCAGGAACAAGCCATCCGCATGCAAGGGCATCGCCAGCTCCACCAGGAGGAGGACGAGGAGGAGGAGAAG
GAGGAGAGGAAGAGGGGGCCCATGGAGACCTTTGAGGACCTGTGGCAGCGGCATCTAGAGAATGGAGGGG
ACCTCCAGAAGCGGGTGGCAGAGAAGGCCAGTGACAAAGAGACGGCCAGTTCCAGGCAGAGGAGAAGGG
GGTGCGGGTGTGGGCGGGGACCGCAGCCTGTGGCAGGGGGCCGAGAGAGGCGGAGGAGAGAGGCGCGAG
GACTTGCCCCACCACCACCACCACCAGCCAGAGGCTGAGCCCAGGCAGGAGAAGGAGGAGGCTT
CGGAGAGGGAGGAAAAGGAGGTGGAGCAGCTGGAGCACTTGAGAGACGAAGTGAAGAAGGTGACAGAGAC
GCTGGGGGAGCAGCTCAGGAGGGAGGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RC233592 representing NM_001243212
 Red=Cloning site Green=Tags(s)

MPPRGPASELLLLRLLLGAATAAPLAPRPSKEELTRCLAEVVTEVLTVGQVQRGPCTALLHKELCGTEP
 HGCASTEELGLLLDGDFKKQEAGKMRSSQEVREDEEEVEAERTHKSEVQEQAIRMQGHRQLHQEEDDEEEK
 EERKRGPMTFEDLWQRHLENGGDLQKRVAEKASDKETAQFQAEKGVRLGGDRSLWQGAERGGGERRE
 DLPHHHHHHHQPEAEPRQEKEEASEREKEVEQLEHLRDELKKVTETLGEQLRREG

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM 001243212

ORF Size: 798 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001243212.1, NP_001230141.1

RefSeq Size: 1074 bp

RefSeq ORF: 801 bp


Locus ID: 643669

UniProt ID: I3L3R5

Cytogenetics: 19q13.2

MW: 30.8 kDa

Product images:



Circular map for RC233592