

Product datasheet for **RC233606**

ZNF501 (NM_001258280) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	ZNF501
Synonyms:	ZNF; ZNF52
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC233606 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATCCAGCCAAATATCACTCAGAATGAAACATGGGAGAGTTAACATGCAGAAGAAACCTTCAAAGT
GTAGTGAATGTGGGAAGTTCTTTACTCAGAGATCATCTTTACCCAGCACCAGAGGATTCACAGAGGAGA
GAAGCCCTATGTGTGCAGTGAATGTGGAAGTTGTTCCGTAACAGTCAAATCTTACTCAACATCTGAGA
ATTCATACCGGAGAGAAACCTTATAATGTAATGAATGTGAGAAAGCCTTCAAACAAAAGCAATCTTG
TTCAGCATCTGAGAATTCATACTGGAGAGAAACCTATAAATGCAATGAATGTGGAAAAGCCTTTTGTC
GAGCCCATCCCTTATTAACACCAGCGAATTCATACTGGAGAGAAACCATATAAATGTACAGAATGTGGC
AAAGCCTTCAGTCAGAGCATATGCCTTACTCGTCATCAGAGAAGTCATTCTGGAGATAAACCTTTTAAGT
GTAATGAATGTGGAAAAGCCTTAAATCAGAGTGCATGTCTCATGCAGCATCAGAGAATTCATTGAGGAG
GAAGCCCTACACATGCACTGAATGTGTAAGCCTTCACTCAGAACTCTTCCCTTGTGAACATGAAAGG
ACTCAGACTGGAGAGAACTTTATAAGTGTAGTGAGTGTGAAAAAAGCTTCCGCAAAACAAGCACACCTTA
GTGAGCATTACAGAATTCATACTGGAGAAAAACCTTATGAGTGTGTTGGATGTGGGAAATCCTTTAGGCA
CAGTTCAGCACTTCTTCGACATCAGAGGCTTCATGCTGGAGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>RC233606 protein sequence
 Red=Cloning site Green=Tags(s)

MNSSQISLRMKHGRVNMQKKPSKCSECGKFFTRSSLTQHQRIRHGEKPYVCSECGSFRKQSNLTQHLR
 IHTGEKPYKCNECEKAFQTKAILVQHLRIHTGEKPYKCNECGKAFQSPSLIKHQRHTGEKPYKCTECG
 KAFSQSICLTRHQRSHSGDKPFKCNCEGKAFNQSAQLMQHQRHSGEKPYTCTECGKAFQTNSSLVEHER
 THTGEKLYKCSECEKTRFRQAHLSHYRIHTGEKPYECVCGCKSFRHSSALLRHQRLHAGE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk6325_g03.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001258280

ORF Size:

813 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001258280.2](#)

RefSeq Size: 3400 bp

RefSeq ORF: 816 bp

Locus ID: 115560

UniProt ID: [Q96CX3](#)


Cytogenetics: 3p21.31

Protein Families: Transcription Factors

MW: 31.2 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC233606