

## Product datasheet for RC233765

### G protein beta subunit like (MLST8) (NM\_001199175) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	G protein beta subunit like (MLST8) (NM_001199175) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	G protein beta subunit like
Synonyms:	GbetaL; GBL; LST8; POP3; WAT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233765 representing NM_001199175 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAACACCTCCCAGGCACGGTGGGCAGTGACCCGGTCATCCTGGCCACTGCAGGCTACGACCACCCG  
TGCCTTCTGGCAGGCCACAGCGGCATCTGCACCCGGACGGTGCAGCACCAGGACTCCGTGAATGCCTT  
GGAGGTACACCCGACCGCAGCATGATTGCTGCTGCAGGTTACCAGCACATCCGCATGTATGATCTCAAC  
TCCAATAACCCTAACCCATCATCAGCTACGACGGCGTCAACAAGAACATCGCGTCTGTGGGCTTCCACG  
AAGACGGCCGCTGGATGTACACGGCGGGGAGGACTGCACAGCCAGGATCTGGGACCTCAGTCCCAGGAA  
CCTGCAGTGCCAGCGGATCTCCAGGTGAACGCACCCATTAAGTGCCTGCTGCACCCCAACCAGGCA  
GAGCTCATCGTGGGTGACCAGAGCGGGGTATCCACATCTGGGACTTGAAAACAGACCACAACGAGCAGC  
TGATCCCTGAGCCCGAGGTCTCCATCACGTCCGCCACATCGATCCCGACGCCAGCTACATGGCAGCTGT  
CAATAGCACCCGAAACTGCTATGTCTGGAATCTGACGGGGGGCATTGGTGACGAGGTGACCCAGCTCATC  
CCCAAGACTAAGATCCCTGCCACACGCGCTACGCCCTGCAGTGTGCTTTCAGCCCCGACTCCACGCTCC  
TCGCCACCTGCTCGGCTGATCAGACGTGCAAGATCTGGAGGACGTCCAATTCTCCCTGATGACGGAGCT  
GAGCATCAAGAGCGGCAACCCCGGGGAGTCTCCCGCGGCTGGATGTGGGGCTGCCCTTCTCGGGGAC  
TCCCAGTACATCGTCACTGCTTCTCGGACAACCTGGCCCGGCTCTGGTGTGGAGACTGGAGAGATCA  
AGAGAGATATGGCGGCCACCAGAAGGCTGTTGTCTGCCTGGCCTCAATGACAGTGTGCTGGCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC233765 representing NM\_001199175  
 Red=Cloning site Green=Tags(s)

MNTSPGTVGSDPVILATAGYDHTVRFWQAHSIGICRTRVQHQDSVNALEVTPDRSMIAAAGYQHIRMVYDLN  
 SNNPNPIISYDGVNKNIASVGFHEDGRWMTGGEDCTARIWDLRSRNLQCQRIFQVNAPINCVCVCLHPNQ  
 ELIVGDQSGAIHIWDLKTDHNEQLIPEPEVSITSAHIDPDASYMAAVNSTGNCYVWNLTTGGIGDEVTQLI  
 PKTKIPAHTRYALQCRFSPDSTLLATCSADQTCKIWRTSNFSLMTEL SIKSGNPGESSRGMWGCASFSGD  
 SQYIVTASSDNLARLWCVETGEIKREYGGHQKAVVCLAFNDSVLG

TRTRPLEQKLISEEDLAANDILDYKDDDDKLV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_001199175

**ORF Size:** 975 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001199175.3</a>
<b>RefSeq Size:</b>	1928 bp
<b>RefSeq ORF:</b>	978 bp
<b>Locus ID:</b>	64223
<b>UniProt ID:</b>	<a href="#">Q9BVC4</a>
<b>Cytogenetics:</b>	16p13.3
<b>Protein Pathways:</b>	mTOR signaling pathway
<b>MW:</b>	36.2 kDa
<b>Gene Summary:</b>	<p>Subunit of both mTORC1 and mTORC2, which regulates cell growth and survival in response to nutrient and hormonal signals. mTORC1 is activated in response to growth factors or amino acids. Growth factor-stimulated mTORC1 activation involves a AKT1-mediated phosphorylation of TSC1-TSC2, which leads to the activation of the RHEB GTPase that potentially activates the protein kinase activity of mTORC1. Amino acid-signaling to mTORC1 requires its relocalization to the lysosomes mediated by the Ragulator complex and the Rag GTPases. Activated mTORC1 up-regulates protein synthesis by phosphorylating key regulators of mRNA translation and ribosome synthesis. mTORC1 phosphorylates EIF4EBP1 and releases it from inhibiting the elongation initiation factor 4E (eIF4E). mTORC1 phosphorylates and activates S6K1 at 'Thr-389', which then promotes protein synthesis by phosphorylating PDCD4 and targeting it for degradation. Within mTORC1, LST8 interacts directly with MTOR and enhances its kinase activity. In nutrient-poor conditions, stabilizes the MTOR-RPTOR interaction and favors RPTOR-mediated inhibition of MTOR activity. mTORC2 is also activated by growth factors, but seems to be nutrient-insensitive. mTORC2 seems to function upstream of Rho GTPases to regulate the actin cytoskeleton, probably by activating one or more Rho-type guanine nucleotide exchange factors. mTORC2 promotes the serum-induced formation of stress-fibers or F-actin. mTORC2 plays a critical role in AKT1 'Ser-473' phosphorylation, which may facilitate the phosphorylation of the activation loop of AKT1 on 'Thr-308' by PDK1 which is a prerequisite for full activation. mTORC2 regulates the phosphorylation of SGK1 at 'Ser-422'. mTORC2 also modulates the phosphorylation of PRKCA on 'Ser-657'. [UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for RC233765