

Product datasheet for **RC233800**

Dysbindin (DTNBP1) (NM_001271668) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dysbindin (DTNBP1) (NM_001271668) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dysbindin
Synonyms:	BLOC1S8; DBND; HPS7; My031; SDY
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233800 representing NM_001271668 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGGAGACCCTTCGCGAGCGGCTGCTGAGCGTGCAGCAGGATTTACCTCCGGGCTGAAGACTTTAA
GTGACAAGTCAAGAGAAGCAAAGTAAAAGCAAACCCAGGTATGAGGATACATGGGCTGCACCTTCACAG
AAGAGCCAAAGACTGTGCAAGTGTGGAGAGCTGGTGGATAGCGAGGTGGTCATGCTTTCTGCGCACTGG
GAGAAGAAAAAGACAAGCCTCGTGGAGCTGCAAGAGCAGCTCCAGCAGCTCCAGCTTTAATCGCAGACT
TAGAATCCATGACAGCAAATCTGACTCATTTAGAGGCGAGTTTTGAGGAGGTAGAGAACAACCTGCTGCA
TCTGGAAGACTTATGTGGCAGTGTGAATTAGAAAAGATGCAAACATATGCAGTCCCAGCAACTGGAGAAT
TACAAGAAAAATAGAGGAAGGAACCTTCAAAGCTGAACTAGATGCAGAGCACGCCCAGAAAGG
TCCTGGAATGGAGCACACCCAGCAAATGAAGCTGAAGGAGCGGCAGAAAGTTTTTGGAGGAAGCCTTCCA
GCAGGACATGGAGCAGTACCTGTCCACTGGCTACCTGCAGATTGCAGAGCGGCGAGAGCCATAGGCAGC
ATGTCATCCATGGAAGTGAACGTGGACATGCTGGAGCAGATGGACCTGATGGACATATCGGACCAGGAGG
CCCTGGACGTCTTCTGAACTCTGGAGGAGAAGAGAACACTGTGCTGTCCCCGCCTTAGGGCTGAATC
CAGTACCTGTGAGAATGAGATTACCCTCCAGGTTCCAAATCCCTCAGAATTAAGAGCCAAGCCACCTTCT
TCTTCTCCACCTGCACCGACTCGGCCACCCGGGACATCAGTGAGGGTGGGGAGTCCCCGTTGTTTCAGT
CCGATGAGGAGGAAGTTCAGGTGGACACTGCCCTGGCCACATCACACACTGACAGAGAGGCCACTCCGGA
TGGTGGTGGAGACAGCGACTCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC233800 representing NM_001271668
 Red=Cloning site Green=Tags(s)


MLETLRERLLSVQQDFTSGLKTLSDKSREAKVKSKPRYEDTWAALHRRAKDCASAGELVDSEVVMLSAHW
 EKKKTSLVELQEQLQLPALIADLESMTANLTHLEASFEEVENNLLHLEDLCGQCELERCKHMQSQQLEN
 YKKNKRKELETFKAELDAEHAQKVLMEHTQMKLKERQKFFEEAFQDMEQYLSTGYLQIAERREPIGS
 MSSMEVNVMDLEQMDLMDISDQEALDVFNLNSGGEENTVLSPALGPESSTCQNEITLQVPNPSELRAKPPS
 SSSTCTDSATRDISEGGESPVVQSDDEEVQVDTALATSHTDREATPDGGEDSDS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001271668

ORF Size: 1002 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271668.1](#), [NP_001258597.1](#)

RefSeq Size: 1373 bp

RefSeq ORF: 1005 bp

Locus ID: 84062

UniProt ID: [Q96EV8](#)


Cytogenetics: 6p22.3

Protein Families: Druggable Genome

MW: 38.1 kDa

Gene Summary: This gene encodes a protein that may play a role in organelle biogenesis associated with melanosomes, platelet dense granules, and lysosomes. A similar protein in mouse is a component of a protein complex termed biogenesis of lysosome-related organelles complex 1 (BLOC-1), and binds to alpha- and beta-dystrobrevins, which are components of the dystrophin-associated protein complex (DPC). Mutations in this gene are associated with Hermansky-Pudlak syndrome type 7. This gene may also be associated with schizophrenia. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC233800