

Product datasheet for RC233816

SPOCK3 (NM_001204355) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPOCK3 (NM_001204355) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SPOCK3
Synonyms:	HSAJ1454; TES-3; TICN3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233816 representing NM_001204355 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATAACAATGGCTCACCAATCTCTCAGTATGACAAGGAAGTCGGACAGTGGAAACAAATCCGAGA
CGATGAAAGAAGCAGGAGTAGACCATAGGCAGTGGAGGGTCCCATTATCCACCTGCAAGCAGTGCCC
AGTGGTCTATCCCAGCCCTGTTTGGTTCAGATGGTCATACCTACTCTTTTCAGTGCAAACTAGAATAT
CAGGCATGTGTCTTAGGAAAACAGATCTCAGTCAAATGTGAAGGACATTGCCCATGTCCTTCAGATAAGC
CCACCAGTACAAGCAGAAATGTTAAGAGAGCATGCAGTGACCTGGAGTTCAGGGAAGTGGCAAACAGATT
GCGGGACTGGTTCAAGGCCCTTCATGAAAGTGGAAAGTCAAACAAGAAGCAAAAACATTGCTGAGGCC
GAGAGAAGCAGATTTCGATACCAGCATCTTGCCAAATTTGCAAGGACTCACTTGGCTGGATGTTTAAACAGAC
TTGATACAAACTATGACCTGCTATTGGACCAGTCAGAGCTCAGAAGCATTACCTTGATAAGAATGAACA
GTGTACCAAGGCATTCTCAATTCTGTGACACATACAAGGACAGTTAATATCTAATAATGAGTGGTGC
TACTGCTCCAGAGACAGCAAGACCCACCTTGCCAGACTGAGCTCAGCAATTCAGAAGCGCAAGGGG
TAAAGAAGCTCCTAGGACAGTATATCCCCTGTGTGATGAAGATGGTTACTACAAGCCAACAATGTCA
TGGCAGTGTGGACAGTGCTGGTGTGTGACAGATATGGAATGAAGTCATGGGATCCAGAATAAATGGT
GTTGCAGATTGTCTATAGATTTTGGATCTCCGGAGATTTTGCTAGTGGCGATTTTCATGAATGGACTG
ATGATGAGGATGATGAAGACGATATTATGAATGATGAAGATGAAATTGAAGATGATGATGAAGATGAAGG
GGATGATGATGATGGTGGTATGACCATGATGTATACATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC233816 representing NM_001204355
Red=Cloning site Green=Tags(s)

MINNGSPQSLSMTRKSDSGTNSMKEAGVDHRQWRGPILSTCKQCPVVYSPVCGSDGHTYSFQCKLEY
 QACVLGKQISVKCEGHPCPSDKPTSTSRNVKRACSDLEFREVANRLRDWFKALHESGSQNKTKTLLRP
 ERSRFDTSILPICKDSLWGMFNRLDTNYDLLLDQSELRSIYLDKNEQCTKAFFNSCDTYKDSLISNNEW
 YCFQRQDPPCQTEL SNIQKRQGVKLLGQYIPLCDEDEGYKPTQCHGSVGCWCVDRYGNEVMGSRING
 VADCAIDFEISGDFASGDFHEWTDDEDEDDIMNDEDEIEDDEDEGDDDDGGDDHDVYI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001204355

ORF Size: 1020 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204355.1](#), [NP_001191284.1](#)

RefSeq Size: 2819 bp

RefSeq ORF: 1023 bp

Locus ID: 50859

UniProt ID: [Q9BQ16](#)


Cytogenetics: 4q32.3

Protein Families: Secreted Protein, Transmembrane

MW: 39.1 kDa

Gene Summary: This gene encodes a member of a novel family of calcium-binding proteoglycan proteins that contain thyroglobulin type-1 and Kazal-like domains. The encoded protein and may play a role in adult T-cell leukemia by inhibiting the activity of membrane-type matrix metalloproteinases. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Oct 2011]

Product images:



Circular map for RC233816