

## Product datasheet for RC233843

### TMEM185B (NM\_024121) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM185B (NM_024121) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TMEM185B
Synonyms:	FAM11B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC233843 representing NM_024121 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**

ATGAACCCAGGGCCTGTTCCAGGACTTCAACCCAGTAAGTTTCTCATCTACACCTGCCTGCTGCTCT  
TCTCGGTGCTGCTGCCCCCGCCTGGACGGCATCATCCAATGGAGCTACTGGCCGCTTTGCCCCAT  
ATGGCTGTGGAAGCTTCTAGTCGTCGACGGCGCCTCCGTGGCGCGGGCCTTTGGCCCGCAACCCCTCGC  
TACCGCACCGAGGAGAGGCCTGTGTGGAGTTCAAAGCCATGCTGATCGCTGTGGGCATCCACCTGCTGC  
TGCTCATGTTCGAAGTCTGGTCTGCGACAGGGTGGAGAGGGGCACCCACTTCTGGCTGCTGGTCTTCAT  
GCCTCTCTTCTCGTGTCCCCGTGTCCTGGCTGCCTGCGTCTGGGGCTTTCGACACGATAGGTCGCTG  
GAGCTGGAGATCCTGTGCTCGGTCAACATCCTGCAGTTTCATCTTCATCGCCCTAAAGCTGGACAGGATTA  
TTCAGTGGCCGTGGCTGGTGGTGTGTTGTCGCCCTGTGGATCCTCATGTCGTTCTTTGCCTGGTCTGCT  
CTATTACATCGTCTGGTCCCTCCTGTTCCCTGCGGTCCCTGGATGTGGTTGCCGAGCAGCGGAGAACACAC  
GTGACCATGGCTATCAGTTGGATAACGATTGTCGTCGCTCTGCTCACTTTTGAGGTCCTGCTGGTTCACA  
GATTGGATGGCCACAATACATTCTCCTACGTCTCCATATTTGCCCCCTTTGGCTTTCCTACTAAGTTT  
AATGGCCACAACATTTAGGCGAAAGGGGGCAATCATTGGTGGTTGGCATTCCGAGAGACTTCTGTCAG  
TTCTGCTTGAATTTTCCATTTTAAAGAGAATATGGGAACATTCATATGATCCTCCATCACAAGATA  
GTGAAGATGCTGAAGAAACATCAGTTCAGAAAGCTCCGAAAATTGCTCCAATATTTGGAAGAAGGCCAG  
AGTAGTTATAACCCAGAGCCCTGGGAAATACGTTCCCCCCCCTCCAAGTTAAATATTGATATGCCAGAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_024121.3](#)

**RefSeq Size:** 5938 bp

**RefSeq ORF:** 1053 bp


**Locus ID:** 79134

**UniProt ID:** [Q9H7F4](#)

**Cytogenetics:** 2q14.2

**MW:** 41 kDa

**Product images:**



Circular map for RC233843