

Product datasheet for **RC234119**

GPR172A (SLC52A2) (NM_001253815) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GPR172A (SLC52A2) (NM_001253815) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	GPR172A
Synonyms:	BVLS2; D15Ert747e; GPCR41; GPR172A; hRFT3; PAR1; RFT3; RFVT2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC234119 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCAGCACCCACGCCCGCCGTCGGTGTGACCCACCTGCTGGTGGCTCTCTTCGGCATGGGCTCCT
GGGCTGCGGTCAATGGGATCTGGGTGGAGCTACCTGTGGTGGTCAAAGAGCTTCCAGAGGGTTGGAGCCT
CCCCTTTACGTCTCTGTGCTTGTGGCTCTGGGAACCTGGGTCTGCTGGTGGTACCCTCTGGAGGAGG
CTGGCCCCAGAAAGGACGAGCAGGTCCCCATCCGGTGGTGCAGGTGCTGGGCATGGTGGGCACAGCCC
TGCTGGCCTCTCTGTGGCACCATGTGGCCCCAGTGGCAGGACAGTTGCATTCTGTGGCCTTCTTAGCACT
GGCCTTTGTGCTGGCACTGGCATGCTGTGCCTCGAATGCACTTTCTGCCCTTCTTGAGCCACCTGCCA
CCTCGTCTTACGGTATTCTTCTGGGTCAAGGCCTGAGTGCCTGCTGCCCTGCGTGTGGCCCTAG
TGCAGGGTGTGGCCGCCTCGAGTGCCCGCCAGCCCCATCAACGGCACCCCTGGCCCCCGCTCGACTT
CCTTGAGCGTTTTCCCGCCAGCACCTTCTTCTGGGCACTGACTGCCCTTCTGGTGCCTTCAGCTGTGCC
TTCCAGGGTCTTCTGCTGCTGTTGCCGCCACCACCATCTGTACCCACAGGGGAGTTAGGATCAGGCCTCC
AGGTGGGAGCCCCAGGAGCAGAGGAAGAGTGGAGAGTCTCACCACTGCAAGAGCCACCAAGCCAGGC
AGCAGGCACCACCCTGGTCCAGACCCTAAGGCCTATCAGCTTCTATCAGCCCGCAGTGCCTGCCTGCTG
GGCCTGTTGGCCGCCACCAACGCGCTGACCAATGGCGTGTGCCTGCCGTGCAGAGCTTTTCTGCTTAC
CCTACGGGCGTCTGGCCTACCACCTGGCTGTGGTGTGGGCGAGTGTGCCAATCCCCTGGCCTGCTTCT
GGCCATGGGTGTGCTGTGCAGGTCTTGGCAGGGCTGGGCGGCCCTCTCTGCTGGGCGTGTCTGTGGG
GGCTACCTGATGGCGCTGGCAGTCTGAGCCCCTGCCCGCCCTGGTGGGCACCTCGGCGGGGTGGTCC
TCGTGGTGTGCTGTGGGTGCTGTGCTTGGCGTGTCTCCTACGTGAAGGTGGCAGCCAGCTCCCTGCT
GCATGGCGGGGGCCGCGGCATTGCTGGCAGCCGCGTGGCCATCCAGGTGGGCTCTCTGCTCGGCGCT
GTTGCTATGTTCCCCCGACCAGCATCTATCACGTGTTCCACAGCAGAAAGGACTGTGCAGACCCTGTG
ACTCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC234119 protein sequence
Red=Cloning site Green=Tags(s)

MAAPTPARPVLTHLLVALFGMGSWAAVNGIWWELPVVVKELPEGWSLPSYVSVLVALGNLGLLVVTLWRR
LAPGKDEQVPIRVVQVLGMVGTALLASLWHHVAVVAGQLHSVAFALAFVLALACCASNVTFLPFLSHLP
PRFLRSFFLQGLSALLPCVLALVQGVGRLECPPAPINGTPGPPLDFLERFPASTFFWALTALLVASAAA
FQGLLLLLPPPSVPTGELGSLQVQVGAEEVEESSPLQEPPSQAAGTTPGDPKAYQLLSARSACLL
GLLAATNALTNGVLPVQSFSLPYGRLAYHLAVVLGSAANPLACFLAMGVLCRSLAGLGLSLLGVFCG
GYLMALAVLSPCPPLVGTSAVVVLSWVLCVGFVSYVVAASSLLHGGRPALLAAGVAIQVGSLLGA
VAMFPPTSIYHVFHSRKDCADPCDS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6085_d12.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_001253815

ORF Size: 1335 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001253815.2](#)

RefSeq Size: 1806 bp

RefSeq ORF: 1338 bp


Locus ID: 79581

UniProt ID: [Q9HAB3](#)

Cytogenetics: 8q24.3

Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	45.8 kDa
Gene Summary:	This gene encodes a membrane protein which belongs to the riboflavin transporter family. In humans, riboflavin must be obtained by intestinal absorption because it cannot be synthesized by the body. The water-soluble vitamin riboflavin is processed to the coenzymes flavin mononucleotide (FMN) and flavin adenine dinucleotide (FAD) which then act as intermediaries in many cellular metabolic reactions. Paralogous members of the riboflavin transporter gene family are located on chromosomes 17 and 20. Unlike other members of this family, this gene has higher expression in brain tissue than small intestine. Alternative splicing of this gene results in multiple transcript variants encoding the same protein. Mutations in this gene have been associated with Brown-Vialetto-Van Laere syndrome 2 - an autosomal recessive progressive neurologic disorder characterized by deafness, bulbar dysfunction, and axial and limb hypotonia. [provided by RefSeq, Jul 2012]

Product images:



Circular map for RC234119