

Product datasheet for **RC234127**

Glucocorticoid Receptor (NR3C1) (NM_001204263) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Glucocorticoid Receptor (NR3C1) (NM_001204263) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Glucocorticoid Receptor
Synonyms:	GCCR; GCR; GCRST; GR; GRL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234127 representing NM_001204263
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTACCACTATGACATGAATACAGCATCCCTTTCTCAACAGCAGGATCAGAAGCCTATTTTAAATGTCA
 TTCACCAATTCCCGTTGGTTCGAAAATTGGAATAGGTGCCAAGGATCTGGAGATGACAACTTGACTTC
 TCTGGGACTCTGAACTCCCTGGTCAACAGTTTTTTCTAATGGCTATTCAAGCCCAGCATGAGACCA
 GATGTAAGCTCTCCTCCATCCAGCTCTCAACAGCAACAACAGGACCACCTCCAAACTCTGCCTGGTGT
 GCTCTGATGAAGCTTCAGGATGTCATTATGGAGTCTTAACTTGTGGAAGCTGTAAAGTTTTCTCAAAG
 AGCAGTGAAGGACAGCACAATTACCTATGTGCTGGAAGGAATGATTGCATCATCGATAAAATTCGAAGA
 AAAAACTGCCAGCATGCCGCTATCGAAAATGTCTTCAGGCTGGAATGAACCTGGAAGCTCGAAAAACA
 AAAAAAATAAAAGGAATTCAGCAGGCCACTACAGGAGTCTACAAGAACTCTGAAAATCCTGGTAA
 CAAAACAATAGTTCCTGCAACGTTACCACAACCTACCCCTACCTGGTGTACTGTTGGAGTTATTGAA
 CCTGAAGTGTATATGCAGGATATGATAGCTCTGTTCCAGACTCAACTGGAGGATCATGACTACGCTCA
 ACATGTTAGGAGGGCGCAAGTATTGCAGCAGTAAATGGGCAAGGCAATACCAGTTTCAGGAACTT
 ACACCTGGATGACCAATGACCCACTGTCAGTACTCCTGGATGTTTCTTATGGCATTGCTCTGGGGTGG
 AGATCATATAGACAATCAAGTGCAAACTGCTGTGTTTTGCTCCTGATCTGATTATTAATGAGCAGAGAA
 TGACTCTACCCTGCATGTACGACCAATGTAACACATGCTGTATGTTTCTCTGAGTTACACAGGCTTCA
 GGTATCTTATGAAGATATCTCTGTATGAAAACCTTACTGCTTCTCTCTCAGTTCTTAAGGACGGTCTG
 AAGAGCCAAGAGCTATTTGATGAAATAGAATGACCTACATCAAAGAGCTAGGAAAAGCCATTGTCAAGA
 GGAAGGAACTCCAGCCAGAAGTGGCAGCGGTTTTATCAACTGACAAAACCTTTGGATTCTATGCATGA
 AGTGGTTGAAAATCTCCTTAACTATTGCTTCAAACATTTTTGGATAAGACCATGAGATTGAATCCCC
 GAGATGTTAGCTGAAATCATACCAATCAGATACCAAAATATTCAAATGGAAATATCAAAAACTTCTGT
 TTCATCAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC234127 representing NM_001204263
 Red=Cloning site Green=Tags(s)

MYHYDMNTASLSQQDQKPIFNVIPPIVGSSENWNRQSGDDNLTSLGTLNFPGRVTFVSNYSSPSMRP
 DVSSPPSSSTATTGPPPKLCLVCSDEASGCHYGVLTCGSKVFFKRAVEGQHNYLCAGRNDICIIDKIRR
 KNCPACRYRKCLQAGMNLARKTKKKIKIGIQATTGVSQETSENPKNKTIVPATLPQLTPTLVSLLEVIE
 PEVLAYGYDSSVPDSTWRIMTTLNMLGGRQVIAAVKWAKAIPGFRNLHLDQMTLLQYSWMFLMAFALGW
 RSYRQSSANLLCFAPDLIINEQRMTLPCMYDQCKHMLYVSELHRLQVSYEYLCMKTLSSVPKDGL
 KSQELFDEIRMTYIKELGKAIKREGNSSQNWQRFYQLTKLLDSMHEVVENLLNYCFQFLDKTMSIEFP
 EMLAEIITNQIPKYSNGNIKKLLFHQK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001204263

ORF Size: 1341 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204263.2](#)

RefSeq Size: 6801 bp


RefSeq ORF: 1344 bp

Locus ID: 2908

UniProt ID: [P04150](#)

Cytogenetics:	5q31.3
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Protein Pathways:	Neuroactive ligand-receptor interaction
MW:	51 kDa
Gene Summary:	This gene encodes glucocorticoid receptor, which can function both as a transcription factor that binds to glucocorticoid response elements in the promoters of glucocorticoid responsive genes to activate their transcription, and as a regulator of other transcription factors. This receptor is typically found in the cytoplasm, but upon ligand binding, is transported into the nucleus. It is involved in inflammatory responses, cellular proliferation, and differentiation in target tissues. Mutations in this gene are associated with generalized glucocorticoid resistance. Alternative splicing of this gene results in transcript variants encoding either the same or different isoforms. Additional isoforms resulting from the use of alternate in-frame translation initiation sites have also been described, and shown to be functional, displaying diverse cytoplasm-to-nucleus trafficking patterns and distinct transcriptional activities (PMID:15866175). [provided by RefSeq, Feb 2011]

Product images:



Circular map for RC234127