

Product datasheet for **RC234204**

RNF14 (NM_001201365) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RNF14 (NM_001201365) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	RNF14
Synonyms:	ARA54; HFB30; HRIHFB2038; TRIAD2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234204 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

RCATGTCGTCAGAAGATCGAGAAGCTCAGGAGGATGAATTGCTGGCCCTGGCAAGTATTTACGATGGAGA
 TGAATTTAGAAAAGCAGAGTCTGTCCAAGTGGAGAAACCAGGATCTATTTGGATTTGCCACAGAATTTCC
 AAGATATTTGTGAGCGGCAATCAAAATGAGTGTCTCCAGAATAGTGGCTTTGAATACACCATTTGCTTTT
 TGCTCCACTTGTGCTGAACTTTGAACTGCCACCAGATTATCCATCCTCTTCCCCACCTTCACTCACACT
 TAGTGGCAAATGGCTGTCACTCACTCAGCTATCTGCTCTATGCAAGCACTTAGACAACCTATGGGAAGAA
 CACCGTGGCAGCGTGGTCTGTTGCTGGATGCAATTTCTTAAGGAAGAGACCCTAGCATACTTGAATA
 TTGTCTCTCTTTGAGCTCAAGATTGGTTCTCAGAAAAAGTGCAGAGAAGGACAGCTCAAGCTTCTCC
 CAACACAGAGCTAGATTTGGAGGAGCTGCTGGATCTGATGTAGACCAAGAGGAAATTTGGATGAGAGA
 GCAGTGCAGGATGTGAATCACTGTCAAATCTGATCCAGGAAATCTTGGACTTTGATCAAGCTCAGCAGA
 TAAATGCTTTAATAGTAAATGTTCTGTGCAAGTATCTGTTTCTGTGAGAAGCTGGGTAGTGAATGCAT
 GTACTTCTGGAGTGCAGGCATGTGACTGCAAAGCCTGTCTGAAGGACTACTTTGAAATCCAGATCAGA
 GATGGCCAGGTTCAATGCCTCAACTGCCAGAACCAAAGTGCCCTTCGGTGGCCACTCCTGGTCAGGTCA
 AAGAGTTAGTGAAGCAGAGTTATTTGCCCGTTATGACCGCTTCTCCTCCAGTCTCCTTGGACCTGAT
 GGCAGATGTGGTACTGCCCGCGCTGCTGCCAGTGCCTGTGATGCAGGAACCTGGCTGCACCATG
 GGTATCTGCTCCAGTCAATTTGCTTCTGTACTTTGTGAGGTTGACCTACCATGGGGTCTCCCAT
 GTAAGGTGACTGCAGAGAAATTAATGGACTTACGAAATGAATACCTGCAAGCGGATGAGGCTAATAAAG
 ACTTTTGGATCAAAGGTATGGTAAGAGAGTGATTCAGAAGGCACCTGGAAGAGATGGAAGTAAGGAGTGG
 CTAGAGAAGAACTCAAAGAGCTGCCATGTTGTGGAATCCCATAGAGAAATTAGACGGATGTAACAAGA
 TGACATGTACTGGCTGTATGCAATATTTCTGTTGATTTGCATGGTTCTCTCTAGAGCAAACCTTA
 CAAACATTTCAATGACCTGTTTACCATGTTTTAACCGCTGTTTTATGCTGTGGATGTTGACGACGAT
 ATTTGGGAAGATGAGGTAGAAGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC234204 protein sequence
 Red=Cloning site Green=Tags(s)

XCRQKIEKLRMNCWPWQVFTMEMNLEKQSLSKVEKPGSIWICHIRISRYL*AAIQMSVSRIVALNTPFAF
 CLHLC*TLNCHQIIHPLPHLHSHLVANGCHQLSYLLYAST*TTYGKNTVAAWSCLPGCNFLRKRP*HT*I
 LSLLLSSRLVLRKKCREGQLKLLPTQS*ILEELLDLM*TKRKLWMREQCRMWNHCQI*SRKSWTLIKLSR
 *NALIVNCSCAVSVSVRSWVNACTSWSAGMCTAKPV*RTTLKSRSEMARFNASTAQNSALRWPLLVRS
 KS*WKQSYLPVMTAFSSPPWT*WQMWCTAPGRAASCL*CRNLAAPWVSAPAAILPSVLCAG*PTMGSPH
 VR*LQRN*WTYEMNTCKRMRLIKDFWIKGMVRE*FRRHWKRWKVRSG*RRTQRAAHVVELP*RN*TDVTR
 *HVLAVNCISVGF AFWLSLEQTL TNISMTLVHHVLTGCFMLWMLTTIFGKMR*K

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mg2876_e12.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001201365

ORF Size: 1422 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001201365.1](#), [NP_001188294.1](#)

RefSeq Size: 4181 bp

RefSeq ORF: 1425 bp

Locus ID: 9604

UniProt ID: [Q9UBS8](#)


Cytogenetics: 5q31.3

Protein Families: Druggable Genome, Transcription Factors

MW: 53.8 kDa

Gene Summary: The protein encoded by this gene contains a RING zinc finger, a motif known to be involved in protein-protein interactions. This protein interacts with androgen receptor (AR) and may function as a coactivator that induces AR target gene expression in prostate. A dominant negative mutant of this gene has been demonstrated to inhibit the AR-mediated growth of prostate cancer. This protein also interacts with class III ubiquitin-conjugating enzymes (E2s) and may act as a ubiquitin-ligase (E3) in the ubiquitination of certain nuclear proteins. Six alternatively spliced transcript variants encoding two distinct isoforms have been reported. [provided by RefSeq, Jan 2011]

Product images:



Circular map for RC234204