

Product datasheet for **RC234295**

p73 (TP73) (NM_001204184) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	p73 (TP73) (NM_001204184) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	p73
Synonyms:	P73
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234295 representing NM_001204184
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCCAGTCCACCGCCACCTCCCTGATGGGGCACCACGTTTGGACCTCTGGAGCTCTCTGGAAC
 CAGACAGCACCTACTTCGACCTTCCCAGTCAAGCCGGGGAATAATGAGGTGGTGGGCGGAACGGATTC
 CAGCATGGACGTCTTCCACCTGGAGGGCATGACTACATCTGTTCATGGCCAGTTCAATCTGCTGAGCAGC
 ACCATGGACCAGATGAGCAGCCGCGCGCTCGGCCAGCCCTACACCCAGAGCAGCCGCCAGCGTGC
 CCACCCACTCGCCCTACGCACAACCCAGCTCCACCTTCGACACCATGTGCGCCGGCGCTGTATCCCTC
 CAACACCGACTACCCCGACCCACCACCTTTGAGGTCACTTTCCAGCAGTCCAGCAGGCCAAGTCAGCC
 ACCTGGACGTACTCCCGCTCTGAAGAACTCTACTGCCAGATCGCAAGACATGCCCATCCAGATCA
 AGGTGTCCACCCGCCACCCAGGACCGCCATCCGGGCCATGCCTGTTTACAAGAAAGCGGAGCAGT
 GACCGACGTCTGAAACGCTGCCCAACCAGAGCTCGGGAGGGACTTCAACGAAGGACAGTCTGTCCA
 GCCAGCCACCTCATCCGGTGGAAAGCAATAATCTCTCGCAGTATGTGGATGACCCTGTACCCGGCAGGC
 AGAGCGTCTGGTGCCTATGAGCCACCACAGGTGGGGACGGAATTCACCACCATCCTGTACAATTCAT
 GTGTAACAGCAGCTGTGTAGGGGCATGAACCGCGGCCATCCTCATCATCACCCCTGGAGATGCGG
 GATGGGACAGTCTGGCCGCGGTCTTTGAGGGCCGATCTGCGCCTGTCTGGCCGCGACCGAAAAG
 CTGATGAGGACCACTACCGGAGCAGCAGGCCCTGAACGAGAGCTCCGCCAAGAAGCGGGCCGACGAA
 GCGTGCCTTAAGCAGAGCCCCCTGCCGTCCCGCCCTTGGTCCCGGTGTGAAGAAGCGGGCCGATGGA
 GACGAGGACACGTACTACCTCAGGTGCGAGGCCGGGAGAAGTTGAGATCCTGATGAAGCTGAAAGAGA
 GCCTGGAGCTGATGGAGTTGGTCCCGCAGCCACTGGTGGACTCCTATCGGCAGCAGCAGCTCCTACA
 GAGGCCGAGTCACTACAGCCCCCTCCTACGGGCCGGTCTCTCGCCATGAACAAGGTGCACGGGGGC
 ATGAACAAGTGCCTCCGTCAACCAGCTGGTGGCCAGCCTCCCCGCACAGTTCGGCAGCTACACCCA
 ACCTGGGGCCGTGGGCCCGGATGCTCAACAACCATGGCCACGCAGTCCAGCCAACGGCGAGATGAG
 CAGCAGCCACAGCGCCAGTCCATGGTCTCGGGTCCCACTGCACTCCGCCACCCCTACCACGCCGAC
 CCCAGCCTCGTCAGGACCTGGGGCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC234295 representing NM_001204184
 Red=Cloning site Green=Tags(s)

MAQSTATSPDGGTTFEHLWSSLEPDSTYFDLPQSSRGNNVVGTDSSMDVHLEGMTTSVMAQFNLLSS
 TMDQMSSRAASAPYTPPEHAASVPTHSPYAQPSSFTDTMSPAPVIPSNTDYPGPHHFEVTFQSSSTAKSA
 TWTYSPLLKLYCQIAKTCPIQIKVSTPPPGTAIRAMPVYKAEHVTDVVKRCPNHELGRDFNEGQSAP
 ASHLIRVEGNNLSQYVDDPVTGRQSVVVPYEPQVGFTEFTTILYNFMCNSSCVGGMNRRPILIIITLEM
 DGQVLRGRSFEGRICACPRDRKADEDHYREQALNESSAKNGAASKRAFQSPPAVPALGAGVKKRRHG
 DEDTYYLQVRGRENFEILMKLKESELELMELVPQPLVDSYRQQQLLQRPSTLQPPSYGPVLSMKNVHGG
 MNKLPVSNQLVGQPPPHSSAATPNLGPVGPMLNNHGHAVPANGEMSSSHSAQSMVSGSHCTPPPPYHAD
 PSLVRTWGP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001204184

ORF Size: 1497 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204184.2](#)

RefSeq Size: 5060 bp

RefSeq ORF: 1500 bp

Locus ID: 7161

UniProt ID: [O15350](#)

Cytogenetics: 1p36.32


Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Neurotrophin signaling pathway, p53 signaling pathway

MW: 54.8 kDa

Gene Summary: This gene encodes a member of the p53 family of transcription factors involved in cellular responses to stress and development. It maps to a region on chromosome 1p36 that is frequently deleted in neuroblastoma and other tumors, and thought to contain multiple tumor suppressor genes. The demonstration that this gene is monoallelically expressed (likely from the maternal allele), supports the notion that it is a candidate gene for neuroblastoma. Many transcript variants resulting from alternative splicing and/or use of alternate promoters have been found for this gene, but the biological validity and the full-length nature of some variants have not been determined. [provided by RefSeq, Feb 2011]

Product images:



Circular map for RC234295