

Product datasheet for **RC234296**

PUF60 (NM_001271100) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PUF60 (NM_001271100) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PUF60
Synonyms:	FIR; RoBPI; SIAHBP1; VRJS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234296 representing NM_001271100
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGAACGGGCAGAGCACAGCCGCCAAGCTGGGGCTGCCTCCCCTGACGCCCGAGCAGCAGGAGGCC
 TTCAGAAGGCCAAGAAGTACGCCATGGAGCAGAGCATCAAGAGTGTGCTGGTGAAGCAGACCATCGCGCA
 CCAGCAGCAGCAGCTACCAACCTGCAGATGGCGGCTCAGCGGCAGCGGGCGCTGGCCATCATGTGCCGC
 GTCTACGTGGGCTCTATCTACTATGAGCTGGGGGAGGACACCATCCGCCAGGCCTTTGCCCTTTGGCC
 CCATCAAGAGCATCGACATGTCTGGGACTCCGTACCATGAAGCACAAGGGCTTTGCCTTCGTGGAGTA
 TGAGGTCCCCGAAGCTGCACAGCTGGCCTTGGAGCAGATGAACTCGGTGATGCTGGGGGGCAGGAACATC
 AAGGTGGCAGACCCAGCAACATAGGGCAGGCCAGCCATCATAGACCAGTTGGCTGAGGAGGCACGGG
 CCTTCAACCGCATCTACGTGGCCTCTGTGCACCAGGACCTCTCAGACGATGACATCAAGAGCGTGTTTGA
 GGCTTTGGCAAGATCAAGTCTGCACACTGGCCCGGGACCCACAACCTGGCAAGCACAAGGGCTACGGC
 TTCATTGAGTACGAGAAGGCCAGTCGTCCCAAGATGCTGTGTCTCCATGAACCTCTTTGACCTGGGTG
 GCCAGTACTTGCGGGTGGCAAGGCTGTACACCGCCCATGCCCTACTCACACCAGCCACGCCTGGAGG
 CCTCCCACCTGCCGCTGCTGTGGCAGCTGCTGCAGCCACTGCCAAGATCACAGCTCAGGAAGCAGTGGCC
 GGAGCAGCGGTGCTGGTACCCTGGGCACACCTGGACTGGTGTCCCAGCACTGACCCTGGCCAGCCCC
 TGGGCACTTTGCCAGGCTGTATGGCTGCCAGGCACCTGGAGTATCACAGGTGTGACCCAGCCCCG
 TCCTCCTATCCGGTACCATCCCCTCGGTGGGAGTGGTGAACCCATCCTGGCCAGCCCTCAACGCTG
 GGTCTCCTGGAGCCCAAGAAGGAGAAGGAAGAAGAGGAGCTGTTCCCGAGTCAGAGCGCCAGAGATGC
 TGAGCGAGCAGGAGCAGATGAGCATCTCGGGCAGTAGCGCCCGACACATGGTGATGCAGAAGCTCTCCG
 CAAGCAGGAGTCTACAGTATGGTTCTGCGCAACATGGTGGACCCCAAGGACATCGATGATGACCTGGAA
 GGGGAGGTGACAGAGGAGTGTGGCAAGTTCGGGGCCGTGAACCGCGTCATCATCTACCAAGAAAACAAG
 GCGAGGAGGAGGATGCAGAAATCATTGTCAAGATCTTTGTGGAGTTTTCCATAGCCTCTGAGACTCATAA
 GGCCATCCAGGCCCTCAATGGCCGCTGGTTTGTGGCCGCAAGGTGGTGGCTGAAGTGTACGACCAGGAG
 CGTTTTGATAACAGTGACCTCTCTGCG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC234296 representing NM_001271100
 Red=Cloning site Green=Tags(s)

MENQSTAAKLG LPLTPEQEQEALQKAKKYAMEQSIKSVLVKQTI AHQQQQLTNLQMAAQRQRALAIMCR
 VYVGSIIYELGEDTIRQAFAPFGPIKSIDMSWDSVTMKHKGF AFVEYEVPEAAQLALEQMNSVMLGGRNI
 KVGGRPSNIGQAQPIIDQLAEERAFNR IYVASVHQDLSDDDIKSVFEAFGKIKSCTLRDPTTGKHKGYG
 FIEYEKAQSSQDAVSSMNLFDLGGQYLRV GKAVTPPMLLTPATPGGLPAAA AVAAAAATAKITAQEAVA
 GAAVLGTLGTPGLVSPALTLAQLGLTPQAVMAA QAPGVITGVTPARPPIPVTIPSVGVVNPILASPPTL
 LLELPKKEKEEEEELFPESERPEMLSEQEHMSISGSSARHMVMQKLLRKQESTVMVLRNMVDPKDIDDDLE
 GEVTEECGKFGAVNRVIIYQEKQGEEDAEIIVKIFVEFSIASETHKAIQALNGRWFAGRKVVAEYVDQE
 RFDNSDL SA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001271100

ORF Size: 1497 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271100.2](#)

RefSeq Size: 1916 bp

RefSeq ORF: 1500 bp

Locus ID: 22827

UniProt ID: [Q9UHX1](#)


Cytogenetics: 8q24.3

Protein Pathways: Spliceosome

MW: 54.5 kDa

Gene Summary: This gene encodes a nucleic acid-binding protein that plays a role in a variety of nuclear processes, including pre-mRNA splicing and transcriptional regulation. The encoded protein forms a complex with the far upstream DNA element (FUSE) and FUSE-binding protein at the myelocytomatosis oncogene (MYC) promoter. This complex represses MYC transcription through the core-TFIID basal transcription factor. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2012]

Product images:



Circular map for RC234296