

Product datasheet for **RC234476**

CFHR4 (NM_001201550) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CFHR4 (NM_001201550) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CFHR4
Synonyms:	CFHL4; FHR-4; FHR4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234476 representing NM_001201550
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTTGTACTAATCAATGTCATTCTGACCTTGTGGGTTTCTGTGCTAATGGACAAGAAGTGAACCTT
 GTGATTTTCCAGAAATCAACATGGAGGTCTATATTATAAGAGTTTTCGTAGACTATACTTTCCAGCAGC
 TGCAGGACAATCTTATTCTTACTGTGATCAAAATTTGTGACTCCTTCAGGAAGTTACTGGGATTAC
 ATTCATTGCACACAAGATGGTTGGTACCAACGGTCCCATGCCTCAGAACATGCTCAAAATCAGATGTAG
 AAATTGAAAATGGATTCAATTTCTGAATCTTCTCTATTATATTTTAAATGAAGAAACACAATATAATTG
 TAAACCAGGATATGCAACAGCAGAGGGAAATCTTCAGGATCAATTACATGTTTGCAAAATGGATGGTCA
 ACACAACCAATTTGCATTAATTTTGTGATATGCCTGTTTTTGAAGATCCAGAGCCAAGAGTAATGGCA
 TGTGGTTAAGCTCCATGACACATTGGACTATGAATGCTATGATGGATATGAAAGCAGTTATGGAACAC
 CACAGATTCATAGTGTGGTGAAGATGGCTGGTCCCATTTGCCAACATGCTATAATTCTTCAGAAAAC
 TGTGGGCTCCTCCACCTATTAGCAATGGAGATACCACGTCCTTCCCGCAAAAAGTGTATCTGCCATGGT
 CAAGAGTCGAGTACCAGTGCCAGTCTACTATGAACCTCAGGGTTCTAAATATGTAACATGTAGTAATGG
 AGACTGGTCAGAACCAAGATGCATATCAATGAAACCTTGTGAGTTTCCAGAAATCAACATGGACAT
 CTATATTATGAGAATACGCGTAGACCATACTTTCCAGTAGCTACAGGACAATCTTACTCTATTACTGTG
 ACCAAAATTTGTGACTCCTTCAGGAAGTTACTGGGATTACATTCAGTGCACACAAGATGGGTGGTTGCC
 AACAGTCCCATGCCTCAGAACATGCTCAAAATCAGATATAGAAATGAAAATGGATTCATTTCTGAATCT
 TCCTCTATTTATATTTTAAATAAAGAAATACAATATAAATGTAACCAGGATATGCAACAGCAGATGGAA
 ATTTCTCAGGTTCAATTACATGTTTGCAAAATGGATGGTCAGCACAACCAATTTGCATTAATTTTGTGA
 TATGCCTGTTTTTGAAGATTCAGAGCCAAGAGTAATGGCATGCGGTTTAAAGCTCCATGACACATTGGAC
 TACGAATGCTACGATGGATATGAAATCAGTTATGGAAACACCACAGTTCCATAGTGTGGTGAAGATG
 GGTGGTCCCATTCCCAACATGTTATAATTCTTCAGAAAAGTGTGGGCTCCTCCACCTATTAGCAATGG
 TGATACCCTCCTTTCTACTAAAAGTGTATGTGCCACAGTCAAGAGTCGAGTACCAATGCCAGTCTAC
 TATGAACCTCAGGTTCTAATTATGTAACATGTAGTAATGGAGAGTGGTGGAAACCAAGATGCATAC
 ATCCATGTATAAATACTGAAGAAAACATGAATAAAAATAACATACAGTTAAAAGGAAAAAGTGACATAAA
 ATATTATGCAAAAACAGGGGATACCATTGAATTTATGTGTAATTTGGGATATAATGCCAATACATCAGTT
 CTATCATTCAAGCAGTGTGTAGGAAGGCATAGTGAATACCCAGATGCGAA

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTAA

Protein Sequence:

>RC234476 representing NM_001201550
 Red=Cloning site Green=Tags(s)

MLLLINVILTLWSCANGQEVKPCDFPEIQHGGLYYKSLRRLYFPAAAGQSYSYCDQNFVTPSGSYWDY
 IHCTQDQWSPVPLRTRCSKSDVEIENGFISESSSIYILNEETQYNCKPGYATAEGNSSGSITCLQNGWS
 TQPICIKFCDMPVFNENSRKSNMWFKLHDTLDYECYDGYESSYGNTTDSIVCGEDGWSHLPTCYNSSEN
 CGPPPPISNGDTSFPPQKVVLPWSRVEYQCQSYELQGSKYVTCNSGDWSEPPRCISMKPCEFPFIQHG
 LYYENTRRPYFPVATGQSYSYCDQNFVTPSGSYWDYIHCTQDQWLPVPLRTRCSKSDIEIENGFISES
 SSIYILNKEIQYCKPGYATADGNSSGSITCLQNGWSAQPICIKFCDMPVFNENSRKSNMRFKLHDTLD
 YECYDGYEISYNTTGSIVCGEDGWSHFPTCYNSEKCGPPPPISNGDTSFLLKVVYVQSRVEYQCQSY
 YELQGSNYVTCNSGEWSEPPRCIHPCTIITENMNKNNIQLKGSIDIKYYAKTGDTIEFMCKLGYNANTSV
 LSFQAVCREGIVEYPRCE

SGP**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-RsrII

Cloning Scheme:


ACCN: NM_001201550

ORF Size: 1734 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001201550.3](#)

RefSeq Size: 2181 bp

RefSeq ORF: 1737 bp

Locus ID: 10877

UniProt ID: [Q92496](#)


Cytogenetics: 1q31.3

Protein Families: Secreted Protein

MW: 65.8 kDa

Gene Summary: This gene is a member of the complement factor H (CFH) gene family, and encodes one of the 5 CFH-related (CFHR) proteins. These 5 genes are closely linked to the CFH gene on chromosome 1q31-q32. The CFHRs are secreted plasma proteins synthesized primarily by the hepatocytes, and composed of highly-related short consensus repeats (SCRs). This protein enhances the cofactor activity of CFH, and is involved in complement regulation. It can associate with lipoproteins and may play a role in lipid metabolism. Alternatively spliced transcript variants encoding different isoforms (varying in the number of SCRs) have been described for this gene. [provided by RefSeq, Jan 2011]

Product images:



Circular map for RC234476