

Product datasheet for **RC234494**

ZNF195 (NM_001242842) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF195 (NM_001242842) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF195
Synonyms:	HRF1; ZNFP104
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234494 representing NM_001242842
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCTGGTGCTCAGACTCTGTTGACGTTCAAGGATGTGGCCATAGAATTCTCCCTGGAGGAGTGGAAAT
 GCCTGGACCTCGCTCAGCAGAATTTGTACAGGGATGTGATGTTGGAGAACACAGAACTTGTCTCCGT
 TGGTCTCACTGTCTGTAAGCCAGGCCTGATCACCTGCCTGGAGCAACGAAAAGAGCCCTGGAATGTGAAG
 AGACAGGAGGCAGCAGACGGACATCCAGCACTGTGCTCACCTGGGGTGTGCAAACCTGTGAAATGGTTTT
 TAGAGTTCAGATGCATCTTCTCACTAGCTATGTCTTCTCATTTCACCAAGACCTTCTGCCAGAGCAGGG
 CATACAAGATGCATTCACAAAAAGAATACTGAGAGGATATGGAAATTTGGCCTTGATAATTTATATTTA
 AGGAAAGACTGGGAAAGTTAGATGAGTGAAGTTGCAAAAAGATTATAATGGACTTAACCAATGTTTCAT
 CAACTACCCATAGCAAAATCTTTCAATATAATAAATATGTTAAAATCTTTGATAACTTTTCAAATTTACA
 TAGACGTAATATAAGTAATACTGGAGAGAAACCTTTCAAATGTCAAGAATGTGGCAAATCCTTTCAAATG
 CTCTCATTCTAACTGAACATCAGAAAAATTCACACTGGAAAAAATTCAAAAATGTGGAGAATGTGGCA
 AAACCTTTATCCAGTGCTCACACTTTACTGAACCTGAGAACATTGACACTGGAGAGAAAACCTTACAAGTG
 TCAAGAATGTAACAACGTCATTAACCTTGTCTCAGTCCTTACTAAAAATAGAATTTATGCCGGAGGGGAA
 CATTACAGATGTGAAGAAATTTGGCAAAGTATTTAACAGTGCTCCACCTTACTGAACATGAGCATGGTA
 CTGAGGAAAAACCTGCAAAATGAAGAGTGCAGCAGTGTCTTTATCTCTTGTCAAGCCTTTCTAATCA
 ACAGATGATTCTGTGGAGAGAAGCTCTCAAATGTGAAACATGGTACAAAGTTTTAACACAGCCCA
 AATCCTTCAAACACCAGAGAAATGAGATTGGAGGAAACCTTTCAAATGTGAGGAATGTGACAGCATCT
 TCAAGTGGTCTCAGACCTTACTAAACATAAGAGAATTCACACTGGTGAGAAACCATACAATGTGACGA
 ATGTGGGAAAGCCTATACACAGTCTCACACTCAGTGAACACAGGAGGATTCACACCGGAGAGAAACCC
 TACCAATGTGAAGAATGTGGGAAGTCTTCAGAACTTGCTCAAGCCTTTTCAACATAAGAGAACTCATT
 CTGAAGAAAAACCTACACGTGTGAAGAATGTGGCAACATCTTTAAGCAGTTATCAGACCTCACTAAGCA
 TAAGAAAACCCATACTGGAGAGAAGCCCTACAAATGTGACGAATGTGAAAAAACTTTACCCAGTCTCC
 AACCTTATTGTACATAAGAGAATTCATACTGGAGAGAAACCTACAAGTGTGAAGAATGTGGCAGAGTCT
 TCATGTGGTCTCAGACATTACCAACATAAGAAAACCCATACTGGAGAGAAACCTACAATGTGACGA
 ATGTGAAAAAACTTTACCCAGTCTCAAACCTTATTGTACATAAGAGAATTCATACTGGAGAGAAACCC
 TACAAGTGTGAAAAGTGTGGCAAAGCCTTACCCAGTTCTCACACCTGACTGTACATGAAAGCATTACATA
 CT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC234494 representing NM_001242842
 Red=Cloning site Green=Tags(s)

MAGAQTLTFRDVAIEFSLEEWKCLDLAQNLRYDVMLENYRNLFVSVGLTVCKPGLITCLEQRKEPWNVK
 RQEAADGHPALCSPGVLQTVKWFLEFRCIFSLAMSSHFTQDLLPEQGIQDAFPKRILRGYNGCLDNLVYL
 RKDWESLDECKLQKDYNLNQCSSTTHSKIFQYNKYVKIFDNFNLHRRNISNTGEKPFKQCQECGKSFQM
 LSFLTEHQKIHTGKKFQKCGEGKTFIQCSHFTEPENIDTGEKPYKQCENNVIKTCSVLTKNRIYAGGE
 HYRCEEFGKVFNQCSHLTEHEHGTEEKPKYECSVVISCSLSNQQMILAGEKLSKCEWYKGFNHSP
 NPSKHQRNEIGGKPFKCEECDSIFKWFSDLTKHKRIHTGEKPYKQCDECGKAYTQSSHLSEHRRIHTGEK
 YQCEECGKVFRTCSLSNHKRTHSEKPYTCEECGNIFKQLSDLTKHKKHTHTGEKPYKQCDECGKNFTQSS
 NLI VH KRIHTGEKPYKCEECGRVFMWFSIDITKHKHTHTGEKPYKQCDECGKNFTQSSNLI VH KRIHTGEK
 YKCEKCGKAFQFSHLTVHESIHT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001242842

ORF Size: 1752 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242842.2](#)

RefSeq Size: 3215 bp

RefSeq ORF: 1755 bp

Locus ID: 7748

UniProt ID: [O14628](#)


Cytogenetics: 11p15.4

Protein Families: Transcription Factors

MW: 68.1 kDa

Gene Summary: This gene encodes a protein belonging to the Krueppel C2H2-type zinc-finger protein family. These family members are transcription factors that are implicated in a variety of cellular processes. This gene is located near the centromeric border of chromosome 11p15.5, next to an imprinted domain that is associated with maternal-specific loss of heterozygosity in Wilms' tumors. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2012]

Product images:



Circular map for RC234494