

## Product datasheet for **RC234540**

### **RUNX1T1 (NM\_001198628) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	RUNX1T1 (NM_001198628) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	RUNX1T1
Synonyms:	AML1-MTG8; AML1T1; CBFA2T1; CDR; ETO; MTG8; ZMYND2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RC234540 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

RCATGATATCTGTCAAAGAAACACTTGGAGAGCACTGAGTTTAGTAATAGGTGACTGCCGAAAAAAGG  
 GAACCTTTGAATATTGTCAAGATCGTACTGAGAAGCACTCCACAATGCCAGACTCACCTGTGGATGTGAAG  
 ACGCAATCTAGGCTGACTCCTCCAACAATGCCACCTCCCCAACTACTCAAGGAGCTCCAAGAACCAGTT  
 CATTTACACCGACAACGTTAACTAATGGCAGGACCATTCTCTACAGCCTTGAATGGCGCCCCCTCACC  
 ACCCAATGGCTTCAGCAATGGCCCTTCTTCTTCTCTCTCTCTGCTAATCAACAGCTGCCCA  
 GCCTGTGGTGCCAGGCACTCAGCAAGTGAAAAGTTCTTACTACCTGCAGCAGTTTGGCAATGACA  
 TTTACCCGAGATAGGAGAAAGATTGCGACCCTCGTTCTGGACTAGTGAACCTCACTTTGACAATTGA  
 AGAATTTCACTCAAAGCTCAAGAAGCTACTAACTTCCACTGAGACCTTTTGCATCCCATTTTGAAG  
 GCCAATTGCCCTGCTGCAGCGTGGCTCCTCCACTGCGCAAGACTGGCCAAACAGAACCCTGCCAGT  
 ACCTCGCCAGCATGAACAGCTGCTTCTGGATGCCAGCACCTCACCTGTTGACTCCTCAGAGCTGCT  
 TCTCGATGTGAACGAAAACGGGAAGAGGCGAATCCAGACAGAACCAAAGAAAATGGCTTTGACAGAGAG  
 CCTTTCAGACTCAGAACATCCAAGCAAGCGACCATGCACTATTAGCCAGGCCAGCGGTACAGTCAAATA  
 ACGGTTATCCTACCAGCCCAATGGCCTGCCTCACCTACCCACCTCCACCTCAGCATTACCGTTTGGTA  
 TGATATGGCCATTGCCCACTACAGGACTCCTATCGACACCCAGCCACAGGGACCTCAGGGACAGA  
 AACAGACCTATGGGGTTGCATGGCACACGTCAAGAAGAAATGATTGATCACAGACTAACAGACAGAGAAT  
 GGGCAGAAGAGTGGAAACATCTTGACCATCTGTTAACTGCATAATGGACATGGTAGAAAAACAAGCGG  
 ATCTCTACCGTACTAAGCGGTGTCAAGAAGCAGACCGGGAAGAATTGAATTACTGGATCCGGCGGTAC  
 AGTGACGCCGAGGACTTAAAAAAGGTGGCGGAGTAGCAGCAGCCACTTAGGCAGCAGAGTCCCGTCA  
 ACCCAGACCCAGTTGCACTAGACGCGCATCGGGAATTCCTTACAGGCCTGCGTCTGGATACGTGCCAGA  
 GGAGATCTGGAAGAAAGCTGAGGAGGCCGTAATGAGGTGAAGCGCCAGGCGATGACGGAGCTGAGAAG  
 GCCGTGTCTGAGGCGGAGCGGAAAGCCACGACATGATCACACAGAGAGGCCAAGATGGAGCGCACGG  
 TCGCCGAGGCCAAACGGCAGGCGCGGAGGACGCACTGGCAGTTATCAATCAGCAGGAGGATTCAAGCGA  
 GAGTTGCTGGAATTGTGGCCGTAAGCGAGTGAACCTGCAGTGGCTGTAACACAGCCCGATACTGTGGC  
 TCATTTTCCAGCACAAGACTGGGAGAAGCACCATCACATCTGTGACAGACCTGCAGGCCACGACGC  
 AGGAGACACACCTGCAGTCAGTCTCTGTACGCCAACAGCGGGCTGGGAGCCGATGGACACACC  
 ACCAGCAGCCACTCCGAGGTCAACCACCCGGGAACCCCTTCACCATAGAGACAACCCCTCGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RC234540 protein sequence  
 Red=Cloning site Green=Tags(s)

X\*YLSKETLGEH\*V\*\*\*VTAGKKGTLNIVKIVLRSTPQCQTHLWM\*RRNLG\*LLQQCHLPQLLKELQEPV  
 HLHRQR\*LMARAILLQP\*MAPPHPMASAMGLPLLPPLWLINSCPQPVVPGNSAS\*KGSLLPCSSLAMT  
 FHPR\*EKEFAPSFWD\*\*TPL\*QLKNFIPNCKLLTSH\*DLLSSH\*F\*RPCTCPCSVSSTAQDWPNRTLPS  
 TSPSMNSCFWMPAPPHLLTPQSCFSM\*TKTGRGELQTEPKKMALTESLCTQNIQASDHALLAQASGTVQI  
 TAYPTSPMACTLPHLHLSITVWMIWPLPTTTGTPIDTPATGTSGETDLWGCMAHVKKK\*LITD\*QTEN  
 GQKSGNILTIC\*TA\*WTW\*KKQGDLSPY\*GGVKKQTGKN\*ITGSGGTVPRT\*KKVAVAATLGSRVPS  
 TQTQLH\*TRIGNSFTGLRLDTCQRRSGRKLRRPSMR\*SARR\*RSCRPRCLRRSGKPTT\*SQQRGPRWSAR  
 SPRPNGRRRTHWQLSISRRIQARVAGIVAVKRKPAVAVTQPDVAHFASFKTGRSTITSVDRPCRPS  
 RETHLQSAPLSRPTAGLGARWTHHQQLRGQPPREPLPP\*RQPL

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**


**ACCN:** NM\_001198628

**ORF Size:** 1812 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001198628.1](#), [NP\\_001185557.1](#)

**RefSeq Size:** 7284 bp


**RefSeq ORF:** 1815 bp

**Locus ID:** 862

**UniProt ID:** [Q06455](#)

<b>Cytogenetics:</b>	8q21.3
<b>Protein Families:</b>	Transcription Factors
<b>Protein Pathways:</b>	Acute myeloid leukemia, Pathways in cancer
<b>MW:</b>	67.6 kDa
<b>Gene Summary:</b>	This gene encodes a member of the myeloid translocation gene family which interact with DNA-bound transcription factors and recruit a range of corepressors to facilitate transcriptional repression. The t(8;21)(q22;q22) translocation is one of the most frequent karyotypic abnormalities in acute myeloid leukemia. The translocation produces a chimeric gene made up of the 5'-region of the runt-related transcription factor 1 gene fused to the 3'-region of this gene. The chimeric protein is thought to associate with the nuclear corepressor/histone deacetylase complex to block hematopoietic differentiation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2010]

**Product images:**



Circular map for RC234540