

## Product datasheet for **RC234580**

### **TNIP1 (NM\_001252385) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	TNIP1 (NM_001252385) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TNIP1
Synonyms:	ABIN-1; NAF1; nip40-1; VAN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RC234580 representing NM\_001252385  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGAAGGAGAGACCCTACCGGATCTACGACCTGGGGCAGCGTGCCCTCAGGAGAGGCATCCGCAG  
 CTTTTGAGCGCCTAGTGAAGGAGAATCCCGGCTGAAGGAAAAATGCAAGGGATAAAGATTTAGGGGA  
 GCTTTTGAAGAGTCCCAGATGGAAGCGACCAGGCTCCGGCAGAAGGCAGAGGAGCTAGTGAAGGACAAC  
 GAGCTGCTCCACCACCTTCTCCCTCCTTGGGCTCCTTCGACCCCTGGCTGAGCTCACAGGAAAGGACT  
 CAAATGTCACAGCATCTCCACAGCCCCTGCATGCCCCAGTGACAAGCCAGCACCAGTCCAGAAGCCTCC  
 ATCCAGTGGCACCTCCTCTGAATTTGAAGTGGTCACTCCTGAGGAGCAGAATTCACCAGAGAGCAGCAGC  
 CATGCCAATGCGATGGCGCTGGGCCCTGCCCGTGAGGACGGCAACCTGATGCTGCACCTGCAGCGCC  
 TGGAGACCACGCTGAGTGTGTGCGGAGGAGCCGGACCAGGCCAGCTTTCACCCACCTGGGCCGCAT  
 GGCCCTGGAGTTCACCCGACTGGCATCCAAGGTGCACAAGAATGAGCAGCGCACCTCCATTTGCGAGACC  
 CTGTGTGAGCAGCTTCGGAAGGAGAACGAGGCTCTGAAGGCCAAGTTGGATAAAGGCCTGGAACAGCGGG  
 ATCAGGCTGCCGAGAGGCTGCCGGAGGAAAAATTTGGAGCTCAAGAAGTTGTTGATGAGCAATGGCAACA  
 AGAGGGTGCCTCTGGCGGCCAGGCTCACCGAAGATGGAAGGGACAGGCAAGAAGGCAGTGGCTGGACAG  
 CAGCAGGCTAGTGTGACGGCAGGTAAGGTCCCAGAGGTGGTGGCCTTGGGCGCAGCCGAGAAGAAGTGA  
 AGATGCTGGAGCAGCAGCGCAGTGAAGTGAACAAGCAGTGGGACCAGCATTTCGGTCCAT  
 GAAGCAGCAGTATGAGCAGAAGTCACTGAGCTGCGTCAAGAGCTGGCTGATTTGAGAAGCAGGTGACT  
 GACCTGGAGGCCGAGCGGAGCAGAAGCAGCGTGACTTTGACCGCAAGCTCCTCTGGCCAAAGTCCAAGA  
 TTGAAATGGAGGAGACCACAAGGAGCAGCTGACAGCAGAGGCCAAGGAGCTGCGCCAAAAGTCAAGTA  
 CCTGCAGGATCAGCTGAGCCCACTCACCGACAGCGTGAGTACCAGGAAAAGGAGATCCAGCGGCTCAAC  
 AAGGCCCTGGAGGAAGCACTGAGCATCCAAACCCGCCATCATCTCCACCAACAGCATTGGGAGCCAG  
 AAGGAGCAGGGCCCTCCTAAGGAAACAGGAGCTGGTCACGCAGAATGAGTTGCTGAAACAGCAGGTGAA  
 GATCTTCGAGGAGGACTTCAGAGGGAGCGCAGTATCGTGAAGCAGTGAATGAGGAGAAGGAGAGCTG  
 AAGAAGCAAGTGGAGAAGCTGCAGGCCAGGTACCCCTGTCAAATGCCAGCTAAAAGCATTCAAAGATG  
 AGGAGAAGGCAAGAGAAGCCCTCAGACAGCAGAAGAGGAAAGCAAAGGCCTCAGGAGAGCGTTACCATGT  
 GGAGCCCCACCCAGAATCTCTGCGGGGCTACCCCTACGCCTACCCGCCATGCCAGCCATGGTGCCA  
 CACCATGGCTTCGAGGACTGGTCCAGATCCGCTACCCCTCCCCCATGGCCATGGAGACCCGCCCC  
 CACTCCCAACTCGCGCCTTCTCATCTGCACCAAGTTTTGTGCGAGCCGGAATACACCTGGCGTCTACCC  
 TGTGGAGGGGTTGAAATCCAAATCAGAGCTCCCAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**


>RC234580 representing NM\_001252385  
 Red=Cloning site Green=Tags(s)

MEGRGPYRIYDPGGSVPSGEASAAFERLVKENSRLKEKMQGIKMLGELLEESQMEATRLRQKAEELVKDN  
 ELLPPSPSLGSLDPLAELTGKDSNVTASPTAPACPSDKPAPVQKPPSSGTSSEFEVVTPEEQNSPSSS  
 HANAMALGPLPREDGNLMLHLQRLLETTLSVCAEEDPHGQLFTHLGRMALEFNRLASKVHKNEQRTSILQT  
 LCEQLRKENEALKAKLDKLEQRDQAAERLREENLELKKLLMSNGNKEGASGRPGSPKMEGTGKKAVAGQ  
 QQASVTAGKVPEVVALGAAEKVKMLEQQRSELEVNKQWDQHFRSMKQQYEQKITELRQKLADLQKQVT  
 DLEAEREQKQRDFDRKLLAKSKIEMEETDKEQLTAEAKELRQKVYLDQQLSPLTRQREYQEKIQRNLN  
 KALEEALSIQTPSPPTAFGSPGAGALLRKQELVTQNELLKQVKIFEEFQFRERSDRERMNEEKEEL  
 KKQVEKLQAQVTLNAQLKAFKDEEKAREALRQQRKAKASGERYHVEPHPEHLGAYPYAYPPMPAMVP  
 HHGFEDWSQIRYPMPMAMEHPPPLPNSRLFHLHQFCRSRNTPGVYYPVEGFETIQRAPK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001252385

**ORF Size:** 1857 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001252385.1](#), [NP\\_001239314.1](#)

**RefSeq Size:** 2950 bp

**RefSeq ORF:** 1860 bp

**Locus ID:** 10318


**UniProt ID:** [Q15025](#)

**Cytogenetics:** 5q33.1

**MW:** 70.5 kDa

**Gene Summary:** This gene encodes an A20-binding protein which plays a role in autoimmunity and tissue homeostasis through the regulation of nuclear factor kappa-B activation. Mutations in this gene have been associated with psoriatic arthritis, rheumatoid arthritis, and systemic lupus erythematosus. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]

**Product images:**



Circular map for RC234580