

Product datasheet for **RC234642**

PKC theta (PRKCQ) (NM_001242413) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PKC theta (PRKCQ) (NM_001242413) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PKC theta
Synonyms:	nPKC-theta; PRKCT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC234642 representing NM_001242413
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCGCCATTTCTTCGGATTGGCTTGTCACAACTTACTGCGGGTCTGCCAGTCTTGTACAGGGCGAGG
CTGTTAACCCCTTACTGTGCTGTGCTCGTCAAAGAGTATGTGCAATCAGAGAACGGGCAGATGTATATCCA
GAAAAAGCCTACCATGTACCCACCTGGGACAGCACTTTTGATGCCCATATCAACAAGGGAAGAGTCATG
CAGATCATTGTGAAAGGCAAAAACGTGGACCTCATCTCTGAAACCACCGTGGAGCTCTACTCGTGGCTG
AGAGGTGCAGGAAGAACAACGGGAAGACAGAAATATGGTTAGAGCTGAAACCTCAAGGCCGAATGCTAAT
GAATGCAAGATACTTTCTGAAAATGAGTGACACAAAGGACATGAATGAATTTGAGACGGAAGGCTTCTTT
GCTTTGCATCAGCGCCGGGTGCCATCAAGCAGGCAAAGGTCCACCACGTCAGTGCCACGAGTTTGTCTGGGCTGAACAAACAGGG
CTACCAGTGCCGACAATGCAATGCAGCAATTCACAAGAAGTGTATTGATAAAGTTATAGCAAAGTGCACA
GGATCAGCTATCAATAGCCGAGAAACCATGTTCCACAAGGAGAGATTCAAAATTGACATGCCACACAGAT
TTAAAGTCTACAATTACAAGAGCCCGACCTTCTGTGAACACTGTGGGACCTGCTGTGGGACTGGCAGC
GCAAGGACTCAAGTGTGATGCATGTGGCATGAATGTGCATCATAGATGCCAGACAAAGGTGGCAACCTT
TGTGGCATAAACCAGAAGCTAATGGCTGAAGCGCTGGCCATGATTGAGAGCACTCAACAGGCTCGTGTCT
TAAGAGATACTGAACAGATCTTCAGAGAAGGTCCGGTTGAAATTTGGTCTCCCATGCTCCATCAAAAATGA
AGCAAGGCCGCCATGTTTACCGACCCGGGAAAAGAGAGCCTCAGGGCATTTCCTGGGAGTCTCCGTTG
GATGAGGTGGATAAAATGTGCCATCTTCAGAACCTGAACTGAACAAAGAAAGACCATCTCTGCAGATTA
AACTAAAATTTGAGGATTTTATCTTGCACAAAATGTTGGGAAAAGGAAGTTTGGCAAGGCTTCCCTGGC
AGAATTCAGAAAACCAATCAATTTTTTCGCAATAAAGGCCTTAAAGAAAGATGTGGTCTTGATGGACGAT
GATGTTGAGTGCACGATGGTAGAGAAGAGATTTCTTCTTGGCCTGGGAGCATCCGTTTCTGACGCACA
TGTTTTGTACATTCAGACCAAGGAAAACCTCTTTTTTGTGATGGAGTACCTCAACGGAGGGGACTTAAT
GTACCACATCAAAGCTGCCACAAGTTCGACCTTTCAGAGCGACGTTTTATGCTGCTGAAATCATTCTT
GGTCTGCAGTTCCTTCATTCAAAGGAATAGTCTACAGGGACCTGAAGCTAGATAACATCCTGTTAGACA
AAGATGGACATATCAAGATCGCGGATTTTGAATGTGCAAGGAGAACATGTTAGGAGATGCCAAGACGAA
TACCTTCTGTGGGACACCTGACTACATCGCCCCAGAGCTCTTCGTGCGAGAACCTGAGAAGAGGCTGGC
GTGAGGGGAGACATCCGCCAGCACCTTTGTTTTGGGAGATCAACTGGGAGGAACCTGAACGGAAGGAGA
TTGACCCACCGTTCCGGCCGAAAGTGAATCACCATTTGACTGCAGCAATTTTCGACAAAGAATTCTAAA
CGAGAAGCCCCGGCTGTCAATTTGCCGACAGAGCACTGATCAACAGCATGGACCAGAATATGTTTCAGGAAC
TTTTCTTCATGAACCCCGGATGGAGCGGCTGATATCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234642 representing NM_001242413
Red=Cloning site Green=Tags(s)


MSPFLRIGLSNFDCGSCQSCQGEAVNPYCAVLVKEYVESENGQMYIQKKPTMYPPWDSTFDAHINKGRVM
 QIIYKGNVDLISETTVELYSLAERCNRKNGKTEIWLELKPQGRMLMNARYFLEMSDTKDMNEFETEGFF
 ALHQRRGAIKQAKVHHVKCHEFTATFFPQPTFCVSCHEFVWGLNKQGYQCRQCNAAIHKKCIDKVIKCT
 GSAINSRETMFHKERFKIDMPHRFKVYNYKSPTFCEHCGTLLWGLARQGLKCDACGMNVHHRQCQTKVANL
 CGINQKLMAEALAMIESTQQARCLRDTEQIFREGPVEIGLPCSIKNEARPPCLPTPGKREPQGISWESPL
 DEVDKMCHLPEPELNKERPSLQIKLKIEDFILHKMLGKGSFGKVFLLAEFKTNQFFAIKALKKDVVLMDD
 DVECTMVEKRVLSLAWHEPFLTHMFCFQTKENLFFVMEYLNNGDLMYHIQSCHKFDLSRATFYAAEIIIL
 GLQFLHSKGIVYRDLKLDNILLDKDGHKIADFGMCKENMLGDAKTNTFCGTPDYIAPEL FVREPEKRLG
 VRGDIRQHPLFREINWEELERKEIDPPFRPKVSPFDCSNFDKEFLNEKPRLSFADRALINSMDQNMFRN
 FFSMNPGMERLIS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001242413

ORF Size: 1929 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001242413.2 , NP_001229342.1
RefSeq Size:	3106 bp
RefSeq ORF:	1932 bp
Locus ID:	5588
UniProt ID:	Q04759
Cytogenetics:	10p15.1
Protein Families:	Druggable Genome, Protein Kinase, Transcription Factors
Protein Pathways:	Adipocytokine signaling pathway, T cell receptor signaling pathway, Tight junction, Vascular smooth muscle contraction
MW:	74.7 kDa
Gene Summary:	<p>Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play a distinct role. The protein encoded by this gene is one of the PKC family members. It is a calcium-independent and phospholipid-dependent protein kinase. This kinase is important for T-cell activation. It is required for the activation of the transcription factors NF-kappaB and AP-1, and may link the T cell receptor (TCR) signaling complex to the activation of the transcription factors. [provided by RefSeq, Jul 2008]</p>

Product images:



Circular map for RC234642