

Product datasheet for **RC234652**

CLPB (NM_001258393) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLPB (NM_001258393) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CLPB
Synonyms:	ANKCLB; HSP78; MEGCANN; MGCA7; SKD3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC234652 representing NM_001258393
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGGGGTCCCTGGTGTGAGGAGAAAAGCACTGGCGCCACGGCTACTCCTCCGGCTGCTCAGGTCCC
 CAACGCTCCGGGGCCATGGAGGTGCTTCCGGCCGGAATGTGACTACTGGGAGTCTCGGGGAGCCGAGTG
 GCTGAGGGTAGCCACCGGGGGCGCCCTGGAACATCGCCGGCCTTGTCTCCGGACGTGGGGCAGCCACC
 GGGGGGCGCCAGGGAGGACGCTTCGATACCAAATGCCTCGCGGCTGCCACTTGGGACGCCTTCCTGGTC
 CCGAAGAAACACTCCCAGGACAGGACAGCTGGAACGGGGTCCCCAGCAGGGCCGGACTGGGCATGTGCGC
 CCTGGCCGCGAGCTGGTGGTTCATTGCTACAGCAAGAGTCCGTCCAACAAGGATGCAGCCCTGTTGAA
 GCTGCCCGTGCCAACAATATGCAAGAAGTCAGCAGTGTGGTACAGGTCTGCTTGTGCTGGGGCTGATC
 CAAACCTTGGAGATGATTTAGCAGTGTTCACAAGACTGCCAAGGAACAGGGAATCCATTCTTTGGAAGT
 CCTGATCACCCGAGAGGATGACTTCAACAACAGGCTGAACAACCGCGCCAGTTTCAAGGGCTGCACGGCC
 TTGCACTATGCTGTTCTTGTGCTGACTACCGCACTGTCAAGGAGCTGCTTGATGGAGGAGCCAACCCCC
 TGCAGAGGAATGAAATGGGACACACACCCTTGGATTATGCCCGAGAAGGGGAAGTATGAAGCTTCTGAG
 GACTTCTGAAGCCAAGTACCAAGAGAAGCAGCGGAAGCGTGAGGCTGAGGAGCGCGCCGCTTCCCCCTG
 GAGCAGCGACTAAAGGAGCACATCATTGGCCAGGAGAGCGCCATCGCCACAGTGGGTGCTGCGATCCGGA
 GGAAGGAGAATGGCTGGTACGATGAAGAACACCCTCTGGTCTTCTCTTCTTGGGATCATCTGGAATAGG
 AAAACAGAGCTGGCCAAGCAGACAGCCAAATATATGCACAAAGATGCTAAAAAGGGCTTCATCAGGCTG
 GACATGTCCGAGTTCAGGAGCGACACGAGGTGGCCAAGTTTATTGGGTCTCCACCAGGCTACGTTGGCC
 ATGAGGAGGGTGGCCAGCTGACCAAGAAGTTGAAGCAGTGCCCAATGCTGTGGTCTCTTGTATGAAGT
 AGACAAGGCCCATCCAGATGTGCTCACCATCATGCTGCAGCTGTTTGATGAGGGCCGGCTGACAGATGGA
 AAAGGGAAGACCATTTGATTGCAAGGACGCCATTTTCATCATGACCTCCAATGTGGCCAGCGACGAGATCG
 CACAGCACGCGCTGCAGCTGAGGCAGGAAGCTTTGGAGATGAGCCGTAACCGTATTGCCGAAAACCTGGG
 GGATGTCCAGATAAGTGACAAGATCACCATCTCAAAGAACTTCAAGGAGAATGTATTGCGCCTATCCTG
 AAAGCTCACTTCCGAGGGATGAGTTTCTGGGACGGATCAATGAGATCGTCTACTTCTCCCTTCTGCGC
 ACTCGGAGCTCATCAACTCGTCAACAAGGAATAAACTTCTGGGCCAAGAGAGCCAAGCAAAGGCACAA
 CATCACGCTGCTCTGGGACCGCGAGGTGGCAGATGTGCTGGTCGACGGCTACAATGTGCACTATGGCGCC
 CGCTCCATCAAACATGAGGTAGAACGCCGTGTGGTGAACCAGCTGGCAGCAGCCTATGAGCAGGACCTGC
 TGCCAGGGGGCTGACTTTGCGCATCACGGTGGAGGACTCAGACAAGCAGCTACTCAAAGCCCAGAACT
 GCCCTACCCCAGGCTGAGAAGCGCCTCCCAAGCTGCGTCTGGAGATCATCGACAAGGACAGCAAGACT
 CGCAGACTGGACATCCGGGCACCACTGCACCCTGAGAAGGTGTGCAACACCATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC234652 representing NM_001258393
Red=Cloning site Green=Tags(s)

MLGSLVLRKALAPRLLRLRSPRTLRGHGGASGRNVTTGSLGEPQWLRVATGGRPGTSPALFSGRGAAT
GGRQGGFRFDTKLAAATWGRLPGPEETLPGQDSWNGVPSRAGLGMCAALAAVHCHYCKSPSNKDAALLE
AARANNMQEVSSVVQVLLAAGADPNLGDDFSSVYKTAKEQGIHSLLEVLITREDDFNRLNLRASFKGCTA
LHYAVLADDYRTVKELLDGGANPLQRNEMGHPTLDYAREGEVMKLLRTSEAKYQEQRKREAEEERRRFP
EQLRKEHIIGQESAIATVGAAIRRKENGWYDEEHLVFLFLGSSGIGKTELAKQTAKYMHKDAKKGFI
DMSEFQERHEVAKFIGSPPGYVGHEEGQLTKKLKQCPNAVVLDFEVDKAHPDVLTIMLQLFDEGRLTDG
KGKTIDCKDAIFIMTSNVASDEIAQHALQLRQEALMSNRNRIENLGDVQISDKITISKNFKENVIRPIL
KAHFRRDEFLGRINEIVYFLPFCHSEL IQLVNKELNFWAKRAKQRHNITLLWDREVADVLVDGYNVHYGA
RSIKHEVERRVYNQLAAAYEQDLLPGGCTLRITVEDSDKQLLKSPELPSPQAEKRLPKLRLEIIDKDSKT
RRLDIRAPLHPEKVCNTI

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM 001258393

ORF Size: 1944 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001258393.3](#)

RefSeq Size: 3095 bp

RefSeq ORF: 1947 bp

Locus ID: 81570


UniProt ID: [Q9H078](#)

Cytogenetics: 11q13.4

MW: 72.7 kDa

Gene Summary: This gene belongs to the ATP-ases associated with diverse cellular activities (AAA+) superfamily. Members of this superfamily form ring-shaped homo-hexamers and have highly conserved ATPase domains that are involved in various processes including DNA replication, protein degradation and reactivation of misfolded proteins. All members of this family hydrolyze ATP through their AAA+ domains and use the energy generated through ATP hydrolysis to exert mechanical force on their substrates. In addition to an AAA+ domain, the protein encoded by this gene contains a C-terminal D2 domain, which is characteristic of the AAA+ subfamily of Caseinolytic peptidases to which this protein belongs. It cooperates with Hsp70 in the disaggregation of protein aggregates. Allelic variants of this gene are associated with 3-methylglutaconic aciduria, which causes cataracts and neutropenia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015]

Product images:



Circular map for RC234652