

Product datasheet for **RC234655**

ZNF816A (ZNF816) (NM_001202456) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | ZNF816A (ZNF816) (NM_001202456) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | ZNF816A |
| Synonyms: | ZNF816A |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

ORF Nucleotide
Sequence:

>RC234655 representing NM_001202456
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTTACGTGAGGAAGCCACCAAGAAGAGCAAAGAAAAGGAGCCAGGGATGGCTCTTCTCAGGGACGCT
TGACTTTCAGGGATGTGGCTATAGAGTTCTTTGGAGGAGTGGAAATGCCTGAACCTGCACAGAGGGC
TTTATACAGGGCTGTGATGTTGGAGAACTACAGGAACCTGGAGTTTGTGGATAGCTCTTTAAAATCCATG
ATGGAGTTCATCAACCAGGCACAGTATTACAGGAGAAGTGATCCACACAGGGACGTTGCAAAGACATA
AAAGTCATCACATTGGAGATTTTGTCTCCAGAAATGAAGAAAGATATTCATCACTTTGAGTTTCAGTG
GCAAGAAGTTGAAAGAAATGGCCATGAAGCACCCATGACAAAAATCAAAAAGTTGACTGGTAGTACAGAC
CGAAGTGATCACAGGCATGCTGGAACAAGCCTATTAAGATCAGCTTGGATTAAGCTTTCATTCGCATC
TGCCCTGAACTCCACATGTTTCAGACTAAAGGGAAAATTAGTAACCAATTGGACAAGTCTATCGGTGCTTC
CTCAGCTTCAGAATCCCAAAGAATTTCTGTAGGCTCAAACTCATATTTCTAATAAGTATGGGAAGAAT
TTCTCCATTCTTCACTCACAAAAATACAGGAAATATGCATGAGAGAAAAACCTTGCCAAAGTAATGAGT
GTGGCAAAGCCTTTAATTATAGCTCACTCTTAAGGAGACACCACATAACCCATTCAAGAGAGAGAGAATA
TAAATGTGATGTATGTGGCAAGATCTTTAATCAGAAGCAATACATTGTATATCATCATAGATGTCACACT
GGTGAGAAAACCTTACAAGTGAATGAGTGTGGGAAGACCTTCACTCAGATGTCATCCCTTGTATGCCATC
GTAGACTTCATACTGGAGAGAAACCTTACAAGTGAATGAGTGTGGCAAGACCTTCAAGTGAATGAGTGTGGCAAGACTTTT
GGTCGAAATTCAGCCCTTGAATTCATAAGGCAATTCATACTGGAGAGAAACCTTACAAGTGAATGAGT
GTGGCAAGACCTTCAGTCAGAAATCATCCCTTCAATGCCATCATATACTTCACACTGGAGAGAAACCTTA
CAAATGTGAAGAATGTGACAATGTTTACATTCGCAGATCACACCTTGAAAGACATAGGAAGATTCATACT
GGAGAGGGATCATACAAATGTAAGGTTTGTGACAAGGTTTTCCGGAGTGATTCATACCTTGCAGAACATC
AGAGAGTTCATACTGGAGAGAAACCATAACAAGTGAATAAATGTGGCAGGAGTTTCAGTCGGAAGTCATC
CCTTCAATACCATCATACACTTCACACTGGAGAGAAACCTTACACATGAATGAATGTGGCAAGGTTTTT
AGTCGAAGAGAAAACCTTGCACGTCATCATAGACTTCATGCTGGAGAGAAACCTTACAAATGTGAAGAAT
GTGACAAAAGTTTTAGTCGCAGATCACACCTTGAAGACATAGGAGAATTCATACTGGGGAAAAACCATA
CAAATGTAAGGTTTGTGACAAGGCTTCCGGAGTGATTCATGCCTTGCAAACCATACGAGAGTTCATACT
GGAGAGAAACCTTACAAGTGAATAAATGTGCGAAGGTTTTAATCAAAAAGGAATCCTTGCACAACATC
AGAGAGTTCATACTGGAGAGAAACCTTACAAGTGAATGAATGTGGCAAGGTTTTAATCAAAAAGCAAG
CCTTGCAAACATCAGAGAGTTCATACTGCAGAGAAACCTTACAAGTGAATGAGTGTGGCAAAGCCTTT
ACTGGACAGTCAACACTTATTCACCATCAAGCAATCCATGGGTGTAGGAAAACCTTACAAATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234655 representing NM_001202456
Red=Cloning site Green=Tags(s)


```
MLREEATKKSKEKEPGMALPQGRLTFRDVAIEFSLEEWKCLNPAQRALYRAVMLENYRNLEFVDSLSKSM
MEFSSTRHSITGEVIHTGLQRHKSHHIGDFCFPEMKKDIHHEFWQVEVRNGHEAPMTKIKKLTGSTD
RSDHRHAGNKPIKDQLGLSFHSHLPELHMFQTKGKISNQLDKSIGASSASESQRISCRKLTHTISNKYKGN
FLHSSFTQIQEICMREKPCQSNCECGKAFNYSSLLRRHHIHTSREREYKCDVCGKIFNQKQYIVYHHRCHT
GEKTYKNECGKTFQMSSLVCHRRLLHTGEKPYKNECGKTFSEKSSLRCHRRLLHTGEKPYKNECGKTF
GRNSALVIHKAihtGEKPYKNECGKTFsQKSSLQCHHILHTGEKPYKCEECDNVYIRRSHLERHRKIHT
GEGSYKCKVCDKVFRRSDSYLAHQVHTGEKPYKCNKCGRSFSRKSSLQYHHTLHTGEKPYKNECGKVF
SRRENLARHRLHAGEKPYKCEECDKVFRRSHLERHRIHTGEKPYKCKVCDKAFRRSDSLANHTRVHT
GEKPYKCNKCAKVFNQKILAQHQVHTGEKPYKNECGKVFNQKASLAKHQVHTAEKPYKNECGKAF
TGQSTLIHHQAIHGCRETLMQ
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8120_c10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001202456

ORF Size: 1953 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001202456.2](#)

RefSeq Size: 2619 bp

RefSeq ORF: 1956 bp

Locus ID: 125893


UniProt ID: [Q0VGE8](#)

Cytogenetics: 19q13.41

MW: 75.7 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC234655