

Product datasheet for **RC234669**

ANKRD6 (NM_001242814) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKRD6 (NM_001242814) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ANKRD6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC234669 representing NM_001242814
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCCAGCAAGATGCGGTGCTGCACTTTCAGAGCGCCTTCTCGTAGCTGCGTACAAAGGCCAAACAG
AGAATGTGGTTTCAGCTCATCAACAAGGGCGCCAGGGTAGCGGTTACCAAGCATGGCCGGACTCCCTGCA
TCTTGCTGCCAATAAGGGCCATCTTCTGTGGTCCAGATCTTGCTGAAGGCTGGCTGCGACCTTGATGTC
CAGGATGATGGGGACCAGACCGCCTTGACCGGGCCACAGTGGTGGGAAACACGGAGATCATCGCGGCGC
TCATCCACGAAGGGTGTGCCCTGGACAGACAAGACAAGGCGGGAAACACAGCTCTGCACCTGGCCTGCCA
GAACAGCCACTCCAGAGCACGCGCTCCTCTGCTGGCCGGTCCCAGCTGACCTCAAAAATAATGCA
GGAGACACCTGTTTGCAGTGTGTCGCGCTATAATCACTTGTCCATCATTAGGCTCCTCCTCACTGCTT
TCTGTTCTGTCCATGAAAAGAACCAGGCTGGAGACACAGCACTTACGTTGCTGCTGCCATAATCACAA
GAAGGTGGCCAAAATCTTACTGGAAGCCGAGCAGATACGACATTGTTAAACAATGTCTTGCGCTTCACT
CGTGGGCGAAGCCTGAGGAAAAAGAGAGAGAGGCTCAAGGAAGAGAGGAGAGCCAGTCTGTGCCAAGAG
ATGAGGTGGCCAAAGCAAGGGCAGTGTCTCAGCAGGAGACACCCCCAGCAGTGAACAGGCTGTGGCCAG
AAAAGAAGAAGCCAGAGAAGAGTTCTGTGACGCTCCCCAGAACCAGAGCAAAGGATGACAGGAGGAGA
AAGTCAAGGCCAAAGGTGTGACGATTTTCTGACCCACCCACCAGCCGACCAACAGCCTGGACACCAGA
AGAACCTGCATGCTCATAATCACCTAAAAAGAGGAACAGGCATCGGTGTTTATCCCCACCCCAACCCCA
TGAGTTCAGGGCGTATCAGCTCTACACATTGTACCGGGCAAGGATGGGAAAGTATGACAGCCACCAATA
AATGGTTGTCGATGTGAACCTCTAATCAACAAGCTGGAGAATCAGTTGGAGGCTACTGTGGAGGAGATAA
AAGCAGAGCTGGGATCGTTTCAGGACAAAATGAATACAAAGCTGGGCGAGATGGAGAATAAGACCCAGCA
CCAAATGCGTGTGGGACAAGCTGATGGTTGAGCGACTTTCTGCAGAGAGGACGGAGTGCCTGAACCGC
CTGCAACAGCACTCAGACACAGAGAAGCATGAGGGGGAGAAACGACAGATATCCTTGGTGGATGAATTAA
AAACCTGGTGCATGTTAAAGATTCAGAATCTGGAGCAGAAGCTTTCTGGAGATTCTAGGGCTGCAGAGC
TAAATCCACACCATCTACTTGTGTGGACCAATTAGTGGTACTGCAGGTCCAGCAGCAGCTTCCGACAGC
TCCCTCCAGTGGTTAGGCCAAAGAGAAGGCCCTCAACTCCACTGTACCCAGAGACTCCAGCAGGAGC
TGTCGTTCTGACTGTACAGGCTCCCGACTGAGAAACGTCAAGGTCCAGACAGCCTTGCTACCCATGAA
TGAGGCAGCCAGATCTGATCAGCAGGCTGGGCCCTGCGTCAACAGAGGCACTCAAATAAGAAAGTCTGGG
AAGAGTGGCCAAACAAGGCATCGTGCCAGCAACCCGACCCAGCAGCACCTGTGGGCGAGCCGCCACCAG
CCACAGGCAGCGAGCAGACTGGCCCTCACATTCGGGACACCTCCAAGCTCTGGAGCTTACCCAGTATTT
TTTTGAGGCTGTTTCTACCCAGATGAAAAGTGGTATGAAAGGAAGATTGAAGAAGCACGAGCCAAAGCC
AATCAGAAAGCCAGCAAGATAAGGCTACATTGAAGGAACACATTAAGTTTAGAAGAGGAACTTGCCA
AACTAAGGACTAGGGTGCAGAAGGAAAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234669 representing NM_001242814
Red=Cloning site Green=Tags(s)

```
MSQQDAVAALSERLLVAAYKGQ TENVVQLINKGARVAVTKHGRTPLHLAANKGHLPVVQILLKAGCDLDV
QDDGDQTLHRATVVGNT EIIAALIHEGCALDRQDKAGNTALHLACQNSHSQSTRVLLLAGSRADLKNNNA
GDTCLHVAARYNHL SIIRLLL TAFCSVHEKNQAGDTALHVAALNHKKVAKILLEAGADTTIVNNVLRFS
RGRSLRKKRERLKEERRAQSVPRDEVAQSKGSVSAGDTPSSEQAVARKEEAREEFLSASPEPRAKDDRRR
KSRPKVSAFSDPTPPADQQPGHQKLNHANNHPKRNHRCSPPPPHEFRAYQLYTLYRGKDGKVMQAPI
NGRCRCEPLINKLENQLEATVEE IKAELGSVQDKMNTKLGQMENKTQHQMRLDKLMVERLSAERTECLNR
LQQHSDTEKHEGEKRQISLVDELKTWMLKIQNLEQKLSGDSRACRAKSTPSTCVDQLVVTAGPAAASDS
SPPVVRPEKALNSTATQRLQELSSSDCTGSRLRNKVQQTALLPMNEAARSDQQAGPCVNRGTQTKKSG
KSGPTRHRAQQPAASSTCGPPPATGSEQTGPHIRDTSQALELTQYFFAVSTQMEKWYERKIEEARSQA
NQKAQQDKATLKEHIKSLEELAKLRTRVQKEN
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001242814

ORF Size: 1989 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242814.1](#), [NP_001229743.1](#)

RefSeq Size: 4961 bp

RefSeq ORF: 1992 bp

Locus ID: 22881


UniProt ID: [Q9Y2G4](#)

Cytogenetics: 6q15

MW: 73.7 kDa

Gene Summary: Recruits CKI-epsilon to the beta-catenin degradation complex that consists of AXN1 or AXN2 and GSK3-beta and allows efficient phosphorylation of beta-catenin, thereby inhibiting beta-catenin/Tcf signals.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC234669