

Product datasheet for **RC234689**

CDHH (CDH13) (NM_001220489) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDHH (CDH13) (NM_001220489) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CDHH
Synonyms:	CDHH; P105
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC234689 representing NM_001220489
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCAGCCGAGAACTCCGCTCGTTCTGTGCGTTCTCCTGTCCCAGGTGCTGCTAACATCTGCAGAAG
ATTTGGACTGCACTCCTGGATTTTCAGCAGAAAGTGTTCATATCAATCAGCCAGCTGAATTCATTGAGGA
CCAGTCAATTCTAAACTTGACCTTCAGTACTGTAAAGGAAACGACAAGCTACGCTATGAGGTCTCGAGC
CCATACTTCAAGGTGAACAGCGATGGCGGCTTAGTTGCTCTGAGAAACATAACTGCAAGTGGGCAAAACTC
TGTTCTGCTCCATGCACGGACCCCCATGCGGAAGATATGGCAGAACTCGTGATTGTCGGGGGAAAGACAT
CCAGGGCTCCTTGCAAGTAGTCGATAGTGACAGGCCAGAAAGTCCAAGTCCGGCTCACTGAAAAGGGA
GTGGATCAAGAGCCTAAAGGAATTTTCAGAATCAATGAGAACACAGGGAGCGTCTCCGTGACACGGACCT
TGGACAGAGAAGTAATCGCTGTTTATCAACTATTTGTGGAGACCACTGATGTCAATGGCAAAACTCTCGA
GGGGCCGGTGCCTCTGGAAGTCATTGTGATTGATCAGAATGACAACCGACCGATCTTTCCGGGAAGCCCC
TACATCGGCCACGTGATGGAAGGGTCACCCACAGGCACACAGTATGCGGATGACAGCCTTTGATGCGAG
ATGACCCAGCCACCGATAATGCCCTCCTGCGGTATAATATCCGTGACGAGACGCCTGACAAGCCATCTCC
CAACATGTTTACATCGATCCTGAGAAAGGAGACATTGTCACTGTTGTGTACCTGCGCTGCTGGACCGA
GAGACTCTGAAAATCCCAAGTATGAACTGATCATCGAGGCTCAAGATATGGCTGGACTGGATGTTGGAT
TAACAGGCACGGCCACAGCCACGATCATGATCGATGACAAAAATGATCACTCACCAAAATTCACCAAGAA
AGAGTTTCAAGCCACAGTCGAGGAAGGAGCTGTGGGAGTTATTGTCAATTTGACAGTTGAAGATAAGGAT
GACCCACACAGGTGCATGGAGGGCTGCCTACCCATCATCAACGAAACCCCGGGCAGAGCTTTGAAA
TCCACACCAACCCTCAAACCAACGAAGGGATGCTTTCTGTTGTCAAACCATGGACTATGAAATTTCTGC
CTTCCACACCCTGCTGATCAAAGTGAAAAATGAAGACCCACTCGTACCCGACGTCTCCTACGGCCCCAGC
TCCACAGCCACCGTCCACATCACTGTCTGGATGTCAACGAGGGCCAGTCTTCTACCCAGACCCCATGA
TGGTGACCAGGCAGGAGGACCTCTCTGTGGCAGCGTGCTGCTGACAGTGAATGCCACGGACCCCGACTC
CCTGCAGCATCAAACCATCAGGTATTCTGTTTACAAGGACCCAGCAGGTTGGCTGAATATTAACCCCATC
AATGGGACTGTTGACACCACAGCTGTGCTGGACCGTGAGTCCCCATTTGTCGACAACAGCGGTGACTG
CTCTCTTCTGGCAATTGACAGTGGCAACCCTCCCGCTACGGGCACTGGGACTTTGCTGATAACCCTGGA
GGACGTGAATGACAATGCCCGTTCATTTACCCACAGTAGCTGAAGTCTGTGATGATGCCAAAAACCTC
AGTGTAGTCATTTGGGAGCATCAGATAAGGATCTTACCCGAATACAGATCCTTTCAAATTTGAAATCC
ACAACAAGCTGTTCTGATAAAGTCTGGAAGATCTCAAGATCAACAATACACACGCCCTGGTAAGCCT
TCTTCAAATCTGAACAAGCAAACCTACAACCTGCCATCATGGTGACAGATTCAGGGAAACCCCATG
ACGAATATCACAGATCTCAGGGTACAAGTGTGCTCCTGCAGGAATTCAAAGTGGACTGCAACGCGGCAG
GGCCCTGCGCTTACGCTGCCCTCAGTCTGCTCCTCAGCCTTTCAGCTTAGCTTGTCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234689 representing NM_001220489
Red=Cloning site Green=Tags(s)


MQPRTPLVLCVLLSQVLLLTSAEDLDCTPGFQQKVFHINQPAEFIEDQSILNLTFSDCKGNDKLRYEVSS
 PYFKVNSDGGVALRNITAVGKTLFVHARTPHAEDMAELVIVGGKDIQGSLOVVDSDRPERSKFRLTGKG
 VDQEPKGI FRINENTGSVSVTRTLDREVI AVYQLFVETTDVNGKLEGPVPLEVIDQNDNRPIFRE
 YIGHVMEGSPGTGTTVMRMTAFDADDPATDNALLRYNIRQQTPDKPSPNMFYIDPEKGDIVTVVSPALLDR
 ETLENPKYELIIEAQDMAGLDVGLTGTATATIMIDDKNDHSPKFTKKEFQATVEEGAVGVI VNLTVEDKD
 DPTTGAWRAAYTIINGNPGQSFEIHTNPQTNEGMLSVVKPLDYEISAFHTLLIKVENEDPLVPDVSYPGS
 STATVHITVLDVNEGPVFYDPMMVTRQEDLSVGSVLLTVNATDPDSLQHQITIRYSVYKDPAGWLNINPI
 NGTVDTTAVLDRESFV DNSVYALFLAIDSGNPPATGTGTLITLEDVNDNAPFIYPTVAEVCDDAKNL
 SVVILGASDKDLHPNTDPFKFEIHKQAVPDKVWKISKINNTHALVSLQLNLNKANYNLPIMVTDSGKPPM
 TNITDLRVQVCSRNKSKVDCNAAGALRFSLPVLLLSLFLACL

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001220489

ORF Size: 2022 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001220489.2](#)

RefSeq Size: 3911 bp

RefSeq ORF: 2025 bp

Locus ID: 1012


UniProt ID: [P55290](#)

Cytogenetics: 16q23.3

MW: 74.2 kDa

Gene Summary: This gene encodes a member of the cadherin superfamily. The encoded protein is localized to the surface of the cell membrane and is anchored by a GPI moiety, rather than by a transmembrane domain. The protein lacks the cytoplasmic domain characteristic of other cadherins, and so is not thought to be a cell-cell adhesion glycoprotein. This protein acts as a negative regulator of axon growth during neural differentiation. It also protects vascular endothelial cells from apoptosis due to oxidative stress, and is associated with resistance to atherosclerosis. The gene is hypermethylated in many types of cancer. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, May 2011]

Product images:



Circular map for RC234689