

Product datasheet for **RC234734**

EZH2 (NM_001203249) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EZH2 (NM_001203249) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	EZH2
Synonyms:	ENX-1; ENX1; EZH2b; KMT6; KMT6A; WVS; WVS2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC234734 representing NM_001203249
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGCCAGACTGGGAAGAAATCTGAGAAGGGACCAGTTTGTGGCGGAAGCGTGTAAAATCAGAGTACA
TGCAGACTGAGACAGCTCAAGAGTTTCAGACGAGCTGATGAAGTAAAGAGTATGTTTAGTTCCAATCGTCA
GAAAATTTTGGAAAGAACGGAAATCTTAAACCAAGAATGGAAACAGCGAAGGATACAGCCTGTGCACATC
CTGACTTCTTGTTCGGTGACCAGTGACTTGGATTTTCCAACACAAGTCATCCATTAAAGACTCTGAATG
CAGTTGCTTCAGTACCCATAATGTATTCTTGGTCTCCCTACAGCAGAATTTTATGGTGGAAAGTAAAC
TGTTTTACATAACATTCTTATATGGGAGATGAAGTTTTAGATCAGGATGGTACTTTTCATTGAAGAACTA
ATAAAAAATTATGATGGGAAAGTACACGGGGATAGAGAATGTGGGTTTATAATGATGAAATTTTTGTGG
AGTTGGTGAATGCCCTTGGTCAATATAATGATGATGACGATGATGATGATGGAGACGATCCTGAAGAAAG
AGAAGAAAAGCAGAAAGATCTGGAGGATCACCGAGATGATAAAGAAAGCCGCCACCTCGGAAATTTCT
TCTGATAAAATTTTTGAAGCCATTTCTCAATGTTCCAGATAAGGGCACAGCAGAAGAACTAAAGGAAA
AATATAAAGAAGCTCACCGAACAGCAGCTCCAGGGCAGCTTCTCCTGAATGTACCCCAACATAGATGG
ACCAAAATGCTAAATCTGTTTCAGAGAGAGCAAAGCTTACACTCCTTTCATACGCTTTTCTGTAGGCGATG
TTTAAATATGACTGCTTCTACATCCTTTTCATGCAACACCCAACACTTATAAGCGGAAGAACACAGAAA
CAGCTCTAGACAACAACCTTGTGGACCACAGTGTACCAGCATTTGGAGGGAGCAAAGGAGTTTGTCTGC
TGCTCTCACCGCTGAGCGGATAAAGACCCCAAAAACGTCCAGGAGGCCGAGAAGAGGACGGCTTCCC
AATAACAGTAGCAGGCCAGCACCCCAACATTAATGTGCTGGAATCAAAGGATACAGACAGTATAGGG
AAGCAGGGACTGAAACGGGGGAGAGAACAATGATAAAGAAGAAGAAGAGAAGAAAGATGAACTTCGAG
CTCCTCTGAAGCAAATTCGGTGTCAAACACCAATAAAGATGAAGCAAATATTGAACCTCCTGAGAAT
GTGGAGTGGAGTGGTGAAGCCTCAATGTTTAGAGTCCTCATTGGCACTTACTATGACAATTTCTGTG
CCATTGCTAGGTTAATTGGGACCAAAACATGTAGACAGGTGATGAGTTTAGAGTCAAAGAATCTAGCAT
CATAGCTCCAGCTCCCGTGAGGATGTGGATACTCCTCAAGGAAAAAGAAGAGGAAACACCGGTTGTGG
GCTGCACACTGCAGAAAGATACAGCTGAAAAAGGGTCAAACCGCTTCCGGGATGCCGCTGCAAAGCAC
AGTGCAACACCAAGCAGTCCCGTGCTACCTGGCTGTCGAGAGTGTGACCTGACCTCTGTCTTACTTG
TGGAGCCGCTGACCATTGGGACAGTAAAAATGTGCTCCTGCAAGAAGTGCAGTATTCAGCGGGCTCCAAA
AAGCATCTATTGCTGGCACCATCTGACGTGGCAGGCTGGGGGATTTTATCAAAGATCCTGTGCAGAAA
ATGAATTCATCTCAGAATACTGTGGAGAGATTATTTCTCAAGATGAAGCTGACAGAAGAGGGAAAGTGTA
TGATAAATACATGTGCAGCTTCTGTTCAACTGAACAATGATTTTGTGGTGGATGCAACCCGCAAGGGT
AACAAAATTCGTTTTGCAAAATCATTGCGTAAATCCAAACTGCTATGCAAAAGTTATGATGGTTAACGGTG
ATCACAGGATAGGTATTTTGGCAAGAGAGCCATCCAGACTGGCGAAGAGCTGTTTTTGTATTACAGATA
CAGCCAGGCTGATGCCCTGAAGTATGTCGGCATCGAAAGAGAAATGGAAATCCCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234734 representing NM_001203249
Red=Cloning site Green=Tags(s)

MGQTGKKSEKGPVCWRKRVKSEYMRLRQLKRFRADEVKSMFSSNRQKILERTEILNQEWKQRRIQPVHI
 LTSCSVTSDLDFPTQVIPLKTLNAVASVPIMYSWSPLQQNFMVEDETVLHNIPYMGDEVLDQDGTFFI EEL
 IKNYDGKLVHGDRECGFINDEIFVELVNALGQYND DDDDDDDDDPEEREKQKDL EDHRDDKESRPPRFKFP
 SDKIFEAISSMFPDKGTAEELKEKYKELTEQQLPGALPPECTPNIDGPNKSVQREQLHSFHTLFCRRC
 FKYDCFLHPFHATPNTYKRKNTETALDNKPCGPQCYQHLEGAKEFAAALT AERIKTPPKRPGRRRRLP
 NNSSRPSTPTINVLESKDTSDREAGTETGGENNDKEEEEKDETSSSSEANSRCQTPIKMKPNIEPPEN
 VEWSGAEASMFVRLIGTYDNFCAIARLIGTKTRQVYEFVYKESIIAPAPAEDVTPPRKKRKHRLW
 AAHCRKIQLKKGQNRFPGCRCKAQCNKQPCPYLAVRECDPLCLTCGAADHWDSKNVSKNCISIQRGSK
 KHL L LAPSDVAGWGI FIKDPVQKNEFISEYCGEII SQDEADRRGKVYDKYMCSFLFNLNDFVVDATRKG
 NKIRFANHSVNPNCYAKVMMVNGDHRIGIFAKRAIQTGEELFFDYRYSQADALKYVGIEREMEIP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001203249

ORF Size: 2085 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001203249.2](#)

RefSeq Size: 2682 bp

RefSeq ORF: 2088 bp

Locus ID: 2146

UniProt ID: [Q15910](#)


Cytogenetics: 7q36.1

Protein Families: Druggable Genome, Transcription Factors

MW: 80.1 kDa

Gene Summary: This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric protein complexes, which are involved in maintaining the transcriptional repressive state of genes over successive cell generations. This protein associates with the embryonic ectoderm development protein, the VAV1 oncoprotein, and the X-linked nuclear protein. This protein may play a role in the hematopoietic and central nervous systems. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Feb 2011]

Product images:



Circular map for RC234734