

Product datasheet for **RC234742**

EPB41L1 (NM_001258330) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EPB41L1 (NM_001258330) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	EPB41L1
Synonyms:	4.1N; MRD11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234742 representing NM_001258330
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTTTTCTGGGGAGGATCAATGAGGTAGAACCTGCAAAGGGCTTAGCAGAGAGCCTAGCACCTACTG
 AGCGCTCAGTAAAGAGCCTAGACATGGAGGAGAAGGACTACAGTGAGCCGATGGCCTTTTCGGAGAGGAC
 CACGCCAGCAAGGCCAGAAATCGCCCCAGAAGATTGCCAAGAAATACAAGAGTGCCATCTGCCGGGTC
 ACTCTGCTTATGCTCGGAGTATGAGTGTGAGGTGGAGAAACATGGCCGGGGCCAGGTGCTGTTTGACC
 TGGTCTGTGAACACCTCAACCTCTAGAGAAGGACTACTTCGGCCTGACCTTCTGTGATGCTGACAGCCA
 GAAGAAGTGGCTGGACCCCTCCAAGGAGATCAAGAAGCAGATCCGGAGTAGCCCCGGAATTTGCCTTC
 ACAGTCAAGTTCTACCCGCTGATCCTGCCAGCTGACAGAAGACATCACAAGATACTACCTGTGCCTGC
 AGCTGCGGGCAGACATCATCACGGGCCGCTGCCATGCTCCTTTGTACGCATGCCCTACTGGGCTCCTA
 CGCTGTGCAGGCTGAGCTGGGTGACTATGATGCTGAGGAGCATGTGGCAACTATGTCAGCGAGCTCCGC
 TTCGCCCTAACCAGACCCGGGAGCTGGAGGAGAGGATCATGGAGCTGCATAAGACATATAGGGGGATGA
 CCCCAGGAGAAGCAGAAATCCACTTCTTAGAGAATGCCAAGAAGCTTTCATGTACGGAGTAGACCTGCA
 CCATGCCAAGGACTCTGAGGGCATCGACATCATGTTAGGCGTTTGTGCCAATGGCCTGCTCATCTACCGG
 GACCGGCTGAGAAATCAACCGCTTTGCTGGCCCAAGATCCTCAAGATCTCCTACAAGAGGAGTAACCTCT
 ATATCAAGATCCGGCCTGGGAGTATGAGCAATTTGAGAGCACAATTGGCTTTAAGCTCCCAAACCCCG
 GTCAGCCAAGAGACTGTGGAAGGTCTGCATCGAGCATCATACATTCTCCGGCTGGTGTCCCCTGAGCCC
 CCACCAAGGGCTTCTGGTGTGGCTCCAAGTTCGGTACAGTGGGAGGACCCAGGCACAGACTCGCC
 AGGCCAGCGCCCTCATTGACCGGCTGCACCCTCTTTGAGCGTCTTCCAGCAAACGGTACACCATGTC
 CCGCAGCCTTGATGGAGCAGAGTTCTCCCGCCAGCCTCGGTACGCGAGAACCATGATGCAGGGCCTGAC
 GGTGACAAGCGGGATGAGGATGGCGAGTCTGGGGGCAACGGTACAGGGCTGAGGAGGGAGAGGTCAGGA
 CTCCAACCAAGATCAAGGAGCTAAAGTCTTAGACAAGCCAGAAGATGTCTTGCTGAAGCACCAGGCCAG
 CATCAATGAGCTCAAAAGGACCTGAAGGAGCCCAACAGCAAACCTCATCCACCGGATCGAGACTGGGAA
 CGGGAGCGCAGGCTGCCCTCCTCCCCGCCTCCCCCTCCCCAAGGGCACCCCTGAGAAAGCCAATGAGC
 CCGTGAAAACAGAAACCATGACTGTCAGCAGTCTGGCCATTAGAAAGAAGATTGAGCCGGAGGCGTACT
 GCAGACCAGAGTCTCCGCTATGGATAACACCCAGCAGTTGATGGGAGTGCCTCAGTGGGAGGGAGTTC
 ATAGCAACCACTCCCTCCATCACACGGAGACCATATCGACCACCATGGAGAACAGTCTCAAGTCCGGGA
 AGGGGGCAGCTGCCATGATCCCAGGCCACAGACGGTGGCCACGGAAATCCGTTCTTTCTCCGATCAT
 CGGGAAAGATGTCTCACCAGCACCTACGGCGCCACTGCGGAAACCCCTCTCAACCTCCACCACCACCCAT
 GTCACCAAACTGTGAAAGGAGGGTTTTCTGAGACAAGGATCGAGAAGCGAATCATCATTACTGGGGATG
 AAGATGTCGATCAAGACCAGGCCCTGGCTTTGGCCATCAAGGAGGCCAACTGCAGCATCCTGATATGCT
 GGTAACCAAGCTGTCGTATACAGAGAAACAGACCCATCCCCAGAGGAGAGGGACAAGAAGCCACAGGAA
 TCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234742 representing NM_001258330
Red=Cloning site Green=Tags(s)


```
MVFLGRINEVEPAKGLAESLAPTERSVKSLDMEEKDYSEADGLSERTTPSKAQKSPQKI AKKYKSAICRV
TLLDASEYECEVEKHGRGQVLFDLVCEHLNLEKDYFGLTFCDADSQKNWLDPSKEIKKQIRSSPWNFAF
TVKFYPPDPAQLTEDITRYYLCLQLRADIIITGRLPCSFVTHALLGSYAVQAE LGDYDAEEHVGNVSELR
FAPNQTRLEERIMELHKTYRGMTPGEAEIHFL ENAKKLSMYGVDLHHA KDSEGIDIMLGVCANGLLIYR
DRLRINRF AWPKILKISYKRSNFYIKIRPGEYEQFESTIGFKLPNHRS AKRLWKVCIEHHTFFRLVSEPEP
PPKGF LVMGSKFRYSGR TQAQTRQASALIDRPAPFFERSSSSKRYTMSRSLDGA EFSRPASVSENHDAGPD
GDKRDEGEGSGQRSEAE EGEVRTPTKIKELKFLDKPEDVLLKHQASINELKRTLKEPN SKLIHRDRDWE
RERRLPSSPASPSPKGTPEKANEPVKTETMTVSSLAIRKKIEPEAVLQTRVSAMDNTQQVDGSASV GREF
IATTPSITTETISTTMENSLKSGKGAAMIPGPQTVATEIRSLSPIIGKDVLTSTYGATAETLSTSTTTT
VTKTVKGGFSETRIEKRIIITGDEVDQDQALALAIKEAKLQHPDMLVTKAVVYRETDPSPEERDKK PQE
S
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001258330

ORF Size: 2103 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001258330.1](#), [NP_001245259.1](#)

RefSeq Size: 5745 bp

RefSeq ORF: 2106 bp

Locus ID: 2036

UniProt ID: [Q9H4G0](#)

Cytogenetics: 20q11.23


Protein Families: Druggable Genome

Protein Pathways: Tight junction

MW: 79.3 kDa

Gene Summary: Erythrocyte membrane protein band 4.1 (EPB41) is a multifunctional protein that mediates interactions between the erythrocyte cytoskeleton and the overlying plasma membrane. The encoded protein binds and stabilizes D2 and D3 dopamine receptors at the neuronal plasma membrane. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2015]

Product images:



Circular map for RC234742